

MISSOURI DEPARTMENT OF HEALTH AND SENIOR SERVICES

STATE PUBLIC HEALTH LABORATORY

BREATH ALCOHOL PROGRAM

INTOX EC/IR II MAINTENANC			REPORT #3
Complete this report at the time of the regular monthly preventive maintenance check (not to exceed 35			
days). Complete this report whenever the instrument is serviced or repaired and whenever it is placed			
into service. Retain the original and send a copy w	within 15 days to the		am, DHSS.
INTOX EC/IR II SN NAME OF AGENCY		DATE OF INSPECTION	
12682 Jefferson Cou	unty Sheriff	08/07/2020	
LOCATION OF INSTRUMENT (STREET AND CITY)		TIME OF INSPECTION	
34 Dillon Plaza High Ridge, MO 63049 NZ		18:35 CDT	
CHECKLIST: Place a mark in the box by each item if found to be satisfactory or is operating within			
established limits. (Write in observed values where determined). Unmarked items must be corrected			
before using instrument.			
X DIAGNOSTIC RECORD			
X BLANK CHECK X CO2 CHECK			
X FC 1 TEMP X FLOW CHECK			
X SRC TEMP X FCB CHECK			
X DET TEMP X CRC COMP CHECK			
X BT TEMP X CRC CAL CHECK			
X STD 2 TEMP X PRINT TEST			
X ETH CHECK			
BREATH ANALYZER ACCURACY STANDARDS			
SIMULATOR SOLUTION			
X STANDARD SUPPLIER Intoximeters LOT# AG912204		EXP. DATE 05/02/2021	
SIMULATOR TEMP (34°C ±0.2°C) SIM	I. SN	SIM. NIST EXP DAT	ſE
—			
CALIBRATION CHECK - (ONLY ONE STANDARD IS TO BE USED PER MAINTENANCE REPORT) Run three tests using a standard solution. All three tests must be within ±5% of the standard value and must have a spread of .005 or less. Mark the box corresponding to the standard solution being used. X 0.10% STANDARD - MUST READ BETWEEN 0.095% AND 0.105% INCLUSIVE 0.08% STANDARD - MUST READ BETWEEN 0.076% AND 0.084% INCLUSIVE 0.04% STANDARD - MUST READ BETWEEN 0.038% AND 0.042% INCLUSIVE			
TEST 1 * 0.099 g/210L TEST 2 * 0.0	99 g/210L	TEST 3 🔗 0.099 g	J/210L
INDICATE THE NUMBER OF BREATH TESTS IN THE FOLLOWING RANGES SINCE THE LAST MAINTENANCE REPORT:			JANCE REPORT:
REFUSALS 1 004 0 .0509 0	.1014 0	.1519 0 C	OVER .19 0
LIST ANY NEW PARTS AND DESCRIBE ANY ALTERATION OR MODIFICATION THAT WAS MADE TO RESTORE THE INSTRUMENT TO OPERATE SATISFACTORILY AND WITHIN ESTABLISHED LIMITS (USE OTHER SIDE IF NECESSARY).			
INSPECTING OFFICER			
SIGNATURE	PRINT FULL NAME		
► Der 1 ~ 7 26	Der Control 26 ALEXANDER KAUSLER		
TYPE II PERMIT NUMBER EXPIRATION DATE TELEPHONE NUMBER			
200129 03/09/2022	(636)797-5000)	
RETURN COMPLETED REPORT TO THE:			
Breath Alcohol Program, Missouri Department of Health and Senior Services,			
	IL OL HEALLII dIQ	Sentor Services,	
by mail, fax, or e-mail			