



MISSOURI DEPARTMENT OF HEALTH AND SENIOR SERVICES
STATE PUBLIC HEALTH LABORATORY

ALCO-SENSOR IV WITH PRINTER MAINTENANCE REPORT

REPORT #7

Complete this report in duplicate at the time of the regular monthly preventative maintenance check, and whenever instrument is repaired. Send copy to Department of Health and Senior Services; retain original in department file.

ALCO SENSOR IV SN 087970	PRINTER SN 08C.3527.187	DATE OF INSPECTION 06/24/2020
-----------------------------	----------------------------	----------------------------------

LOCATION OF INSTRUMENT (STREET AND CITY) 9701 Marion Park Drive, Kanas City, MO 64137	TIME OF INSPECTION 10:35 am
--	--------------------------------

CHECKLIST: Place a mark in the box by each item if found to be satisfactory or if operating within established limits. (Write in observed values where determined.) Unmarked items must be corrected before using instrument.

- DIGITAL READOUT (ALL ELEMENTS OPERATIONAL)
- TEMPERATURE OF ALCO SENSOR (10°C - 40°C)
- PRINTER WORKING PROPERLY
- TIME AND DATE DISPLAYING PROPERLY

BREATH ALCOHOL ACCURACY STANDARDS

<input type="checkbox"/> SIMULATOR SOLUTION	<input checked="" type="checkbox"/> COMPRESSED ETHANOL-GAS MIXTURE
<input checked="" type="checkbox"/> STANDARD SUPPLIER <u>Intoximateres</u> LOT # <u>AG905605</u> EXP. DATE <u>02/25/2021</u>	
<input type="checkbox"/> SIMULATOR TEMPERATURE (34°C ± 0.2°C) _____ SIMULATOR SN _____ SIMULATOR EXP DATE _____	

CALIBRATION CHECK – (ONLY ONE STANDARD IS TO BE USED PER MAINTENANCE REPORT)
Run three tests using a standard solution. All three tests must be within ±5% of the standard value and must have a spread of .005 or less. Check the box corresponding to the standard solution being used. (PRINTOUT ATTACHED)

- 0.100% STANDARD - MUST READ BETWEEN 0.095% and 0.105% INCLUSIVE
- 0.080% STANDARD - MUST READ BETWEEN 0.076% and 0.084% INCLUSIVE
- 0.040% STANDARD - MUST READ BETWEEN 0.038% and 0.042% INCLUSIVE

TEST 1 0.079	TEST 2 0.079	TEST 3 0.079
---------------	---------------	---------------

RFI DETECTOR OPERATING

INDICATE THE NUMBER OF BREATH TESTS IN THE FOLLOWING RANGES SINCE THE LAST MAINTENANCE REPORT: (DO NOT INCLUDE SELF-ADMINISTERED TESTS)

REFUSALS	(0-.04)	(.05-.09)	(.10-.14)	(.15-.19)	(OVER .19)
----------	---------	-----------	-----------	-----------	------------

List any new parts and describe any alteration or modification that was made to restore the instrument to operate satisfactorily and within established limits (use other side if necessary).

INSPECTING OFFICER	
SIGNATURE 	PRINT NAME Douglas Davidson 5646
TYPE II PERMIT NUMBER/EXPIRATION DATE 290087/04-22-2021	TELEPHONE NUMBER (816) 234-5000

Return completed report to the: Breath Alcohol Program, MO Department of Health and Senior Services, Southeast District Office
2875 James Boulevard
Poplar Bluff, MO 63901

..

AS IV Serial no: 087970
Version no: 532B

TEST RECORD 00297

Temp	Date	Time	s/ 210L
Air Blank:			
	06/24/20	22:40	.000
Calibration Check:			
	30 06/24/20	22:40	.079

Air Blank:
06/24/20 22:40 .000
Calibration Check:
30 06/24/20 22:40 .079

Subject Name

Test 1

Subject I.D.

Operator Name, I.D.

Davidson 5646

Location

290087 04/22/2021

..

AS IV Serial no: 087970
Version no: 532B

TEST RECORD 00298

Temp	Date	Time	s/ 210L
Air Blank:			
	06/24/20	22:42	.000
Calibration Check:			
	31 06/24/20	22:42	.079

Air Blank:
06/24/20 22:42 .000
Calibration Check:
31 06/24/20 22:42 .079

Subject Name

Test 2

Subject I.D.

Operator Name, I.D.

Davidson 5646

Location

290087 04/22/2021

..

AS IV Serial no: 087970
Version no: 532B

TEST RECORD 00299

Temp	Date	Time	s/ 210L
Air Blank:			
	06/24/20	22:44	.000
Calibration Check:			
	31 06/24/20	22:44	.079

Air Blank:
06/24/20 22:44 .000
Calibration Check:
31 06/24/20 22:44 .079

Subject Name

Test 3

Subject I.D.

Operator Name, I.D.

Davidson 5646

Location

290087 04/22/2021

..

AS IV Serial no: 087970
Version no: 532B

TEST RECORD 00300

Temp	Date	Time	s/ 210L
VOID: RFI			
	12 06/24/20	22:46	

VOID: RFI

12 06/24/20 22:46

Subject Name

RFI Test

Subject I.D.

Operator Name, I.D.

Davidson 5646

Location

290087 04/22/2021

MISSOURI DEPARTMENT OF HEALTH AND SENIOR SERVICES
**BLOOD ALCOHOL TEST REPORT - ALCO-SENSOR IV
 WITH PRINTER**

Case Number: _____

FORM #8

SUBJECT'S NAME	DATE OF TEST
----------------	--------------

OPERATIONAL CHECKLIST: ALCO-SENSOR IV WITH PRINTER

ALCO-SENSOR SERIAL NO. 08C.3527.187	PRINTER SERIAL NO. 0807970	LOCATION OF INSTRUMENT
---	--------------------------------------	------------------------

- 1. Examination of mouth conducted. If any substance is observed or indicated to be present, the substance observed or indicated must be removed prior to starting the 15 minute observation period.
- 2. Subject observed for at least 15 minutes by PO Davidson 5646.
No smoking, oral intake or vomiting during this time; if vomiting occurs, start over with 15 minute observation period.
- 3. Make sure printer is connected to Alco-Sensor IV.
- 4. Turn printer on.
- 5. Insert mouthpiece into Alco-Sensor IV.
- 6. Observe temperature display, make sure temperature reading is between 10°C and 40°C.
- 7. When "BLNK" is displayed on Alco-Sensor IV, air blank is taken.
- 8. When "TEST" is displayed on Alco-Sensor IV, take subject breath sample.
- 9. When "SET" is displayed on Alco-Sensor IV, press SET button.
- 10. When printer has completed printing test result, tear off tape and fill in subject and officer information.
- 11. Press red button to eject mouthpiece.
- 12. Turn printer off.
- 13. Attach printout to this report.

CERTIFICATION BY OPERATOR	BAC
----------------------------------	-----

As set forth in the rules promulgated by the Department of Health and Senior Services related to the determination of blood alcohol by breath analysis, I certify that:

- 1. There was no deviation from the procedure approved by the department.
- 2. To the best of my knowledge the instrument was functioning properly.
- 3. I am authorized to operate the instrument.
- 4. No radio transmission occurred inside the room where and when this was being conducted.

NAME OF OPERATOR Davidson	PERMIT NO. 290087	EXPIRATION DATE 04/22/2021
WITNESS (IF ANY)		DATE





STATE OF MISSOURI
DEPARTMENT OF HEALTH AND SENIOR SERVICES
BREATH ALCOHOL PROGRAM

PERMIT
TYPE II

DOUGLAS DAVIDSON

is hereby authorized to instruct and supervise operators, train instructors, inspect, calibrate, perform field service and repairs, and operate the following breath analyzer(s):

ALCO-SENSOR IV WITH PRINTER, INTOXILYZER 8000

for the determination of the alcoholic content of blood from a sample of expired air. Permit issued under the provisions of sections 577.020 through 577.041, RSMo and 306.111 through 306.119 RSMo.

DATE 4/22/2019
NUMBER 290087
EXPIRES 4/22/2021

Davidson
DIRECTOR OF STATE PUBLIC HEALTH LABORATORY

Davidson
DIRECTOR OF DEPARTMENT OF HEALTH AND SENIOR SERVICES
LAB-4 (RE-10)

STATE OF MISSOURI
DEPARTMENT OF HEALTH AND SENIOR SERVICES
BREATH ALCOHOL PROGRAM
INSTRUMENT OPERATOR CARD

This named operator is authorized to operate an alcohol breath alcohol instrument for the determination of the alcoholic content in breath (over a specified air volume).

Operator: **DAVIDSON, DOUGLAS**
Permit No. **290087**
Date Issued **4/22/2019** Date Expires **4/22/2021**

2



7 Eastgate Dr. • P.O. Box 790 • Jacksonville, IL 62651-0790
217-245-2183 • Fax: 217-243-7634 • www.ilmoproducts.com

Certificate of Analysis

Certificate ID: 11721
Part #: BAC105L080T
Cylinder Size: 105L
Lot Number: 09919880A1
Expiration: 3/5/2021

0.080 BAC (For the calibration of instruments used to determine breath alcohol concentration)

Contents: 105 Liters @ 1000 psig 70°F (21°C)
Analytical Accuracy (U, L): 1/-0.002 BAC (6/24/01) NDIR
Method: NDIR
Distributed by: Cyl Inc.
316 East Ninth Street
Owensboro, KY 42303
Phone 866-835-0690
www.alcoholtest.com

*NIST Traceable Reference Material
Cylinder No. CC274523 / Job No. 09160306
Certified 362.2 μmol/mol Ethanol in Nitrogen

Store in dry area, away from sources of heat,
ignition and direct sunlight. Do not allow storage
area to exceed 52 °C (125 °F).

Specialty Gas Lab Tech
Specialty Gas Lab Tech
Date 01-24-19



The calibration results within this certificate were obtained using equipment and methods certified by the National Institute of Standards and Technology (NIST) and apply only to the items contained on this certificate. ILMO Products Company makes no warranty or representation as to the accuracy of the information provided for any particular purpose. The information is at the sole discretion and risk of the user. Liability shall be limited to established replacement cost of this material or service.

ISO/IEC 17025:2005 Accredited Laboratory