

MISSOURI DEPARTMENT OF HEALTH AND SENIOR SERVICES STATE PUBLIC HEALTH LABORATORY ALCO-SENSOR IV WITH PRINTER MAINTENANCE REPORT

RECEIVED By Stephen Wilson at 9:34 am, Apr 13, 2020 #7

Complete this report in d Send copy to Department	•				ever instrument is repaired.	
ALCO SENSOR IV SN		PRINTER SN		DATE O	F INSPECTION	
LOCATION OF INSTRUMENT (TIME OF	INSPECTION		
CHECKLIST: Place a ma ues where determined.)				ing within established lin	nits. (Write in observed val-	
DIGITAL READOUT (ALL ELEMENTS OPERATIONAL)						
TEMPERATURE OF ALCO SENSOR (10°C - 40°C)						
BREATH ALCOHOL AC	CURACY STANDA	RDS				
	ATOR SOLUTION					
STANDARD SUPPLIER		L	LOT # EXP. DATE			
SIMULATOR TEMPERATURE (34°C ± 0		0.2°C) SIM	ULATOR SN	SIMULATOR EXP DATE		
Run three tests using a standard solution. All three tests must be within ±5% of the standard value and must have a spread of .005 or less. Check the box corresponding to the standard solution being used. (PRINTOUT ATTACHED) 0.100% STANDARD - MUST READ BETWEEN 0.095% and 0.105% INCLUSIVE 0.080% STANDARD - MUST READ BETWEEN 0.076% and 0.084% INCLUSIVE 0.040% STANDARD - MUST READ BETWEEN 0.038% and 0.042% INCLUSIVE						
TEST 1 🖝		TEST 2 🖝		TEST 3 🖝		
INDICATE THE NUMBER OF BREATH TESTS IN THE FOLLOWING RANGES SINCE THE LAST MAINTENANCE REPORT: (DO NOT INCLUDE SELF-ADMINISTERED TESTS)						
REFUSALS	(004)	(.0509)	(.1014)	(.1519)	(OVER .19)	
List any new parts and describe any alteration or modification that was made to restore the instrument to operate satisfactorily and within established limits (use other side if necessary).						
INSPECTING OFFICER						
SIGNATURE Patn. Edu	die Lee 311			PRINT NAME		
TYPE II PERMIT NUMBER/EXPIRAT				TELEPHONE NUMBER		
Return completed report to the: Breath Alcohol Program, MO Department of Health and Senior Services, Southeast District Office 2875 James Boulevard Poplar Bluff, MO 63901						



CERTIFICATE OF ANALYSIS

Certified Alcohol Reference Solution for Simulator

Random Samples of Lot Number 19160 of Alcohol Reference Solution for Simulator were analyzed by gas chromatography on July 10, 2019, using a Perkin Elmer Gas Chromatograph Autosystem XL S/N: 610N9030209, and found to contain 0.1202% (w/vol) ethyl alcohol. The expiration date for this lot number is July 9, 2021 at 11:59 PM.

When used in a calibrated Simulator, operating at $34^{\circ}C$ +/- .2°C, this solution will give a breath alcohol analysis instrument reading of 0.100 g/210L +/- 3%.

The alcohol and water used in this solution were free of test interfering substances.

Ted L. Pauley, President GUTH LABORATORIES, INC.

NIST Traceability:

Testing was conducted using Cerilliant Reference Standard lot number FN04271602 whose values are traceable to NIST. All balances are calibrated annually by an outside agency using NIST traceable weights. Calibration verification is done prior to each use utilizing NIST traceable weights.

Deeratov Name, I.D. Subject

Location

Location

Or or Contor Name: S. Ň) eend B Secol Secol Secol

Subject I. D.

Subject Name

6 ş

Subject Have

VUID: RFI 12 04/12/20 21:50

Temp

naj e

Time 210L

Air Blank: 04/12/20 21:49 .000 Calibration Check: 23 04/12/20 21:49 .101

Temp s/ Date Time 210L

TIST RECORD 00228

AS IV Serial no: 111640 Version no: 532B

TIST RECORD 00229

AS IU Serial no: 111640 Version no: 532B

subject , Operator Names Location)ard 5 1987 1987 1987 5 I. D. ~

Subject Name

Air Blank: 04/12/20 21:44 .000 Calibration Check: 20 04/12/20 21:44 .102

Temp Date Time 2191

TEST RECORD 00225 Ŵ.,

AS IU Serial no: 111640 Version no: 532B

Subject I.D. 2

Operator Name, 1

(]_D=

Location

Subject Name

Air Blank: 04/12/20 21:48 .000 Calibration Check: 22 04/12/20 21:48 .101

Temp Date Time 210L

TEST RECORD 00227

AS IU Serial no: 111640 Version no: 532B



STATE OF MISSOURI DEPARTMENT OF HEALTH AND SENIOR SERVICES BREATH ALCOHOL PROGRAM



PERMIT

EDDIE LEE

is hereby authorized to instruct and supervise operators, train instructors, inspect, calibrate, perform field service and repairs, and operate the following breath analyzer(s):

ALCO-SENSOR IV WITH PRINTER

for the determination of the alcoholic content of blood from a sample of expired air. Permit issued under the provisions of sections 577.020 through 577.041, RSMo and 306.111 through 306.119 RSMo.

DATE 5/22/2018

NUMBER 280192

EXPIRES 5/22/2020

MO 580-0771 (6-10)

113 horas

DIRECTOR OF STATE PUBLIC HEALTH LABORATORY

UUM

DIRECTOR OF DEPARTMENT OF HEALTH AND SENIOR SERVICES LAB-4 (R6-10)

