

MISSOURI DEPARTMENT OF HEALTH AND SENIOR SERVICES STATE PUBLIC HEALTH LABORATORY

RECEIVED By Brian Lutmer at 9:57 am, Nov 04, 2024

ALCO-SENSOR IV WITH PRINTER MAINTENANCE REPORT

REPORT #7

Complete this report in duplicate at the time of Send copy to Department of Health and Senior				never instrument is repaired.	
ALCO SENSOR IV_SN 097415	NAME OF AGENCY Missouri State Hi	ghway Patrol	il.	of Inspection - 62 - 2024	
LOCATION OF INSTRUMENT (STREET AND CITY) Zone 4/15 Office 275-B Southwest outer Road, Branson				9 am	
CHECKLIST: Place a mark in the box by each item if found to be satisfactory or if operating within established limits. (Write in observed values					
where determined.) Unmarked items must be corrected before using instrument.					
DIGITAL READOUT (ALL ELEMENTS OPERATIONAL)					
✓ TEMPERATURE OF ALCO SENSOR (10°C - 40°C)					
✓ PRINTER WORKING PROPERLY					
✓ TIME AND DATE DISPLAYING PROPERLY					
BREATH ALCOHOL ACCURACY STANDARD	5				
✓ SIMULATOR SOLUTION ☐ COMPRESSED ETHANOL-GAS MIXTURE					
STANDARD SUPPLIER Guth Laboratories LOT # 23390 EXP. DATE 10-17-2025					
☑ SIMULATOR TEMPERATURE (34°C ± 0.2°	°C) 33.99 SIM	M. SNMP2422	SIM. NIST E	EXP DATE 12-05-2024	
 ✓ CALIBRATION CHECK – (ONLY ONE STANDARD IS TO BE USED PER MAINTENANCE REPORT) Run three tests using a standard solution. All three tests must be within ±5% of the standard value and must have a spread of .005 or less. Check the box corresponding to the standard solution being used. (PRINTOUT ATTACHED) ✓ 0.100% STANDARD - MUST READ BETWEEN 0.095% and 0.105% INCLUSIVE ○ 0.080% STANDARD - MUST READ BETWEEN 0.076% and 0.084% INCLUSIVE ○ 0.040% STANDARD - MUST READ BETWEEN 0.038% and 0.042% INCLUSIVE 					
TEST 1 .099 TE	TEST 2 ▼ .099		TEST 3 .098 .		
✓ RFI DETECTOR OPERATING					
INDICATE THE NUMBER OF BREATH TESTS IN THE FOLLOWING RANGES SINCE THE LAST MAINTENANCE REPORT: (DO NOT INCLUDE SELF-ADMINISTERED TESTS)					
REFUSALS 0 (004) 0	(.0509) 0	(.1014) 0	(.1519) 0	(OVER .19) 0	
List any new parts and describe any alteration established limits (use other side if necessary).	or modification that w	vas made to restore	the instrument to ope	erate satisfactorily and within	
INSPECTING OFFICER					
SIGNATURE .			Brandon Hall		
230011 / 01-24-2025			(417) 895-6868		
Return completed report to the: Breath Alco		partment of Health ar	nd Senior Services, S	outheast District Office	

05 IV Serial no: 857415 TEST RECORD 00991 9/ Temp Date Time 210L Air Blank: 11/02/24 10:29 .000 Calibration Check: 19 11/02/24 10:29 .099 Subject Name Acc check Subject I.D. Brandon Hall 347 Operator Name, I.D. Zone 4/15 office Location Branson

AS IV Serial no: 097415 Version no: 532B TEST RECORD 00992 Time 210L Date Temp Air Blank: 11/02/24 10:32 .000 Calibration Check: 20 11/02/24 10:32 .099 Subject Name Acc Check Subject I.D. Brandon Hall 347 Operator Name, I.D. Zone 4/15 office Location Branson

AS IV Serial no: 097415 Version no: 532B TEST RECORD 00993 Temp Time 210L Date Air Blank: 11/02/24 10:35 .000 Calibration Check: 21 11/02/24 10:35 .098 Subject Name Acc check Subject I.D. Brandon Hall 347 Operator Name, I.D. Zone 4/15 office Location

Branson

Version no: 532B

TEST RECORD 00994
9/
Temp Date Time 210L

VOID: RFI
12 11/02/24 10:39

Subject Name
Acc check

Subject I.D.

Grandon Hall 347

Operator Name, I.D.

Zone 4//5 office

Location

Branson

AS IV Serial no: 097415



CERTIFICATE OF ANALYSIS

Certified Alcohol Reference Solution for Simulator

Random Samples of Lot Number 23390 of Alcohol Reference Solution for Simulator were analyzed by gas chromatography on October 18, 2023, using a Perkin Elmer Gas Chromatograph Autosystem XL S/N: 610N9030209, and found to contain 0.1207% (w/vol) ethyl alcohol. The expiration date for this lot number is October 17, 2025 at 11:59 PM.

When used in a calibrated Simulator, operating at 34°C +/- .2°C, this solution will give a breath alcohol analysis instrument reading of 0.100 g/210L +/- 3%.

The alcohol and water used in this solution were free of test interfering substances.

Ted L. Pauley, President GUTH LABORATORIES, INC.

NIST Traceability:

Testing was conducted using Cerilliant Reference Standard lot number FN03072301 whose values are traceable to NIST.

All balances are calibrated annually by an outside agency using NIST traceable weights. Calibration verification is done prior to each use utilizing NIST traceable weights.



STATE OF MISSOURI

DEPARTMENT OF HEALTH AND SENIOR SERVICES
BREATH ALCOHOL PROGRAM



PERMIT TYPE II

BRANDON A. HALL

is hereby authorized to instruct and supervise operators, train instructors, inspect, calibrate, perform field service and repairs, and operate the following breath analyzer(s):

ALCO-SENSOR IV WITH PRINTER, INTOX DMT

for the determination of the alcoholic content of blood from a sample of expired air. Permit issued under the provisions of sections 577.020 through 577.041, RSMo and 306.111 through 306.119 RSMo.

DATE	1/24/2023	Mike Massim		
5/112		DIRECTOR OF STATE PUBLIC HEALTH LABORATORY		
NUMBER	230011			
EXPIRES	1/24/2025	Davla I. Nichelson		
		DIRECTOR OF DEPARTMENT OF HEALTH AND SENIOR SERVICES		

MO 580-0771 (6-10)

LAB-4 (R6-10)



The named cardholder is authorized to operate an evidential breath alcohol instrument for the determination of the alcoholic content in breath form of expired air in Missouri.

Operator HALL, BRANDON

Permit No 230011

