

MISSOURI DEPARTMENT OF HEALTH AND SENIOR SERVICES STATE PUBLIC HEALTH LABORATORY

RECEIVED By Brian Lutmer at 9:57 am, Nov 04, 2024

ALCO-SENSOR IV WITH PRINTER MAINTENANCE REPORT

REPORT #7

Complete this report in duplicate at the time of Send copy to Department of Health and Senior				ce check, an	d whene	ever instrument is	s repaired.		
ALCO SENSOR IV SN 105447	NAME OF AGENCY Missouri State Highway Patrol				DATE OF INSPECTION 11/01/2024				
LOCATION OF INSTRUMENT (STREET AND CITY) 1000 Osage Beach Parkway, Osage Beach, Mo					TIME OF INSPECTION 11:32 pm				
CHECKLIST: Place a mark in the box by each item if found to be satisfactory or if operating within established limits. (Write in observed values									
where determined.) Unmarked items must be corrected before using instrument.									
DIGITAL READOUT (ALL ELEMENTS OPERATIONAL)									
✓ TEMPERATURE OF ALCO SENSOR (10°C - 40°C)									
PRINTER WORKING PROPERLY									
✓ TIME AND DATE DISPLAYING PROPERLY									
BREATH ALCOHOL ACCURACY STANDAR	RDS								
☑ SÍMULATOR SOLUTION		COMP	RESSED	ETHANOL-G	GAS MIX	TURE			
STANDARD SUPPLIER Guth Laboratories LOT # 23390 EXP. DATE 10/17/2025									
✓ SIMULATOR TEMPERATURE (34°C ± 0	.2°C) 34.00 s	IM. SN	MP2315	SIM.	NIST EX	(P DATE 11/20/	2024		
Run three tests using a standard solution less. Check the box corresponding to the . 0.100% STANDARD - MUST READ 0.080% STANDARD - MUST READ 0.040% STANDARD - MUST READ	standard solution bein BETWEEN 0.095% ar BETWEEN 0.076% är	g used. (PRIN nd 0.105% INC nd 0.084% INC	TOUT AT LUSIVE LUSIVE						
TEST 1 ▼ .098	TEST 2 .097			TEST 3 .097					
☑ RFI DETECTOR OPERATING					,				
INDICATE THE NUMBER OF BREATH TES (DO NOT INCLUDE SELF-ADMINISTERED		NG RANGES			INTENA	NCE REPORT:			
REFUSALS 0 (004) 0	(.0509)	(.1014)	0	(.1519)	0	(OVER .19)	0		
List any new parts and describe any alteration established limits (use other side if necessary All okay		was made to	restore th	e instrument	to opera	ate satisfactorily	and within		
INSPECTING OFFICER		22900		RINT NAME					
SIGNATURE					J.C. Simpson				
TYPE I PERMIT NUMBER/EXPIRATION DATE 230,135 - 07/06/2025	1	(573) 751-1000							
Return completed report to the: Breath Al by mail, f	cohol Program, MO D	epartment of H	lealth and	Senior Serv	ices, Soi	utheast District C	Office		

AS IV Serial no: 185447 Version no: 532B TEST RECORD 01134

Temp

Date

Time 210L

Calibration Check: 21 11/01/24 23:38 .098 Air Blank: 11/01/24 23:38 .000

Subject Name

1651

Subject I.D.

Operator Name, I.D.

Locat ion IL. SIMBON BRIBE

HOO CITY PKWY,

AS IV Serial no: 105447 Version no: 532B

TEST RECORD 01135

Temp Date Time 210L

Air Blank:

11/01/24 23:42 .000 Calibration Check: 22 11/01/24 23:42 .097

Subject Name

Subject I.D. 133

Operator Name, I.D.

400 Locat ion Tic. 5 Fm Brn, 230135

> Subject Name Calibration Check: 22 11/01/24 23:45 .097 Temp Air Blank: AS IV Serial no: 185447 Version no: 532B TEST RECORD 01136 11/01/24 23:45 .000 Date Time 210L

Subject I.D.

Opérator Name, I.D.

S.C. SIMBON OSVISS

Location libr COTY PICLY

BSAGE SCHU M

AS IV Serial no: 185447 Version no: 532B

UOID: RFI 12 11/01/24 23:47

Temp

Date

Time 210L

TEST RECORD 01137

Subject Name

Subject I.D.

Operator Name, I.D.

J.C. SIMPSON 93/135 Location

poor cary bend

OSAGE BEARD MA



CERTIFICATE OF ANALYSIS

Certified Alcohol Reference Solution for Simulator

Random Samples of Lot Number 23390 of Alcohol Reference Solution for Simulator were analyzed by gas chromatography on October 18, 2023, using a Perkin Elmer Gas Chromatograph Autosystem XL S/N: 610N9030209, and found to contain 0.1207% (w/vol) ethyl alcohol. The expiration date for this lot number is October 17, 2025 at 11:59 PM.

When used in a calibrated Simulator, operating at 34°C +/- .2°C, this solution will give a breath alcohol analysis instrument reading of 0.100 g/210L +/- 3%.

The alcohol and water used in this solution were free of test interfering substances.

Ted L. Pauley, President GUTH LABORATORIES, INC.

NIST Traceability:

Testing was conducted using Cerilliant Reference Standard lot number FN03072301 whose values are traceable to NIST.

All balances are calibrated annually by an outside agency using NIST traceable weights. Calibration verification is done prior to each use utilizing NIST traceable weights.



STATE OF MISSOURI

DEPARTMENT OF HEALTH AND SENIOR SERVICES
BREATH ALCOHOL PROGRAM



PERMIT TYPE II

JUSTICE SIMPSON

is hereby authorized to instruct and supervise operators, train instructors, inspect, calibrate, perform field service and repairs, and operate the following breath analyzer(s):

ALCO-SENSOR IV WITH PRINTER, INTOX DMT

for the determination of the alcoholic content of blood from a sample of expired air. Permit issued under the provisions of sections 577.020 through 577.041, RSMo and 306.111 through 306.119 RSMo.

DATE	7/6/2023		Mike Massure
DAIL			DIRECTOR OF STATE PUBLIC HEALTH LABORATORY
NUMBER	230135	 	Davla I. Nichelson
EXPIRES	7/6/2025	 	DIRECTOR OF DEPARTMENT OF HEALTH AND SENIOR SERVICES

MO 580-0771 (6-10)

LAB-4 (R6-10)



The named cardholder is authorized to operate an evidential breath alcohol instrument for the determination of the alcoholic content in breath form of expired ai in Missouri.

Operator SIMPSON, JUSTICE

Permit No 230135

Date Issued 7/6/2023 Date Expires 7/6/2025

