

By Tracy Crews at 7:25 am, Dec 17, 2024



MISSOURI DEPARTMENT OF HEALTH AND SENIOR SERVICES STATE PUBLIC HEALTH LABORATORY

ALCO-SENSOR IV WITH PRINTER MAINTENANCE REPORT

REPORT #7

T.	STATE OF THE PARTY											
Cor	nplete thi	is report in d Departmen	uplicate a t of Health	t the time o	of the regular mont or Services; retain	thly prevent original in c	ative maintena lepartment file.	nce check, a	and whenev	er instrument is	repaired.	
	O SENSOR	IV SN			NAME OF AGENC		sior Springs Departme		DATE OF	INSPECTION 12/13/2024		
107977 LOCATION OF INSTRUMENT (STREET AND CITY) 301 S Main St, Excelsior Sprii				ior Sprin					TIME OF INSPECTION 1445			
CHECKLIST: Place a mark in the box by each item if found to be satisfactory or if operating within established limits. (Write in observed values												
whe	where determined.) Unmarked items must be corrected before using instrument.											
X	M DIGITAL READOUT (ALL ELEMENTS OPERATIONAL)											
×	▼ TEMPERATURE OF ALCO SENSOR (10°C - 40°C)											
×	PRINTE	R WORKING	3 PROPE	RLY								
X	TIME AN	ND DATE DIS	SPLAYING	PROPER	LY							
BRI	EATH AL	COHOL AC	CURACY	STANDAR	DS:							
X	SIMULA	TOR SOLU	ION				COMPRESSED	ETHANOL-	GAS MIXT	VRE		
X	STANDA	ARD SUPPLI	ER C	uth		LOT#_	23390	EXP. DA	те10	0/17/2025		
×	SIMULA	TOR TEMPE	RATURE	(34°C ± 0.2	2°C) 34.0 C	SIM. SN	MP3564	SIM	. NIST EXF	DATE 05/20)/2025	
`	CALIBRATION CHECK – (ONLY ONE STANDARD IS TO BE USED PER MAINTENANCE REPORT) Run three tests using a standard solution. All three tests must be within ±5% of the standard value and must have a spread of .005 or less. Check the box corresponding to the standard solution being used. (PRINTOUT ATTACHED) 0.100% STANDARD - MUST READ BETWEEN 0.095% and 0.105% INCLUSIVE 0.080% STANDARD - MUST READ BETWEEN 0.076% and 0.084% INCLUSIVE 0.040% STANDARD - MUST READ BETWEEN 0.038% and 0.042% INCLUSIVE											
TES	T 1 🖝	.100		Т	EST 2 🖝	.099	-	TEST 3 🖝	.098	3		
X	RFI DET	ECTOR OPE	RATING									
	-	HE NUMBER			S IN THE FOLLO ESTS)	WING RAN	IGES SINCE T	HE LAST M	AINTENAN	ICE REPORT:		
REF	USALS	0	(004)	0	(.0509) 0	(.10	14) 0	(.1519)	0	(OVER .19)	0	
List any new parts and describe any alteration or modification that was made to restore the instrument to operate satisfactorily and within established limits (use other side if necessary).												
INS	PECTING	OFFICER		9 3								
SIGNATURE									Brian K Kennedy			
TYPE II PERMIT NUMBER/EXPIRATION DITE 230039 03/09/2025								TELEPHONE NUMBER 816-629-7108				
					× .							
Return completed report to the: Breath Alcohol Program, MO Department of Health and Senior Services, Southeast District Office by mail, fax, or email.												

AS IV Serial no: 107977 Version no: 5328

Temp Date Time 216L TEST RECORD 01460

Calibration Check: 21 12/13/24 14:51 .108 fir Blank: 12/13/24 14:51 .000

Subject Name

Subject I.D.

Deerator Name, I.D.

AS IV Serial no: 107977 Version no: 5328

TEST RECORD 01461

Temp Date Time 210L

12/13/24 14:53 .080 Calibration Check: 22 12/13/24 14:53 .099 Air Blank:

Subject Name

Simulator

Subject I.D.

Menueny 9/1

Location ES/AD

AS IV Serial no: 107977 Version no: 532B TEST RECORD 01462

Temp Date Time 219L

Air Blank: 12/13/24 14:55 .000 Calibration Check: 23 12/13/24 14:55 .098

Subject Name

Subject I.D.

Opera or Mame, I.D. KENNEDY 911

ESMO.

AS IV Serial no: 187977 Version no: 532B

Temp Date Time 210L TEST RECORD 01463

UOID: RFI 12 12/13/24 14:56

Subject Name

Subject I.D.

Mexwell 4/

Location ESMO



CERTIFICATE OF ANALYSIS

Certified Alcohol Reference Solution for Simulator

Random Samples of Lot Number 23390 of Alcohol Reference Solution for Simulator were analyzed by gas chromatography on October 18, 2023, using a Perkin Elmer Gas Chromatograph Autosystem XL S/N: 610N9030209, and found to contain 0.1207% (w/vol) ethyl alcohol. The expiration date for this lot number is October 17, 2025 at 11:59 PM.

When used in a calibrated Simulator, operating at 34°C +/- .2°C, this solution will give a breath alcohol analysis instrument reading of 0.100 g/210L +/- 3%.

The alcohol and water used in this solution were free of test interfering substances.

Ted L. Pauley, President GUTH LABORATORIES, INC.

NIST Traceability:

Testing was conducted using Cerilliant Reference Standard lot number FN03072301 whose values are traceable to NIST.

All balances are calibrated annually by an outside agency using NIST traceable weights. Calibration verification is done prior to each use utilizing NIST traceable weights.



Missouri Department of Health and Senior Services

P.O. Box 570, Jefferson City, MO 65102-0570 Phone: 573-751-6400 FAX: 573-751-6010 RELAY MISSOURI for Hearing and Speech Impaired 1-800-735-2966 VOICE 1-800-735-2466



Paula Nickelson **Acting Director**

Michael L. Parson Governor

SIMULATOR CERTIFICATION REPORT

SIMULATOR INFORMATION

Simulator Serial Number: MP3564

Manufacturer: Guth

Model Number:

12V500

Agency:

EXCELSIOR SPRINGS

Agency Address: 301 S MAIN, EXCELSIOR SPRINGS, MO 64024

NIST THERMOMETER INFORMATION

Serial Number:

17KMM00690

Bias:

0.00

Uncertainty:

0.02

Date of Certification:

10/27/2023

Date of Expiration:

10/27/2024

ENVIRONMENTAL CONDITIONS

The environmental conditions during testing are within the tolerances of DHSS BAP method 3.

VERIFICATION RESULTS

Simulator Average

NIST Average

Combined Uncertainty

33.99

34.00

.03

The combined uncertainty is calculated with a k=2 value.

ADJUSTMENT RESULTS

No adjustment was needed.

Date of testing:

5/20/2024

Certification Expiration:

5/20/2025

Simulator testing technician: R. SCHILDKNECHT

Notes on Condition: none

Deviation(s) from method: none

DHSS BAP Scientist Approving:

BRIANNA MEDRANO

Certification No:

MP3564 5202024

DHSS BAP Scientist Approving

Simulator Calibration Certification

Issued by Lab Manager, DHSS BAP Revision Date: 06/25/2022

Breath Alcohol Program 1903 Northwood Drive, Suite 4 Poplar Bluff, MO 63901

DHSS BAP Document 3.6A Revision 2 Page 1 of **1**



STATE OF MISSOURI

DEPARTMENT OF HEALTH AND SENIOR SERVICES

BREATH ALCOHOL PROGRAM



PERMIT TYPE II

BRIAN K. KENNEDY

is hereby authorized to instruct and supervise operators, train instructors, inspect, calibrate, perform field service and repairs, and operate the following breath analyzer(s):

ALCO-SENSOR IV WITH PRINTER

for the determination of the alcoholic content of blood from a sample of expired air. Permit issued under the provisions of sections

377.020	tillough 377.041, Itoliko and 333.111 tillough 3301110 itolik	Mike Masson
DATE	3/9/2023	DIRECTOR OF STATE PUBLIC HEALTH LABORATORY
NUMBER	230039	Davla I. nielselson
EXPIRES	3/9/2025	DIRECTOR OF DEPARTMENT OF HEALTH AND SENIOR SERVICES

MO 580-0771 (6-10)

LAB-4 (R6-10)



The named cardholder is authorized to operate an evidential breath alcohol instrument for the determination of the alcoholic content in breath form of expired ai in Missouri.

Operator KENNEDY, BRIAN Permit No 230039

577 020 through 577 041 RSMo and 306 111 through 306 119 RSMo

Date Issued 3/9/2023 Date Expires 3/9/2025

