

MISSOURI DEPARTMENT OF HEALTH AND SENIOR SERVICES STATE PUBLIC HEALTH LABORATORY

ALCO-SENSOR IV WITH PRINTER MAINTENANCE REPORT

REPORT #7

7 mm2		40 - 28 66-75 51-515				
Complete this report in du Send copy to Department	plicate at the time of Health and Sen	of the regular mo ior Services; retai	nthly preventative mainten n original in department file	ance check, and v	whenever instrument is repaired.	
ALCO SENSOR IV SN 108271		Missouri	NCY State Highway Patrol		ATE OF INSPECTION 11/29/2024	
LOCATION OF INSTRUMENT (STREET AND CITY) 3131 E Kearney, Springfield, Missouri 65803				1	ME OF INSPECTION 4:26	
CHECKLIST: Place a mar	k in the box by each	n item if found to b	e satisfactory or if operating	within established	limits. (Write in observed values	
where determined.) Unmarked items must be corrected before using instrument.						
✓ DIGITAL READOUT	(ALL ELEMENTS C	PERATIONAL)				
▼ TEMPERATURE OF ALCO SENSOR (10°C - 40°C)						
PRINTER WORKING PROPERLY						
TIME AND DATE DIS	SPLAYING PROPE	RLY				
BREATH ALCOHOL ACC	CURACY STANDA	RDS				
✓ SIMULATOR SOLUTION COMPRESSED ETHANOL-GAS MIXTURE						
✓ STANDARD SUPPLI	ER Guth Labora	atories, INC	LOT # 24310	EXP. DATE _	08/27/2026	
SIMULATOR TEMPE	RATURE (34°C ± 0	_{0.2°C)} 33.99	SIM. SN MP2418	SIM. NI	ST EXP DATE 12/05/2024	
0.100% STANDA	ARD - MUST READ ARD - MUST READ	BETWEEN 0.09 BETWEEN 0.07	n being used. (PRINTOUT) 5% and 0.105% INCLUSIV 6% and 0.084% INCLUSIV 8% and 0.042% INCLUSIV	E E		
TEST 1 ▼ .101		TEST 2 ☞ .099		TEST 3 - .099		
✓RFI DETECTOR OPE	RATING					
			LOWING RANGES SINCE	THE LAST MAIN	TENANCE REPORT:	
(DO NOT INCLUDE SEL	r-ADMINISTERED	1 (2015)	1		1	
REFUSALS	(004)	(.0509)	(.1014)	(.1519)	(OVER .19)	
List any new parts and describe any alteration or modification that was made to restore the instrument to operate satisfactorily and within established limits (use other side if necessary).						
	nge	ed ti	me +	1 n	ninute	
INSPECTING OFFICER SIGNATURE				PRINT NAME D W Henley #727		
TYPE II PERMIT NUMBER/EXP PRATION DATE 240144		4 6120120	6/28/2026		TELEPHONE NUMBER	
				(417) 895 -		
Return completed repo		Alcohol Program, fax, or email.	MO Department of Health a	and Senior Service	es, Southeast District Office	
MO 580-1351 (5-19)		AN EQUAL OPPO	ORTUNITY/AFFIRMATIVE ACTION EMPLOYER		LAB-114	

TEST RECORD 60896
TEST RECORD 60896
Temp Date Time 2181
Air Blank:
11/29/24 14:29 .000
Calibration Check:
18 11/29/24 14:29 .101
Subject Name
TEST 1
Subject I.D.
Operator Name I.D.
Operator Name I.D.
Location
3)3) E KEARNEY
SQUINGFIELD MOS
SQUINGFIELD MOS
SQUINGFIELD MOS

resion no: 188271
Uersion no: 532B

TEST RECORD 00807
TEST RECORD 00807
Air Blank:
Air Blank:
11/29/24 14:30 .000
Calibration Check:
Calibration Check:
19 11/29/24 14:30 .009
Subject Name
TEST
Subject I.D.
Subject I.D.

Operator Hame, I.D.

WHENEY #797

LOCATION

SORTHOFTELDING

SORTHOFTELDING

OPERATOR HAME, I.D.

OPERATOR HAME, I.D.

LOCATION

LOCATION

SORTHOFTELDING

OPERATOR HAME, I.D.

OP

HS IV Serial not 198271
Uersion no: \$32B

TEST RECORD 00808

Temp Date Time 210L

Air Blank:
11/29/24 14:32 .000
Calibration Check:
20 11/29/24 14:32 .009
Subject Name

TEST 3

Subject I.D.

Operator Name, I.D.

UN HENLEY#737

Location

3)3) E KEARNEY

CILH Y MISSONS

AS IV Serial not 188271 Version not 5328

Temp Date Time 210L

UOID: RFI
12 11/29/24 14:33

Subject Name

\(\begin{array}{c} \begin{array}{c} \text{Subject I.D.} \end{array} \]

COLA POTENDA



CERTIFICATE OF ANALYSIS

Certified Alcohol Reference Solution for Simulator

Random Samples of Lot Number 24310 of Alcohol Reference Solution for Simulator were analyzed by gas chromatography on August 28, 2024, using a Perkin Elmer Gas Chromatograph Autosystem XL S/N: 610N9030209, and found to contain 0.1209% (w/vol) ethyl alcohol. The expiration date for this lot number is August 27, 2026 at 11:59 PM.

When used in a calibrated Simulator, operating at 34°C +/- .2°C, this solution will give a breath alcohol analysis instrument reading of 0.100 g/210L +/- 3%.

The alcohol and water used in this solution were free of test interfering substances.

Ted L. Pauley, President GUTH LABORATORIES, INC.

NIST Traceability:

1.

Testing was conducted using Cerilliant Reference Standard lot number FN03072301 whose values are traceable to NIST.

All balances are calibrated annually by an outside agency using NIST traceable weights. -Calibration verification is done prior to each use utilizing NIST traceable weights.



STATE OF MISSOURI

DEPARTMENT OF HEALTH AND SENIOR SERVICES
BREATH ALCOHOL PROGRAM



PERMIT TYPE II

DAVID W. HENLEY, JR.

is hereby authorized to instruct and supervise operators, train instructors, inspect, calibrate, perform field service and repairs, and operate the following breath analyzer(s):

ALCO-SENSOR IV WITH PRINTER, INTOX DMT



The named cardholder is authorized to operate an evidential breath alcohol instrument for the determination of the alcoholic content in breath form of expired a in Missouri.

Operator HENLEY, JR., DAVID

Permit No 240144

Date Issued 6/28/2024 Date Expires 6/28/2026

