

# MISSOURI DEPARTMENT OF HEALTH AND SENIOR SERVICES STATE PUBLIC HEALTH LABORATORY

### ALCO-SENSOR IV WITH PRINTER MAINTENANCE REPORT

REPORT #7

- shridde.x					
Complete this report in duplicate at the time Send copy to Department of Health and Ser	e of the regular monthly nior Services; retain origi	preventative mainten nal in department file	ance check, and where.	never instrument is repaired.	
ALCO SENSOR IV SN 108391	NAME OF AGENCY MSHP Troop F			OF INSPECTION 1/2024	
LOCATION OF INSTRUMENT (STREET AND CITY) 2920 North Shamrock, Jefferson City, Missouri				of inspection 0 pm	
CHECKLIST: Place a mark in the box by each	h item if found to be satis	factory or if operating	within established limi	its. (Write in observed values	
where determined.) Unmarked items must be corrected before using instrument.					
DIGITAL READOUT (ALL ELEMENTS OPERATIONAL)					
TEMPERATURE OF ALCO SENSOR (10°C - 40°C)					
PRINTER WORKING PROPERLY					
TIME AND DATE DISPLAYING PROPERLY					
BREATH ALCOHOL ACCURACY STANDA	RDS				
SIMULATOR SOLUTION   COMPRESSED ETHANOL-GAS MIXTURE					
STANDARD SUPPLIER Guth Laboratories LOT # 23390 EXP. DATE 10/17/2025					
SIMULATOR TEMPERATURE (34°C ±	0.2°C) <u>34.01</u> SI	M. SNMP241	0 SIM. NIST E	XP DATE <u>11/05/2025</u>	
CALIBRATION CHECK – (ONLY ONE S Run three tests using a standard solutio less. Check the box corresponding to the O.100% STANDARD - MUST READ O.080% STANDARD - MUST READ O.040% STANDARD - MUST READ	n. All three tests must be standard solution being BETWEEN 0.095% and BETWEEN 0.076% and	e within ±5% of the s g used. (PRINTOUT / d 0.105% INCLUSIV d 0.084% INCLUSIV	standard value and mu ATTACHED) E E	ist have a spread of .005 or	
TEST 1 • 0.100	TEST 2 • 0.100		TEST 3 • 0.100		
□ RFI DETECTOR OPERATING					
INDICATE THE NUMBER OF BREATH TESTS IN THE FOLLOWING RANGES SINCE THE LAST MAINTENANCE REPORT: (DO NOT INCLUDE SELF-ADMINISTERED TESTS)					
REFUSALS 0 (004) 0	(.0509) 0	(.1014) 0	(.1519) 0	(OVER .19) 0	
List any new parts and describe any alterati established limits (use other side if necessar All checks normal.	on or modification that v	vas made to restore	the instrument to ope	rate satisfactorily and within	
INSPECTING OFFICER					
SIGNATURE MAH FA			PRINT NAME Matthow D. Faston		
TYPE II PERMIT NUMBER/EXPIRATION DATE			Matthew D. Easton TELEPHONE NUMBER		
230127 / June 20, 2025			(573) 751-1000		
Return completed report to the: Breath Alcohol Program, MO Department of Health and Senior Services, Southeast District Office by mail, fax, or email.					

AS IV Serial no: 188391 Version no: 5328 TEST RECORD 00949 9/ Date Time 210L Air Blank: 12/07/24 11:56 .000 Calibration Check: 14 12/07/24 11:56 .100 Subject Name Test 1 Subject I.D. Operator Name, I.D. M.D. Easton 230127 Location

140

Troup F

Location

Troop F HQ

AS IV Serial no: 108391 Version no: 532B TEST RECORD Temp Date Time 210L Air Blank: 12/07/24 12:09 .000 Calibration Check: 18 12/07/24 12:09 .100 Subject Name lest 3 Subject I.D. Operator Name, I.D. M.D. Easton 230127 Location Troop F HQ

AS IV Serial no: 188391
Version no: 5328
TEST RECORD 88952

Temp Date Time 2101
12 12/87/24 12:17
Subject Name
RFI Check
Subject I.D.

Operator Name, I.D.

M.D. Easton
Tooo F H Q



#### CERTIFICATE OF ANALYSIS

Certified Alcohol Reference Solution for Simulator

Random Samples of Lot Number 23390 of Alcohol Reference Solution for Simulator were analyzed by gas chromatography on October 18, 2023, using a Perkin Elmer Gas Chromatograph Autosystem XL S/N: 610N9030209, and found to contain 0.1207% (w/vol) ethyl alcohol. The expiration date for this lot number is October 17, 2025 at 11:59 PM.

When used in a calibrated Simulator, operating at 34°C +/- .2°C, this solution will give a breath alcohol analysis instrument reading of 0.100 g/210L +/- 3%.

The alcohol and water used in this solution were free of test interfering substances.

Ted L. Pauley, President GUTH LABORATORIES, INC.

NIST Traceability:

Testing was conducted using Cerilliant Reference Standard lot number FN03072301 whose values are traceable to NIST.

All balances are calibrated annually by an outside agency using NIST traceable weights. Calibration verification is done prior to each use utilizing NIST traceable weights.



### STATE OF MISSOURI

DEPARTMENT OF HEALTH AND SENIOR SERVICES
BREATH ALCOHOL PROGRAM



# PERMIT TYPE II

## MATTHEW D. EASTON

is hereby authorized to instruct and supervise operators, train instructors, inspect, calibrate, perform field service and repairs, and operate the following breath analyzer(s):

## ALCO-SENSOR IV WITH PRINTER, INTOX DMT

for the determination of the alcoholic content of blood from a sample of expired air. Permit issued under the provisions of sections 577.020 through 577.041, RSMo and 306.111 through 306.119 RSMo.

DATE6/20/2023	Mile Masson
	DIRECTOR OF STATE PUBLIC HEALTH LABORATORY
NUMBER 230127	
EXPIRES 6/20/2025	Daves I. Nechelson
₩C 690.0771 (6.10)	DIRECTOR OF DEPARTMENT OF HEALTH AND SENIOR SERVICES

MO 580-0771 (6-10)

LAB-4 (R6-10)



The named cardholder is authorized to operate an evidential breath alcohol instrument for the determination of the alcoholic content in breath form of expired ai in Missouri

Operator EASTON, MATTHEW

Permit No 230127

