

# MISSOURI DEPARTMENT OF HEALTH AND SENIOR SERVICES STATE PUBLIC HEALTH LABORATORY

### ALCO-SENSOR IV WITH PRINTER MAINTENANCE REPORT

REPORT #7

Complete this report in d Send copy to Departmen							d whene	ver instrument is	repaired.	
ALCO SENSOR IV SN 111333		NAME OF AGENCY Missouri State Highway Patrol					DATE OF INSPECTION 09/13/2024			
LOCATION OF INSTRUMENT (STREET AND CITY) 2201 I-70 Drive Northwest Columbia, Missouri 65202							TIME OF INSPECTION 5:25 pm			
CHECKLIST: Place a mark in the box by each item if found to be satisfactory or if operating within established limits. (Write in observed values where determined.) Unmarked items must be corrected before using instrument.										
✓ DIGITAL READOUT (ALL ELEMENTS OPERATIONAL)										
✓ TEMPERATURE OF ALCO SENSOR (10°C - 40°C)										
✓ PRINTER WORKING PROPERLY										
✓ TIME AND DATE DISPLAYING PROPERLY										
BREATH ALCOHOL ACCURACY STANDARDS										
✓ SIMULATOR SOLUTION ☐ COMPRESSED ETHANOL-GAS MIXTURE										
							NIST EXP DATE _04/03/2025			
<ul> <li>✓ CALIBRATION CHECK – (ONLY ONE STANDARD IS TO BE USED PER MAINTENANCE REPORT)         Run three tests using a standard solution. All three tests must be within ±5% of the standard value and must have a spread of .005 or less. Check the box corresponding to the standard solution being used. (PRINTOUT ATTACHED)         ✓ 0.100% STANDARD - MUST READ BETWEEN 0.095% and 0.105% INCLUSIVE         ○ 0.080% STANDARD - MUST READ BETWEEN 0.076% and 0.084% INCLUSIVE         ○ 0.040% STANDARD - MUST READ BETWEEN 0.038% and 0.042% INCLUSIVE     </li> </ul>										
TEST 1 <b>●</b> 0.104		TEST 2   0.104				TEST 3    0.104				
☑ RFI DETECTOR OPERATING										
INDICATE THE NUMBER OF BREATH TESTS IN THE FOLLOWING RANGES SINCE THE LAST MAINTENANCE REPORT: (DO NOT INCLUDE SELF-ADMINISTERED TESTS)										
REFUSALS 0	(004) 0	(.0509)	0	(.1014)	0	(.1519)	0	(OVER .19)	0	
List any new parts and describe any alteration or modification that was made to restore the instrument to operate satisfactorily and within established limits (use other side if necessary).  Battery Replaced Time Adjusted (Forward one minute)										
INSPECTING OFFICER							PRINT NAME			
SIGNATURE					Trooper D. L. Holliday					
түр <del>е и релміт</del> пимвеп/ехріпаті 240204/ 09-06-2026	85	4 6			TELEPHONE NUMBER (573) 751-1000					
Return completed report to the: Breath Alcohol Program, MO Department of Health and Senior Services, Southeast District Office by mail, fax, or email.										

AS IV Serial no: 111333 Version no: 532B

TEST RECORD 01120

Temp Date Time 210L

Air Blank:

09/13/24 17:48 .000

Calibration Check:

23 09/13/24 17:48 .104

Subject Name

Acc 1

Subject I.D.

Ac 1

Operator Name, I.D.

David Holliday 240 204

Location

2201 I-70 Drive

Worthwest Columbia, mo 65202

> AS IV Serial no: 111333 Version no: 532B

> > TEST RECORD 01122

Temp Date Time 210L

Air Blank:

09/13/24 17:52 .000

Calibration Check:

24 09/13/24 17:52 .104

Subject Name

Acc 3

Subject I.D.

Acc 3

Operator Name, I.D.

David Holliday 240204

Location

2201 I-70 Drau

Northwest Colombia, mo

65202

AS IV Serial no: 111333 Version no: 532B

TEST RECORD 01121

Temp Date Time 210L

Air Blank:

09/13/24 17:50 .000

Calibration Check: 23 09/13/24 17:50 .104

23 09/13/24 17:50 .104

Subject Name

Acc 2

Subject I.D.

Acc 2

Operator Name, I.D.

David Helliday 240204

Location

2201 I-70 Drive

Northwest Colombia no

65202

AS IV Serial no: 111333 Version no: 532B

TEST RECORD 01123

Temp Date Time 210L

UOID: RFI 12 09/13/24 17:56

Subject Name

RFI Test

Subject I.D.

Operator Name, I.D.

David Holl deg 240204

Location

2201 I.70 Dage

Watherst Columbia no



#### CERTIFICATE OF ANALYSIS

Certified Alcohol Reference Solution for Simulator

Random Samples of Lot Number 23180 of Alcohol Reference Solution for Simulator were analyzed by gas chromatography on May 18, 2023, using a Perkin Elmer Gas Chromatograph Autosystem XL S/N: 610N9030209, and found to contain 0.1220% (w/vol) ethyl alcohol. The expiration date for this lot number is May 17, 2025 at 11:59 PM.

When used in a calibrated Simulator, operating at 34°C +/- .2°C, this solution will give a breath alcohol analysis instrument reading of 0.100 g/210L +/- 3%.

The alcohol and water used in this solution were free of test interfering substances.

Ted L. Pauley, President GUTH LABORATORIES, INC.

NIST Traceability:

Testing was conducted using Cerilliant Reference Standard lot number FN11172002 whose values are traceable to NIST.

All balances are calibrated annually by an outside agency using NIST traceable weights. Calibration verification is done prior to each use utilizing NIST traceable weights.



### STATE OF MISSOURI

DEPARTMENT OF HEALTH AND SENIOR SERVICES
BREATH ALCOHOL PROGRAM



# PERMIT TYPE II

## DAVID L. HOLLIDAY

is hereby authorized to instruct and supervise operators, train instructors, inspect, calibrate, perform field service and repairs, and operate the following breath analyzer(s):

## ALCO-SENSOR IV WITH PRINTER, INTOX DMT

for the determination of the alcoholic content of blood from a sample of expired air. Permit issued under the provisions of sections

DATE 9/6/2024

NUMBER 240204

EXPIRES 9/6/2026

DIRECTOR OF DEPARTMENT OF HEALTH AND SENIOR SERVICES

MO 580-0771 (6-10)

LAB-4 (R6-10)



The named cardholder is authorized to operate an evidential breath alcohol instrument for the determination of the alcoholic content in breath form of expired all in Missouri.

Operator HOLLIDAY, DAVID

Permit No 240204

Date Issued 9/6/2024 Date Expires 9/6/2026

