



MISSOURI DEPARTMENT OF HEALTH AND SENIOR SERVICES STATE PUBLIC HEALTH LABORATORY

ALCO-SENSOR IV WITH PRINTER MAINTENANCE REPORT

REPORT #7

- 1858 A.				
Complete this report in duplicate at the time of the regular monthly preventative maintenance check, and whenever instrument is repaired. Send copy to Department of Health and Senior Services; retain original in department file.				
ALCO SENSOR IV SN	NAME OF AGENCY		DATE OF	NSPECTION
111642	[Campbell Po	ce Departme	n+ 112-	3-2024
LOCATION OF INSTRUMENT (STREET AND CITY)			TIME OF I	NSPECTION
204 Grand Ave Car	mpbell, MO 639.	33		
CHECKLIST: Place a mark in the box by each item if found to be satisfactory or if operating within established limits. (Write in observed values				
where determined.) Unmarked items must be corrected before using instrument.				
DIGITAL READOUT (ALL ELEMENTS OPERATIONAL)				
TEMPERATURE OF ALCO SENSOR (10°C - 40°C)				
PRINTER WORKING PROPERLY				
TIME AND DATE DISPLAYING PROPERLY				
BREATH ALCOHOL ACCURACY STANDARDS				
SIMULATOR SOLUTION		COMPRESSED ETH	ANOL-GAS MIXT	URE
STANDARD SUPPLIER Uth		23390 E	(P. DATE <u>10-1</u>	1-2025
SIMULATOR TEMPERATURE (34°C ± 0.2°	c) 34°C sim. sn	SD 1743	SIM. NIST EXF	DATE 12-19-2025
Run three tests using a standard solution. All three tests must be within ±5% of the standard value and must have a spread of .005 or less. Check the box corresponding to the standard solution being used. (PRINTOUT ATTACHED) 0.100% STANDARD - MUST READ BETWEEN 0.095% and 0.105% INCLUSIVE 0.080% STANDARD - MUST READ BETWEEN 0.076% and 0.084% INCLUSIVE 0.040% STANDARD - MUST READ BETWEEN 0.038% and 0.042% INCLUSIVE				
TEST 1 - CO98 TE	ST 2 * .098	TEST	3 ** • 0 9 5	2
TRFI DETECTOR OPERATING				
INDICATE THE NUMBER OF BREATH TESTS IN THE FOLLOWING RANGES SINCE THE LAST MAINTENANCE REPORT: (DO NOT INCLUDE SELF-ADMINISTERED TESTS)				
	.100509)		19)	(OVER .19)
List any new parts and describe any alteration or modification that was made to restore the instrument to operate satisfactorily and within established limits (use other side if necessary).				
INSPECTING OFFICER				
SIGNATURE		PRINTN		
· Serry D. Hararans			cry Don I	largraves
TYPE PERINT NUMBER/EXPIRATION DATE		1 · · · · · · · · · · · · · · · · · · ·	ONE NUMBER)/ j j
230319/ 12-21-2025	7	(573	51 217 - 9	<u>/</u> \d
Return completed report to the: Breath Alcohol Program, MO Department of Health and Senior Services, Southeast District Office by mail, fax, or email.				

AS IV Serial no: 111642 Version no: 532B

TEST RECORD 00607

9/ Time 210L

Air Blank: 12/03/24 13:55 .000

Date

9ubject Test: Man 18 12/03/24 13:55 .098

Subject Name

test#1

Temp

Subject I.D.

Jessy Hagsaves 230319 Operator Name, I.D.

Campbell Police Dept Location

Booking room

AS IV Serial no: 111642 Version no: 532B

TEST RECORD 00608

Temp Date Time 210L

Air Blank: 12/03/24 13:56 .000 Subject Test: Man 19 12/03/24 13:56 .098

Subject Name

+es+ #2

Subject I.D.

Jerry Hagraves 230319 Operator Name, I.D.

Campbell Police Dept Location

Booking room

AS IV Serial no: 111642 Version no: 532B

TEST RECORD 00609

Temp Date Time 218L

Air Blank: 12/03/24 13:58 .000 Subject Test: Man 19 12/03/24 13:58 .098

Subject Name

tes+ # 3

Subject I.D.

<u>Jerry Horgraves 230319</u> Operator Mame, I.D.

Campbell Police Dept Location

Booking room

AS IV Serial no: 111642 Version no: 532B

TEST RECORD 00610

Temp Date Time 210L VOID: RTI

12 12/33/24 14:00

Subject Name

test#4 RFI test

Jerry Hargra ves 230319 Operator Name, I.D.

Campbell Police Dept location

Booking room



Missouri Department of Health and Senior Services

P.O. Box 570, Jefferson City, MO 65102-0570 Phone: 573-751-6400 FAX: 573-751-6010 RELAY MISSOURI for Hearing and Speech Impaired 1-800-735-2966 VOICE 1-800-735-2466



Paula Nickelson **Acting Director**

Michael L. Parson Governor

SIMULATOR CERTIFICATION REPORT

SIMULATOR INFORMATION

Simulator Serial Number: SD1743

Manufacturer: Guth

Model Number:

10-4D

Agency:

CAMPBELL PD

Agency Address: 204 W GRAND, CAMPBELL, MO 63933

NIST THERMOMETER INFORMATION

Serial Number:

17KMM00690

Bias:

0.00

Uncertainty:

0.02

Date of Certification:

10/27/2023

Date of Expiration:

10/27/2024

ENVIRONMENTAL CONDITIONS

The environmental conditions during testing are within the tolerances of DHSS BAP method 3.

VERIFICATION RESULTS

Simulator Average

NIST Average

Combined Uncertainty .06

34.00

34.03 The combined uncertainty is calculated with a k=2 value.

ADJUSTMENT RESULTS

No adjustment was needed.

Date of testing:

12/19/2023

Certification Expiration:

12/19/2024

Simulator testing technician: R. SCHILDKNECHT

Notes on Condition: none

Deviation(s) from method: none

DHSS BAP Scientist Approving:

BRIANNA MEDRANO

Certification No:

SD1743_12192023

DHSS BAP Scientist Approving

Simulator Calibration Certification

Issued by Lab Manager, DHSS BAP Revision Date: 06/25/2022

Breath Alcohol Program 1903 Northwood Drive, Suite 4 Poplar Bluff, MO 63901

DHSS BAP Document 3.6A Revision 2 Page 1 of 1



CERTIFICATE OF ANALYSIS

Certified Alcohol Reference Solution for Simulator

Random Samples of Lot Number 23390 of Alcohol Reference Solution for Simulator were analyzed by gas chromatography on October 18, 2023, using a Perkin Elmer Gas Chromatograph Autosystem XL S/N: 610N9030209, and found to contain 0.1207% (w/vol) ethyl alcohol. The expiration date for this lot number is October 17, 2025 at 11:59 PM.

When used in a calibrated Simulator, operating at 34°C +/- .2°C, this solution will give a breath alcohol analysis instrument reading of 0.100 g/210L +/- 3%.

The alcohol and water used in this solution were free of test interfering substances.

Ted L. Pauley, President GUTH LABORATORIES, INC.

NIST Traceability:

Testing was conducted using Cerilliant Reference Standard lot number FN03072301 whose values are traceable to NIST.

All balances are calibrated annually by an outside agency using NIST traceable weights. Calibration verification is done prior to each use utilizing NIST traceable weights.



STATE OF MISSOURI

DEPARTMENT OF HEALTH AND SENIOR SERVICES **BREATH ALCOHOL PROGRAM**



PERMIT TYPE II JERRY HARGRAVES

is hereby authorized to instruct and supervise operators, train instructors, inspect, calibrate, perform field service and repairs, and operate the following breath analyzer(s):

ALCO-SENSOR IV WITH PRINTER for the determination of the alcoholic content of blood from a sample of expired air. Permit issued under the provisions of sections 577.020 through 577.041, RSMo and 306.111 through 306.119 RSMo. Mile Masson DIRECTOR OF STATE PUBLIC HEALTH LABORATORY DATE 12/21/2023 Daven J. nichelson NUMBER 230319 DIRECTOR OF DEPARTMENT OF HEALTH AND SENIOR SERVICES EXPIRES 12/21/2025 1 AB-4 (PIS-10)

MO 580-0771 (6-10)



STATE OF MISSOURI DEPARTMENT OF HEALTH AND SENIOR SERVICES BREATH ALCOHOL PROGRAM

INSTRUMENT OPERATOR CARD

The named cardholder is authorized to operate an evidential breath alcohol instrument for the determination of the alcoholic content in breath form of expired air

HARGRAVES, JERRY Operator

230319 Permit No

Date Expires 12/21/2025 Date Issued 12/21/2023

