

MISSOURI DEPARTMENT OF HEALTH AND SENIOR SERVICES STATE PUBLIC HEALTH LABORATORY

ALCO-SENSOR IV WITH PRINTER MAINTENANCE REPORT

REPORT #7

THEET							
Complete this report in duplicate at the time Send copy to Department of Health and Sen				whenev	ver instrument is	repaired.	
ALCO SENSOR IV SN 111743	NAME OF AGENCY Missouri State Highway Patrol			DATE OF INSPECTION 08/02/2024			
LOCATION OF INSTRUMENT (STREET AND CITY) 6012 NE Antioch Road, Gladstone, MO 64119				TIME OF INSPECTION 10:50 pm			
CHECKLIST: Place a mark in the box by each			within establishe	ed limits.	. (Write in observe	ed values	
where determined.) Unmarked items must be	e corrected before using	instrument.					
DIGITAL READOUT (ALL ELEMENTS C	PERATIONAL)			n-1472-2017			
✓ TEMPERATURE OF ALCO SENSOR (1)	0°C - 40°C)						
☑ PRINTER WORKING PROPERLY							
☑ TIME AND DATE DISPLAYING PROPE	RLY						
BREATH ALCOHOL ACCURACY STANDA	RDS						
☑ SIMULATOR SOLUTION ☐ COMPRESSED ETHANOL-GAS MIXTURE							
✓ STANDARD SUPPLIER Guth Laboratories LOT # 23180 EXP. DATE 05/17/2025							
SIMULATOR TEMPERATURE (34°C ± 0	0.2°C)34.00 SII	M. SNMP242	4 SIM. N	IST EXF	P DATE 11/27/2	2024	
CALIBRATION CHECK – (ONLY ONE S' Run three tests using a standard solution less. Check the box corresponding to the 0.100% STANDARD - MUST READ 0.080% STANDARD - MUST READ 0.040% STANDARD - MUST READ	n. All three tests must be standard solution being BETWEEN 0.095% and BETWEEN 0.076% and	e within ±5% of the s gused. (PRINTOUT / d 0.105% INCLUSIV d 0.084% INCLUSIV	standard value ar ATTACHED) E E		have a spread of	of .005 or	
TEST 1 ▼ .105	TEST 2 • .105	TEST 3 .105					
☑ RFI DETECTOR OPERATING	Maria de la composición del composición de la co						
INDICATE THE NUMBER OF BREATH TES (DO NOT INCLUDE SELF-ADMINISTERED		G RANGES SINCE	THE LAST MAIN	ITENAN	ICE REPORT:		
REFUSALS 0 (004) 0	(.0509) 0	(.1014) 0	(.1519)	0	(OVER .19)	0	
List any new parts and describe any alterative stablished limits (use other side if necessary		was made to restore	the instrument to	o operat	e satisfactorily a	nd within	
INSPECTING OFFICER							
SIGNATURE			PRINT NAME CPL. B. E. THOMAN #1256				
TYPE II PERMIT NUMBER/EXPIRATION DATE			TELEPHONE NUMBER				
240111 05/21/2026			(816) 622-0800				
	lcohol Program, MO Del fax, or email.	partment of Health a	nd Senior Service	es, Sout	heast District Off	fice	

98 IV Serial no: 111743 Jersion no: 532B

TEST RECORD 99849

Tamp Date Time 2101
Air Blank: 88/02/24 22:54 .000
Calibration Chark:

Subject Name

MAINT TOT 41 Subject I.D.

Operator Names I.D.

BETHORAN 24011

CLAY CO. Z.O.

AS IV Serial r : 111743 Version no: 582E

TEST RECOPD 80850

Temp Date Time 2101 Air Blank:

08/02/24 0 157 gen Calibration #434 23 08/02/24 / 157 mma

Subject Name

MAINT TEST #L Subject I.D.

Operator Name, I.D.

BE THOTIAN 240111 Location CLAY (O. Z.U. AS IV Serial no: 111743

TEST RECORD 00851

Temp Date Time 2101 Air Blank:

88/82/24 33:00 .900 Calibration flool 24 88/02/20 7438 .187

Subject Name

MINTER #3

Operator Mames I.D.

B.E. THOMAN 24011

CLAY (0. E.O -

AS IU Serial no: 111743 version no: 532B

TEST RECORD 60852

Temp Date Time 210L

VOID: RFI -12 08/02/24 23:01

Subject Name

RFITEST Subject I.D.

Operator Name: I.D.

B.E. THOMAN ZYON)

CLAY (0.2.0 -



CERTIFICATE OF ANALYSIS

Certified Alcoho Reference Solution for Simulator

Random Samples of Lot Number 23180 of
Alcohol Reference Solution for Simulator were analyzed by
gas chromatography on May 18, 2023, using a Perkin Elmer Gas Chromatograph
Autosystem XL S/N: 610N9030209, and found to contain 0.1220% (w/vol)
ethyl alcohol. The expiration date for this lot
number is May 17, 2025 at 11:59 PM.

When used in a calibrated Simulator, operating at 34°C +/- .2°C, this solution will give a breath alcohol analysis instrument reading of 0.100 g/210L +/- 3%.

The alcohol and water used in this solution were free of test interfering substances.

Ted L. Pauley, President GUTH LABORATORIES, INC.

NIST Traceability:

Testing was conducted using Cerviliant Reference Standard tot number FN11172002 whose values are traceable to NIST.

All balances are calibrated annually by an outside agency using NIST traceable weights. Calibration verification is done pripr to each use utilizing NIST traceable weights.



STATE OF MISSOURI

DEPARTMENT OF HEALTH AND SENIOR SERVICES **BREATH ALCOHOL PROGRAM**



PERMIT TYPE II

BRUCE E. THOMAN

is hereby authorized to instruct and supervise operators, train instructors, inspect, calibrate, perform field service and repairs, and operate the following breath analyzer(s):

ALCO-SENSOR IV WITH PRINTER, INTOX DMT

for the determination of the alcoholic content of blood from a sample of expired air. Permit issued under the provisions of sections 577.020 through 577.041, RSMo and 306.111 through 306.119 RSMo. MA 1 INA

Mike Massin				
DIRECTOR OF STATE PUBLIC HEALTH LABORATORY				
Davea J. Nichelson				
DIRECTOR OF DEPARTMENT OF HEALTH AND SENIOR SERVICES				

MO 580-0771 (6-10)

