

## MISSOURI DEPARTMENT OF HEALTH AND SENIOR SERVICES STATE PUBLIC HEALTH LABORATORY

### ALCO-SENSOR IV WITH PRINTER MAINTENANCE REPORT

REPORT #7

Complete this report in duplicate at the time Send copy to Department of Health and Sent				rhenever instrument is repaired.	
ALCO SENSOR IV SN 111774				TE OF INSPECTION 5/04/2024	
LOCATION OF INSTRUMENT (STREET AND CITY) 106 S. Folger St, Carrollton			TIN	ME OF INSPECTION :35 am	
CHECKLIST: Place a mark in the box by each item if found to be satisfactory or if operating within established limits. (Write in observed values where determined.) Unmarked items must be corrected before using instrument.					
DIGITAL READOUT (ALL ELEMENTS OPERATIONAL)					
▼ TEMPERATURE OF ALCO SENSOR (10°C - 40°C)					
☑ TIME AND DATE DISPLAYING PROPERLY					
BREATH ALCOHOL ACCURACY STANDARDS					
SIMULATOR SOLUTION   COMPRESSED ETHANOL-GAS MIXTURE					
STANDARD SUPPLIER Guth Labs	LOT # 23390 EXP. DATE 10/17/2025				
Z SIMULATOR TEMPERATURE (34°C $\pm$ 0.2°C) 34.02 SIM. SN MP3879 SIM. NIST EXP DATE $04/22/2025$					
less. Check the box corresponding to the standard solution being used. (PRINTOUT ATTACHED)  ✓ 0.100% STANDARD - MUST READ BETWEEN 0.095% and 0.105% INCLUSIVE  ☐ 0.080% STANDARD - MUST READ BETWEEN 0.076% and 0.084% INCLUSIVE  ☐ 0.040% STANDARD - MUST READ BETWEEN 0.038% and 0.042% INCLUSIVE					
TEST 1 🕶 .103	TEST 2 🖛 .103		TEST 3 🖛 .103		
☑ RFI DETECTOR OPERATING					
INDICATE THE NUMBER OF BREATH TESTS IN THE FOLLOWING RANGES SINCE THE LAST MAINTENANCE REPORT: (DO NOT INCLUDE SELF-ADMINISTERED TESTS)					
REFUSALS 0 (004) 0	(.0509) 0	(.1014) 0	(.1519) 1	(OVER .19) 0	
List any new parts and describe any alteration established limits (use other side if necessary		vas made to restore	the instrument to o	operate satisfactorily and within	
INSPECTING OFFICER			DDINT MAME		
SIGNATURE Brian C Words		-	PRINT NAME Brian C. Woods	3	
TYPE II PERMIT NUMBER/EXPIRATION DATE 220181 07/12/2024			TELEPHONE NUMBER (660) 542-2200		
Return completed report to the: Breath A by mail, 1	lcohol Program, MO Der ax, or email.	partment of Health a	nd Senior Services	, Southeast District Office	

AS:TU Serial not 111774 AS IV Serial not 111774 Version not 5320 AS IV Serial no: 111774 Warsion no: 1928 Version not 5028 yysy pyroed 90511 TEST RECORD 00510 TEST RICORD BUSIS Temp Rate Time 210L Temp Date Time 210L Temp Date Time 210). and the second second second Air Black: Air Blank: Air Mianki 65744724 65137 .006 25/84/24 85:35 .060 880 (817) (817) (817) (818) (819) Calibration Check 21 85/04/20 85:35 .183 Calibration Checks Calling Align Clercki 22 95/241/24 95:32 .183 10 6 786, 21 65:29 .163 Sulport Rage Two sales - Pasmer Subject Mame Test#2 Tes+#1 Test#3 defend 1.0. Subject 1.1. Subject I.D. Oronator Names I.D. Operator Names J.D. Operator Name, I.D. B.C. Woods, 782 B.C. Woods, 782 B. C. Woods, 782

> AS IV Serial not 111774 Version not 1828

> > TEXT RECORD 80513

Tem: Cale Time 2101

VOID: REI 12 85/84/24 05:40

Subject Hame

RFI Test Subject I.D.

Operator Name: I.D.

B.C. Woods, 782

AS IV Serial net 111774 Version set 1928

TEST RACCOU 88514

s/ Tone Hate Time 218L

Air Mlanko 85/MAZZ4 08:32 .000 Subject fost: Muto 24 05/00,24 08:42 .000

lawladreck Harrie

Blank Test SubJool I.D.

Copyrator Mames 1.D.

B. C. Woods, 782



### CERTIFICATE OF ANALYSIS

Certified Alcohol Reference Solution for Simulator

Random Samples of Lot Number 23390 of Alcohol Reference Solution for Simulator were analyzed by gas chromatography on October 18, 2023, using a Perkin Elmer Gas Chromatograph Autosystem XL S/N: 610N9030209, and found to contain 0.1207% (w/vol) ethyl alcohol. The expiration date for this lot number is October 17, 2025 at 11:59 PM.

When used in a calibrated Simulator, operating at 34°C +/- .2°C, this solution will give a breath alcohol analysis instrument reading of 0.100 g/210L +/- 3%.

The alcohol and water used in this solution were free of test interfering substances.

Ted L. Pauley, President GUTH LABORATORIES, INC.

NIST Traceability:

Testing was conducted using Cerilliant Reference Standard lot number FN03072301 whose values are traceable to NIST.

All balances are calibrated annually by an outside agency using NIST traceable weights. Calibration verification is done prior to each use utilizing NIST traceable weights.



### STATE OF MISSOURI

DEPARTMENT OF HEALTH AND SENIOR SERVICES
BREATH ALCOHOL PROGRAM



# PERMIT TYPE II

## **BRIAN C. WOODS**

is hereby authorized to instruct and supervise operators, train instructors, inspect, calibrate, perform field service and repairs, and operate the following breath analyzer(s):

### ALCO-SENSOR IV WITH PRINTER

for the determination of the alcoholic content of blood from a sample of expired air. Permit issued under the provisions of sections 577.020 through 577.041, RSMo and 306.111 through 306.119 RSMo.

DATE 7/12/2022

DIRECTOR OF STATE PUBLIC HEALTH LABORATORY

NUMBER **220181**EXPIRES **7/12/2024** 

O 580-0771 (6-10)

Davis J. Nichelson

DIRECTOR OF DEPARTMENT OF HEALTH AND SENIOR SERVICES

LAB-4 (R6-10)



#### STATE OF MISSOURI

DEPARTMENT OF HEALTH AND SENIOR SERVICES BREATH ALCOHOL PROGRAM

#### **INSTRUMENT OPERATOR CARD**

The named cardholder is authorized to operate an evidential breath alcohol instrument for the determination of the alcoholic content in breath form of expired a in Missourt.

Operator WOODS, BRIAN

Permit No 220181 Date Issued 7/12/2022 Da

22 Date Expires 7/12/2024

