

MISSOURI DEPARTMENT OF HEALTH AND SENIOR SERVICES STATE PUBLIC HEALTH LABORATORY

ALCO-SENSOR IV WITH PRINTER MAINTENANCE REPORT

REPORT #7

- This is							
Complete this report in duplicate at the tin Send copy to Department of Health and S				ck, and whenev	ver instrument is	repaired.	
ALCO SENSOR IV SN 119933					DATE OF INSPECTION 06/14/2024		
LOCATION OF INSTRUMENT (STREET AND CITY) 891 TECHNOLOGY DRIVE, WELDON SPRING (TROOP C HQ) (MOBILE INS				TIME OF I	TIME OF INSPECTION 2:14 pm		
CHECKLIST: Place a mark in the box by ea				•	3333	ed values	
where determined.) Unmarked items must		250			<u> </u>		
✓ DIGITAL READOUT (ALL ELEMENTS OPERATIONAL)							
✓ TEMPERATURE OF ALCO SENSOR (10°C - 40°C)							
☑ TIME AND DATE DISPLAYING PROPERLY							
BREATH ALCOHOL ACCURACY STAND	ARDS						
✓ SIMULATOR SOLUTION ☐ COMPRESSED ETHANOL-GAS MIXTURE							
✓ STANDARD SUPPLIER GUTH LABORATORIES, INC. LOT # 23390 EXP. DATE 10/17/2025							
✓ SIMULATOR TEMPERATURE (34°C ± 0.2°C) 34.00 SIM. SN MP2459 SIM. NIST EXP DATE 01/18/2025							
Run three tests using a standard solution. All three tests must be within ±5% of the standard value and must have a spread of .005 or less. Check the box corresponding to the standard solution being used. (PRINTOUT ATTACHED) 0.100% STANDARD - MUST READ BETWEEN 0.095% and 0.105% INCLUSIVE 0.080% STANDARD - MUST READ BETWEEN 0.076% and 0.084% INCLUSIVE 0.040% STANDARD - MUST READ BETWEEN 0.038% and 0.042% INCLUSIVE							
TEST 1 .096	TEST 2 ▼ .096		TEST 3	.096			
☑ RFI DETECTOR OPERATING							
INDICATE THE NUMBER OF BREATH TESTS IN THE FOLLOWING RANGES SINCE THE LAST MAINTENANCE REPORT: (DO NOT INCLUDE SELF-ADMINISTERED TESTS)							
REFUSALS 0 (004) 0	(.0509)	(.1014)	0 (.151	9) 0	(OVER .19)	0	
List any new parts and describe any alter established limits (use other side if necess		at was made to r	estore the instru	ument to operat	te satisfactorily a	and within	
INSPECTING OFFICER							
SIGNATURE 1 11 1/2				PRINT NAME A. MICHELS			
TYPE II PERMIT NUMBER/EXPIRATION DATE			The second second	TELEPHONE NUMBER			
PERMIT #230300, EXPIRES 12/11/2025			(636) 3	(636) 300-2800			
Return completed report to the: Breath by ma	n Alcohol Program, MO lil, fax, or email.	Department of He	ealth and Senior	Services, Sout	theast District Of	ffice	

AS IV Serial no: 119933
Version no: 532C

TEST RECORD - REPRINT

TEST RECORD 00041

Temp Date Time 210L

Air Blank: 06/14/24 14:13 .000
Calibration Check: 21 08/14/24 14:13 .096

Subject Name

Main +

Subject I.D.

Operator Name, I.D.

A-Michels #2503ccc

Location

Troop CHO

AS IV Serial no: 119933
Version no: 532C

TEST RECORD 00042

Jemp Date Time 210L

Air Blank:
06/14/24 14:17 .000
Calibration Check:
22 06/14/24 14:17 .096

Subject Name

Main +
Subject I.D.

Operator Name, I.D.

A. Michels 4230300
Location

Toop C HG

AS IV Serial no: 119933
Version no: 532C

TEST RECORD 00043

Temp Date Time 210L

Air Blank: 06/14/24 14:19 .000
Calibration Check: 23 06/14/24 14:19 .096

Subject Name

Main †
Subject I.D.

Operator Name, I.D.

A.Michels #250300
Location

Troop CHQ

AS IV Serial no: 119933
Version no: 532C

TEST RECORD 00044

Temp Date Time 210L

VOID: RFI
12 06/14/24 14:20

Subject Name

Main 7

Subject I.D.

Operator Name. I.D.

A. Michels #230500

Location

From C H Q



CERTIFICATE OF ANALYSIS

Certified Alcohol Reference Solution for Simulator

Random Samples of Lot Number 23390 Alcohol Reference Solution for Simulator were analyzed by chromatography on October 18, 2023, using a Perkin Elmer Gas gas Chromatograph Autosystem XL S/N: 610N9030209, and found to contain 0.1207% (w/vol) ethyl alcohol. The expiration date for this lot number is October 17, 2025 at 11:59 PM.

When used in a calibrated Simulator, operating at 34°C +/- .2°C, this solution will give a breath alcohol analysis instrument reading of 0.100 g/210L +/- 3%.

The alcohol and water used in this solution were free of test interfering substances.

> Ted L. Pauley, President GUTH LABORATORIES, INC.

NIST Traceability:

Testing was conducted using Cerilliant Reference Standard lot number FN03072301 whose values are traceable to NIST.

All balances are calibrated annually by an outside agency using NIST traceable weights. Calibration verification is done prior to each use utilizing NIST traceable weights.