



MISSOURI DEPARTMENT OF HEALTH AND SENIOR SERVICES STATE PUBLIC HEALTH LABORATORY

ALCO-SENSOR IV WITH PRINTER MAINTENANCE REPORT

REPORT #7

Weekler .						
Complete this report in duplicate at the time Send copy to Department of Health and Sen	수 없는 경소에게 있다면 하나의 어린다면 가득하면 아이들이 그 물건이 살아 있다면 하나 이 사람이 없다면 하나 하다.			wheneve	er instrument is	repaired.
ALCO SENSOR IV SN 120302	NAME OF AGENCY National Park Service (Arch Grounds)			DATE OF INSPECTION 06/14/2024		
LOCATION OF INSTRUMENT (STREET AND CITY) 891 TECHNOLOGY DRIVE, WELDON SPRING (TROOP C HQ) (MOBILE INS				TIME OF INSPECTION 1:53 pm		
CHECKLIST: Place a mark in the box by eac	n item if found to be satisfa	ctory or if operating	within establishe	d limits. (Write in observe	ed values
where determined.) Unmarked items must b	e corrected before using in	nstrument.				
☑ DIGITAL READOUT (ALL ELEMENTS (OPERATIONAL)					
✓ TEMPERATURE OF ALCO SENSOR (10°C - 40°C)						
✓ PRINTER WORKING PROPERLY						
☑ TIME AND DATE DISPLAYING PROPE						
BREATH ALCOHOL ACCURACY STANDA	RDS					
SIMULATOR SOLUTION COMPRESSED ETHANOL-GAS MIXTURE						
✓ STANDARD SUPPLIER GUTH LABORATORIES, INC. LOT # 23390 EXP. DATE 10/17/2025						
✓ SIMULATOR TEMPERATURE (34°C ±	SIMULATOR TEMPERATURE (34°C ± 0.2°C) 34.00 SIM. SN MP2459 SIM. NIST EXP DATE 01/18/2025					
less. Check the box corresponding to the 0.100% STANDARD - MUST REAL 0.080% STANDARD - MUST REAL 0.040% STANDARD - MUST REAL	BETWEEN 0.095% and BETWEEN 0.076% and	0.105% INCLUSIVE 0.084% INCLUSIVE				
TEST 1 .101	▼.101 TEST 2 ▼ .100		TEST 3 • .099			
☑ RFI DETECTOR OPERATING						
INDICATE THE NUMBER OF BREATH TES (DO NOT INCLUDE SELF-ADMINISTERED		RANGES SINCE 1	THE LAST MAIN	ITENANO	CE REPORT:	
REFUSALS 0 (004) 0	(.0509) 0	(.1014) 0	(.1519)	0	(OVER .19)	0
List any new parts and describe any alterat established limits (use other side if necessar		as made to restore	the instrument to	o operate	satisfactorily a	nd within
INSPECTING OFFICER	e de la composition della comp					
SIGNATURE			PRINT NAME A. MICHELS			
TYPE II PERMIT NUMBER/EXPIRATION DATE			TELEPHONE NUMBER			
PERMIT #230300, EXPIRES 12/11/2025			(636) 300-2800			
Return completed report to the: Breath A by mail,	alcohol Program, MO Depa fax, or email.	artment of Health ar	nd Senior Service	es, South	east District Of	fice

AS IV Serial no: 120302
Version no: 511A

TEST RECORD 00021

Temp Date Time %BAC

Air Blank:
06/14/24 13:53 .000
Calibration Check:
22 06/14/24 13:53 .101

Subject Name
Maint

Subject I.D.

Operator Name, I.D.

A. Michels #250300

Location

7500p CHG

Troop CHQ

AS IV Serial no: 120302
Version no: 511A

TEST RECORD 00023

Temp Date Time %BAC

Air Blank:
06/14/24 13:55 .000
Calibration Check:
23 06/14/24 13:55 .099

Subject Name
Main +
Subject I.D.

Operator Name, I.D.

A. Michel #230300
Location

Troop C HCe

AS IV Serial no: 120302
Version no: 5114

TEST RECORD 00024

Temp Date Time %BAC
Void: RFI
12 06/14/24 13:56

Subject Name

Alloy P
Subject I.D.

Operator Name, I.D.

Alloyels HJ323200
Location

Texy CH&



CERTIFICATE OF ANALYSIS

Certified Alcohol Reference Solution for Simulator

Random Samples of Lot Number 23390 Alcohol Reference Solution for Simulator were analyzed by chromatography on October 18, 2023, using a Perkin Elmer Gas gas Chromatograph Autosystem XL S/N: 610N9030209, and found to contain 0.1207% (w/vol) ethyl alcohol. The expiration date for this lot number is October 17, 2025 at 11:59 PM.

When used in a calibrated Simulator, operating at 34°C +/- .2°C, this solution will give a breath alcohol analysis instrument reading of 0.100 g/210L +/- 3%.

The alcohol and water used in this solution were free of test interfering substances.

> Ted L. Pauley, President GUTH LABORATORIES, INC.

NIST Traceability:

Testing was conducted using Cerilliant Reference Standard lot number FN03072301 whose values are traceable to NIST.

All balances are calibrated annually by an outside agency using NIST traceable weights. Calibration verification is done prior to each use utilizing NIST traceable weights.



STATE OF MISSOURI

DEPARTMENT OF HEALTH AND SENIOR SERVICES
BREATH ALCOHOL PROGRAM



PERMIT TYPE II

ADAM J. MICHELS

is hereby authorized to instruct and supervise operators, train instructors, inspect, calibrate, perform field service and repairs, and operate the following breath analyzer(s):

ALCO-SENSOR IV WITH PRINTER

for the determination of the alcoholic content of blood from a sample of expired air. Permit issued under the provisions of sections 577.020 through 577.041, RSMo and 306.111 through 306.119 RSMo.

DATE	12/11/2023	/ (ike / lassmi
		DIRECTOR OF STATE PUBLIC HEALTH LABORATORY
NUMBER	230300	
EXPIRES 12/11/	12/11/2025	Davla J. Nichelson
		DIRECTOR OF DEPARTMENT OF HEALTH AND SENIOR SERVICES

MO 580-0771 (6-10)

LAB-4 (R6-10)



The named cardholder is authorized to operate an evidential breath alcohol instrument for the determination of the alcoholic content in breath form of expired air in Missouri.

Operator MICHELS, ADAM Permit No 230300

