RECEIVED

By Tracy Crews at 6:29 pm, Dec 29, 2024



MISSOURI DEPARTMENT OF HEALTH AND SENIOR SERVICES STATE PUBLIC HEALTH LABORATORY BREATH ALCOHOL PROGRAM

INTOX DMT MAINTENANCE REPORT

REPORT #1

ANTONIO CONTRACTOR ANTONIO CONTR			
Complete this report at the time of the regular mont Complete this report whenever the instrument is set Retain the original and send a copy within 15 days	rviced or repaired and whenever it i	s placed into service.	
INTOX DMT SN NAME OF AGENCY Peculiar Police	Peculiar Police Department		
LOCATION OF INSTRUMENT (STREET AND CITY) 224 N. Main Peculiar, MO 64078		TIME OF INSPECTION 13:12:10	
CHECKLIST: Place a mark in the box by each item if found to be satisfactory or is operating within established limits. (Write in observed values where determined). Unmarked items must be corrected before using instrument.			
☑ DIAGNOSTIC RECORD			
DATE AND TIME 12/29/2024 13:12:13 ☐ DETECTOR			
☑ PROGRAM	☑ PROGRAM ☑ FILTER 1		
☑ SAMPLE CHAMBER 48.7°C	SAMPLE CHAMBER 48.7°C		
☑ BREATH TUBE 46.7°C ☑ FILTER 3			
☑ PUMP ☑ INTERNAL STANDARD			
BREATH ANALYZER ACCURACY STANDARDS			
☑ SIMULATOR STANDARD	☐ COMPRESSED ETHANOL-GAS MIXTURE		
STANDARD SUPPLIER GUTH	LOT#_23390	EXP. DATE	0/25/2025
	SIM. SN <u>SD2260</u>	SIM. NIST EXP DATE	01/19/2025
 CALIBRATION CHECK - (ONLY ONE STANDARD IS TO BE USED PER MAINTENANCE REPORT) Run three tests using a standard. All three tests must be within ±5% of the standard value and must have a spread of .005 or less. Mark the box corresponding to the standard being used. ☑ 0.10% STANDARD - MUST READ BETWEEN 0.095% AND 0.105% INCLUSIVE ☐ 0.08% STANDARD - MUST READ BETWEEN 0.076% AND 0.084% INCLUSIVE ☐ 0.04% STANDARD - MUST READ BETWEEN 0.038% AND 0.042% INCLUSIVE 			
	EST 2: 0.098	TEST 3: 0.099	
□ PERFORM R.F.I. TEST			
INDICATE THE NUMBER OF BREATH TESTS IN THE FOLLOWING RANGES SINCE THE LAST MAINTENANCE REPORT:			
REFUSALS: 5 004: 41	0509: 1 .1014: 0	.1519: 0	OVER .19: 0
LIST ANY NEW PARTS AND DESCRIBE ANY ALTERATION OR MODIFIC ESTABLISHED LIMITS (USE OTHER SIDE IF NECESSARY)	CATION THAT WAS MADE TO RESTORE THE IN	STRUMENT TO OPERATE SATISFACTORILY AF	ND WITHIN
INSPECTING OFFICER			
PRINT FULL NAME ALEXANDER B GOODMAN			
TYPE II PERMIT WANTER CALL		LEPHONE NUMBER 316-779-5102	
RETURN COMPLETED REPORT TO THE Breath Alcohol Program, Missouri Department of Health and Senior Services by mail, fax, or email			



CERTIFICATE OF ANALYSIS

Certified Alcohol Reference Solution for Simulator

Random Samples of Lot Number 23390 of Alcohol Reference Solution for Simulator were analyzed by gas chromatography on October 18, 2023, using a Perkin Elmer Gas Chromatograph Autosystem XL S/N: 610N9030209, and found to contain 0.1207% (w/vol) ethyl alcohol. The expiration date for this lot number is October 17, 2025 at 11:59 PM.

When used in a calibrated Simulator, operating at 34°C +/- $\chi^{2}^{\circ}\text{C}$, this solution will give a breath alcohol analysis instrument reading of 0.100 g/210L +/- 3%

The alcohol and water used in this solution were free of test interfering substances.

Ted L. Pauley, President GUTH LABORATORIES, INC.



e (* - 11 3£ 2025

STATE OF MISSOURI

DEPARTMENT OF HEALTH AND SENIOR SERVICES
BREATH ALCOHOL PROGRAM



PERMIT TYPE II

ALEXANDER B. GOODMAN

