



MISSOURI DEPARTMENT OF HEALTH AND SENIOR SERVICES
 STATE PUBLIC HEALTH LABORATORY
 BREATH ALCOHOL PROGRAM
DATAMASTER MAINTENANCE REPORT

received 7/14/14-cd ^{REPORT #6}

Complete this report at the time of the regular monthly preventive maintenance check (not to exceed 35 days).
 Complete this report whenever the instrument is serviced or repaired and whenever it is placed into service.
 Retain the original and send a copy within 15 days to the Breath Alcohol Program, DHSS.

REVIEWED
 By Carol Day at 1:39 pm, Jul 28, 2014

DATAMASTER SN 201194	NAME OF AGENCY Desloge Police Department	DATE OF INSPECTION 07/03/2014
LOCATION OF INSTRUMENT (STREET AND CITY) 1000 N. Desloge Dr., Desloge, MO 63601		TIME OF INSPECTION 6:04 pm

CHECKLIST: Place a mark in the box by each item if found to be satisfactory or if operating within established limits. (Write in observed values where determined.) Unmarked items must be corrected before using instrument.

<input checked="" type="checkbox"/> DIAGNOSTIC CHECK (PRINTOUT ATTACHED)	DATE AND TIME (from printout) 07/03/2014 18:04
<input checked="" type="checkbox"/> COMPUTER	<input checked="" type="checkbox"/> DETECTOR
<input checked="" type="checkbox"/> PROGRAM	<input checked="" type="checkbox"/> FILTERS
<input checked="" type="checkbox"/> HEATERS SAMPLE CHAMBER 48°C	<input checked="" type="checkbox"/> QUARTZ STANDARD
<input checked="" type="checkbox"/> FLOW DETECTOR	<input checked="" type="checkbox"/> CALIBRATION
<input checked="" type="checkbox"/> PUMP HIGH SPEED	<input checked="" type="checkbox"/> PRINTER

INDICATOR LIGHTS

SIMULATOR SOLUTION SUPPLIER Guth Laboratories LOT # 14030 EXP. DATE 01/22/2014

SIMULATOR TEMP (34°C ± 0.2°C) 34 °C SIMULATOR SN 093752 EXP. DATE 01/09/2015

CALIBRATION CHECK - (ONLY ONE STANDARD IS TO BE USED PER MAINTENANCE REPORT)

Run three tests using a standard solution. All three tests must be within ±5% of the standard value and must have a spread of .005 or less. Mark the box corresponding to the standard solution being used. (PRINTOUT ATTACHED)

0.100% STANDARD - MUST READ BETWEEN 0.095% AND 0.105% INCLUSIVE
 0.080% STANDARD - MUST READ BETWEEN 0.076% AND 0.084% INCLUSIVE
 0.040% STANDARD - MUST READ BETWEEN 0.038% AND 0.042% INCLUSIVE

TEST 1 \bullet .098	TEST 2 \bullet .098	TEST 3 \bullet .098
-----------------------	-----------------------	-----------------------

PERFORM R.F.I. TEST (PRINTOUT ATTACHED)

INDICATE THE NUMBER OF BREATH TESTS IN THE FOLLOWING RANGES SINCE THE LAST MAINTENANCE REPORT:
 (DO NOT INCLUDE SELF-ADMINISTERED TESTS)

REFUSALS	2	(0-.04)	0	(.05-.09)	0	(.10-.14)	0	(.15-.19)	1	OVER .19	0
----------	---	---------	---	-----------	---	-----------	---	-----------	---	----------	---

LIST ANY NEW PARTS AND DESCRIBE ANY ALTERATION OR MODIFICATION THAT WAS MADE TO RESTORE THE INSTRUMENT TO OPERATE SATISFACTORILY AND WITHIN ESTABLISHED LIMITS (USE OTHER SIDE IF NECESSARY).

This instrument is operating within DOH specs.

INSPECTING OFFICER

SIGNATURE <i>G. Brad Judge</i>	PRINT FULL NAME G. Brad Judge
TYPE II PERMIT NUMBER/EXPIRATION DATE 230003, 01/08/2015	TELEPHONE NUMBER (573) 431-1463

RETURN COMPLETED REPORT TO THE: Breath Alcohol Program, MO Department of Health and Senior Services, Southeast District Office
 2875 James Blvd.
 Poplar Bluff, MO 63901



GUTH LABORATORIES, INC.

690 NORTH 67th STREET • HARRISBURG, PA 17111-4511 • TELEPHONE: 717-564-5470

CERTIFICATE OF ANALYSIS

Certified Alcohol Reference Solution for Simulator

Random Samples of Lot Number 14030 of Alcohol Reference Solution for Simulator were analyzed by gas chromatography on January 22, 2014, using a Perkin Elmer Gas Chromatograph Autosystem XL S/N: 610N9030209, and found to contain 0.1215% (w/vol) ethyl alcohol. The expiration date for this lot number is January 20, 2016 at 11:59 PM.

When used in a calibrated Simulator, operating at 34°C +/- .2°C, this solution will give a breath alcohol analysis instrument reading of 0.100 g/210L +/- 3%.

The alcohol and water used in this solution were free of test interfering substances.

Ted L. Pauley, President
GUTH LABORATORIES, INC.

NIST Traceability:

Testing was conducted using Cerilliant Reference Standard lot number FN122211-02 whose values are traceable to NIST.

All balances are calibrated annually by an outside agency using NIST traceable weights. Calibration verification is done prior to each use utilizing NIST traceable weights.



STATE OF MISSOURI
 DEPARTMENT OF HEALTH AND SENIOR SERVICES
 BREATH ALCOHOL PROGRAM



**PERMIT
 TYPE II**

BRAD JUDGE

In hereby authorized to perform and operate the following breath alcohol program:

DATAMASTER

For the Department of the State Health Laboratory, located at 577.026 through 577.041, 55NCS and 506.111 through 506.115. Renewal date: 01/06/2013.

DATE 01/06/2013

NUMBER 230003

EXPIRES 01/06/2015

MO 690-0771 (8-10)

DIRECTOR OF STATE PUBLIC HEALTH LABORATORY

Adelle Shedd

DIRECTOR OF DEPARTMENT OF HEALTH AND SENIOR SERVICES

LAB-4 (8-10)

FACE THIS SIDE DOWN - THIS EDGE IN FIRST

BAC DataMaster
Evidence Ticket

STATE OF MISSOURI
DESLOGE POLICE DEPARTMENT

BAC DATAMASTER SERIAL NUMBER 201194
07/03/14

TESTING OFFICER:

JUDGE/BRAD
OFFICER I.D.: 406
PERMIT NUMBER: 230003
EXPIRATION DATE: 01/08/15
MISCELLANEOUS DATA:
GUTH .100 SOL,
LOT 14030

--- SUPERVISOR MODE ---

BLANK TEST	.000	18:06
INTERNAL STANDARD	VERIFIED	18:06
EXTERNAL STANDARD	.098	18:07
BLANK TEST	.000	18:07
EXTERNAL STANDARD	.098	18:08
BLANK TEST	.000	18:09
EXTERNAL STANDARD	.098	18:09
BLANK TEST	.000	18:10

N = 3
SIM. = .1
AVG. = .098

FACE THIS SIDE DOWN - THIS EDGE IN FIRST

BAC DataMaster
Evidence Ticket

STATE OF MISSOURI
DESLOGE POLICE DEPARTMENT

BAC DATAMASTER SERIAL NUMBER 201194
07/03/14
18:04

--- DIAGNOSTIC CHECK ---

COMPUTER:	OKAY
PROGRAM (04-07-2009):	OKAY
HEATERS	
SAMPLE CHAMBER:	48c
FLOW DETECTOR:	OKAY
PUMP	
HIGH SPEED:	OKAY
DETECTOR:	OKAY
FILTERS:	OKAY
QUARTZ STANDARD:	OKAY
CALIBRATION:	OKAY

PRINTER TEST

!"#\$%&'()*+,-./0123456789:;<=>?@ABCDEFGHIJ
KLMNOPQRSTUVWXYZ[\]^_`abcdefghijklmnopqr
stuvwxyz{|}~

OPERATOR SIGNATURE

P.D. Judge 406

OPERATOR SIGNATURE

P.D. Judge 406

Stock No.

REORDER ALL SUPPLIES FROM N.P.A.S.
P.O. BOX 1435, MANSFIELD, OH 44901

Card Stock No.
60021

REORDER ALL SUPPLIES FROM N.P.A.S.
P.O. BOX 1435, MANSFIELD, OH 44901

FACE THIS SIDE DOWN - THIS EDGE IN FIRST

BAC DataMaster
Evidence Ticket

STATE OF MISSOURI
DESLOGE POLICE DEPARTMENT

BAC DATAMASTER SERIAL NUMBER 201194
07/03/14

ARREST TIME: 17:00
SUBJECT NAME:
JUDGE/BRADLEY
DOB: 07/29/91 SEX: M
STATE/D.L.: MO/1234567
ARRESTING OFFICER:
JUDGE/BRAD
OFFICER I.D.: 407
TESTING OFFICER:
JUDGE/BRAD
OFFICER I.D.: 407
PERMIT NUMBER: 230003
EXPIRATION DATE: 01/08/15
MISCELLANEOUS DATA:
RFI CHK

--- BREATH ANALYSIS ---

BLANK TEST	.000	18:12
INTERNAL STANDARD	VERIFIED	18:12
RADIO INTERFERENCE		

OPERATOR SIGNATURE

P.O. Judge 406

Card Stock No.
60021

REORDER ALL SUPPLIES FROM N.P.A.S.
P.O. BOX 1435, MANSFIELD, OH 44901