



MISSOURI DEPARTMENT OF HEALTH AND SENIOR SERVICES
 STATE PUBLIC HEALTH LABORATORY
 BREATH ALCOHOL PROGRAM
DATAMASTER MAINTENANCE REPORT

RECEIVED
 By Carol Day at 8:29 am, Sep 10, 2014

REPORT #6

Complete this report at the time of the regular monthly preventive maintenance check (not to exceed 35 days).
 Complete this report whenever the instrument is serviced or repaired and whenever it is placed into service.
 Retain the original and send a copy within 15 days to the Breath Alcohol Program, DHSS.

DATAMASTER SN 201252	NAME OF AGENCY Vinita Park Police Department	DATE OF INSPECTION 09/10/2014
LOCATION OF INSTRUMENT (STREET AND CITY) 8473 Midland Blvd, Vinita Park, MO 63114		TIME OF INSPECTION 3:32 am

CHECKLIST: Place a mark in the box by each item if found to be satisfactory or if operating within established limits. (Write in observed values where determined.) Unmarked items must be corrected before using instrument.

<input checked="" type="checkbox"/> DIAGNOSTIC CHECK (PRINTOUT ATTACHED)	DATE AND TIME (from printout) 09/10/14, 03:32
<input checked="" type="checkbox"/> COMPUTER	<input checked="" type="checkbox"/> DETECTOR
<input checked="" type="checkbox"/> PROGRAM	<input checked="" type="checkbox"/> FILTERS
<input checked="" type="checkbox"/> HEATERS SAMPLE CHAMBER 34 °C	<input checked="" type="checkbox"/> QUARTZ STANDARD
<input checked="" type="checkbox"/> FLOW DETECTOR	<input checked="" type="checkbox"/> CALIBRATION
<input checked="" type="checkbox"/> PUMP HIGH SPEED	<input checked="" type="checkbox"/> PRINTER

<input checked="" type="checkbox"/> INDICATOR LIGHTS
<input checked="" type="checkbox"/> SIMULATOR SOLUTION SUPPLIER <u>Guth Laboratories</u> LCT # <u>14030</u> EXP. DATE <u>01/20/2016</u>
<input checked="" type="checkbox"/> SIMULATOR TEMP (34°C ± 0.2°C) <u>34</u> °C SIMULATOR SN <u>SD2308</u> EXP. DATE <u>02/12/2015</u>

CALIBRATION CHECK - (ONLY ONE STANDARD IS TO BE USED PER MAINTENANCE REPORT)

Run three tests using a standard solution. All three tests must be within ±5% of the standard value and must have a spread of .005 or less. Mark the box corresponding to the standard solution being used. (PRINTCUT ATTACHED)

0.100% STANDARD - MUST READ BETWEEN 0.095% AND 0.105% INCLUSIVE
 0.080% STANDARD - MUST READ BETWEEN 0.076% AND 0.084% INCLUSIVE
 0.040% STANDARD - MUST READ BETWEEN 0.038% AND 0.042% INCLUSIVE

TEST 1 → .100	TEST 2 → .101	TEST 3 → .100
---------------	---------------	---------------

PERFORM R.F.I. TEST (PRINTOUT ATTACHED)

INDICATE THE NUMBER OF BREATH TESTS IN THE FOLLOWING RANGES SINCE THE LAST MAINTENANCE REPORT:
 (DO NOT INCLUDE SELF-ADMINISTERED TESTS)

REFUSALS	1	(0-.04)	0	(.05-.09)	0	(.10-.14)	1	(.15-.19)	1	OVER .19	0
----------	---	---------	---	-----------	---	-----------	---	-----------	---	----------	---

LIST ANY NEW PARTS AND DESCRIBE ANY ALTERATION OR MODIFICATION THAT WAS MADE TO RESTORE THE INSTRUMENT TO OPERATE SATISFACTORILY AND WITHIN ESTABLISHED LIMITS (USE OTHER SIDE IF NECESSARY).

INSPECTING OFFICER	
SIGNATURE <i>Sgt J. Martin 220</i>	PRINT FULL NAME J. Martin DSN: 220
TYPE II PERMIT NUMBER/EXPIRATION DATE 230096, 05/28/2015	TELEPHONE NUMBER (314) 428-7373

RETURN COMPLETED REPORT TO THE: Breath Alcohol Program, MO Department of Health and Senior Services, Southeast District Office
 2875 James Blvd.
 Poplar Bluff, MO 63901



GUTH LABORATORIES, INC.

690 NORTH 67th STREET • HARRISBURG, PA 17111-4511 • TELEPHONE: 717-564-5470

CERTIFICATE OF ANALYSIS

Certified Alcohol Reference Solution for Simulator

Random Samples of Lot Number 14030 of Alcohol Reference Solution for Simulator were analyzed by gas chromatography on January 22, 2014, using a Perkin Elmer Gas Chromatograph Autosystem XL S/N: 610N6030209, and found to contain 0.1215% (w/vol) ethyl alcohol. The expiration date for this lot number is January 20, 2016 at 11:59 PM.

When used in a calibrated Simulator, operating at 34°C +/- .2°C, this solution will give a breath alcohol analysis instrument reading of 0.100 g/210L +/- 3%.

The alcohol and water used in this solution were free of test interfering substances.

Ted L. Pauley, President
GUTH LABORATORIES, INC.

NIST Traceability:

Testing was conducted using Cerilliant Reference Standard lot number FN122211-02 whose values are traceable to NIST.

All balances are calibrated annually by an outside agency using NIST traceable weights. Calibration verification is done prior to each use utilizing NIST traceable weights.

Face This Side Down - This Edge In First

BAC DataMaster Evidence Ticket

STATE OF MISSOURI
DNC DATA MASTER SERIAL NUMBER 281282
09/10/14

TESTING OFFICER:
PARTING, J
OFFICER T.D.A. 2828
PERMIT NUMBER: 280996
EXPIRATION DATE: 05/28/15
MISCELLANEOUS 08101

--- SUPERVISOR NOTE ---

BLANK TEST	.000	03:03
INTERNAL STANDARD	VERIFIED	03:05
EXTERNAL STANDARD	.100	03:35
BLANK TEST	.000	03:36
INTERNAL STANDARD	.101	03:36
BLANK TEST	.000	03:37
EXTERNAL STANDARD	.100	03:37
BLANK TEST	.000	03:38

11-2
281, 7-1
281282

Operator Signature

Sgt J. Martin DSU: 220

Face This Side Down - This Edge In First

BAC DataMaster Evidence Ticket

STATE OF MISSOURI
DNC DATA MASTER SERIAL NUMBER 281282
09/10/14
03:38

--- SUPPLEMENTAL CHECK ---

COMPUTER:	DNA
PROGRAM (CPL-N7-280996)	PKRY
MEASURES	
SAMPLE NUMBER:	434
FLOW DETECTOR:	DRAY
PUMP	
FLOW SPEED:	DRAY
REACTOR:	DRAY
FILTERS:	DRAY
MURKIN STANDARD:	DRAY
CALIBRATION:	DRAY

PRINTER TEST

11-2
281, 7-1
281282

Operator Signature

Sgt J. Martin DSU

Printed on recycled paper

Face This Side Down - This Edge In First

BAC DataMaster Evidence Ticket

STATE OF MISSOURI
BAC DATA MASTER SERIAL NUMBER 801252
09/18/14

ARREST TIME: 08:00

SUBJECT NAME:

DOB: 03/21/66

SEX: M

STATE: MO, MOJESBTRIVE

ARRESTING OFFICER:

MARTIN, J

OFFICER I.D. #: 828

TESTING OFFICER:

MERTON, J

OFFICER I.D. #: 828

PERMIT NUMBER: 800098

EXPIRATION DATE: 06/29/10

MISCELLANEOUS DATA

SET CODE:

--- BREATH ANALYSIS ---

BRAND TEST

INTERNAL STANDARD

DRUG INTERFERENCE

0.000

NEGATIVE

0.000

0.000

Operator Signature

Sgt. J. Martin DSN: 7

Printed on recycled paper with agrifib-based inks

DNISU 2206-00



STATE OF MISSOURI
DEPARTMENT OF HEALTH AND SENIOR SERVICES
BREATH ALCOHOL PROGRAM

2

**PERMIT
TYPE II**

JESSE MARTIN

is hereby authorized to instruct and supervise operators, train instructors, inspect, calibrate, perform field service and repairs, and operate the following breath analyzer(s):

DATAMASTER

for the determination of the alcoholic content of blood from a sample of expired air. Permit issued under the provisions of sections 577.020 through 577.041, RSMo and 306.111 through 306.119 RSMo.

DATE 05/28/2013

NUMBER 230096

EXPIRES 05/28/2015

MOI 500-0771 (8-10)

W's

DIRECTOR OF STATE PUBLIC HEALTH LABORATORY

Gal Vandenbergh
Acting Director

DIRECTOR OF DEPARTMENT OF HEALTH AND SENIOR SERVICES

LAB-4 (88-10)