



MISSOURI DEPARTMENT OF HEALTH AND SENIOR SERVICES
 STATE PUBLIC HEALTH LABORATORY
 BREATH ALCOHOL PROGRAM
DATAMASTER MAINTENANCE REPORT

REPORT #6

Complete this report at the time of the regular monthly preventive maintenance check (not to exceed 30 days).
 Complete this report whenever the instrument is serviced or repaired and whenever it is placed into service.
 Retain the original and send a copy within 15 days to the Breath Alcohol Program, DHSS.

RECEIVED
 By Carol Day at 1:52 pm, Mar 27, 2015

DATAMASTER SN 201211	NAME OF AGENCY Bella Villa P.D.	DATE OF INSPECTION 03/27/2015
-------------------------	------------------------------------	----------------------------------

LOCATION OF INSTRUMENT (STREET AND CITY) 751 Avenue H, Bella Villa, Mo 63125	TIME OF INSPECTION 6:21 am
---	-------------------------------

CHECKLIST: Place a mark in the box by each item if found to be satisfactory or if operating within established limits. (Write in observed values where determined.) Unmarked items must be corrected before using instrument.

<input checked="" type="checkbox"/> DIAGNOSTIC CHECK (PRINTOUT ATTACHED)	DATE AND TIME (from printout) <u>03/27/2015, 06:21 Hours</u>
<input checked="" type="checkbox"/> COMPUTER	<input checked="" type="checkbox"/> DETECTOR
<input checked="" type="checkbox"/> PROGRAM	<input checked="" type="checkbox"/> FILTERS
<input checked="" type="checkbox"/> HEATERS SAMPLE CHAMBER <u>49 °C</u>	<input checked="" type="checkbox"/> QUARTZ STANDARD
<input checked="" type="checkbox"/> FLOW DETECTOR	<input checked="" type="checkbox"/> CALIBRATION
<input checked="" type="checkbox"/> PUMP HIGH SPEED	<input checked="" type="checkbox"/> PRINTER

INDICATOR LIGHTS

SIMULATOR SOLUTION SUPPLIER Guth LOT # 14110 EXP. DATE 05/01/2014

SIMULATOR TEMP (34°C ± 0.2°C) 34.0 °C SIMULATOR SN DS1095 EXP. DATE 04/17/2015

CALIBRATION CHECK - (ONLY ONE STANDARD IS TO BE USED PER MAINTENANCE REPORT)

Run three tests using a standard solution. All three tests must be within ±5% of the standard value and must have a spread of .005 or less. Mark the box corresponding to the standard solution being used. (PRINTOUT ATTACHED)

0.100% STANDARD - MUST READ BETWEEN 0.095% AND 0.105% INCLUSIVE
 0.080% STANDARD - MUST READ BETWEEN 0.076% AND 0.084% INCLUSIVE
 0.040% STANDARD - MUST READ BETWEEN 0.038% AND 0.042% INCLUSIVE

TEST 1 <input checked="" type="checkbox"/> 0.101	TEST 2 <input checked="" type="checkbox"/> 0.101	TEST 3 <input checked="" type="checkbox"/> 0.101
--	--	--

PERFORM R.F.I. TEST (PRINTOUT ATTACHED)

INDICATE THE NUMBER OF BREATH TESTS IN THE FOLLOWING RANGES SINCE THE LAST MAINTENANCE REPORT:
 (DO NOT INCLUDE SELF-ADMINISTERED TESTS)

REFUSALS	0	(0-.04)	0	(.05-.09)	0	(.10-.14)	0	(.15-.19)	0	OVER .19	0
----------	---	---------	---	-----------	---	-----------	---	-----------	---	----------	---

LIST ANY NEW PARTS AND DESCRIBE ANY ALTERATION OR MODIFICATION THAT WAS MADE TO RESTORE THE INSTRUMENT TO OPERATE SATISFACTORILY AND WITHIN ESTABLISHED LIMITS (USE OTHER SIDE IF NECESSARY).

INSPECTING OFFICER

SIGNATURE 	PRINT FULL NAME Peter Palombo
---------------	----------------------------------

TYPE II PERMIT NUMBER/EXPIRATION DATE 240177 03/27/2015	TELEPHONE NUMBER (314) 638-8840
--	------------------------------------

RETURN COMPLETED REPORT TO THE: **Breath Alcohol Program, MO Department of Health and Senior Services, Southeast District Office
 2875 James Blvd.
 Poplar Bluff, MO 63901**



GUTH LABORATORIES, INC.

890 NORTH 87th STREET • HARRISBURG, PA 17111-4611 • TELEPHONE: 717-654-6470

CERTIFICATE OF ANALYSIS

Certified Alcohol Reference Solution for Simulator

Random Samples of Lot Number 14110 of Alcohol Reference Solution for Simulator were analyzed by gas chromatography on May 5, 2014, using a Perkin Elmer Gas Chromatograph Autosystem XL S/N: 610N9030209, and found to contain 0.1206% (w/vol) ethyl alcohol. The expiration date for this lot number is May 1, 2016 at 11:59 PM.

When used in a calibrated Simulator, operating at 34°C +/- .2°C, this solution will give a breath alcohol analysis instrument reading of 0.100 g/210L +/- 3%.

The alcohol and water used in this solution were free of test interfering substances.

Ted L. Pauley, President
GUTH LABORATORIES, INC.

NIST Traceability:

Testing was conducted using Cerilliant Reference Standard lot number FN122211-02 whose values are traceable to NIST.

All balances are calibrated annually by an outside agency using NIST traceable weights. Calibration verification is done prior to each use utilizing NIST traceable weights.

BAC DataMaster Evidence Ticket

BAC DataMaster Evidence Ticket

STATE OF MISSOURI
BAC DATAMASTER SERIAL NUMBER 201211
09/27/15

STATE OF MISSOURI
BAC DATAMASTER SERIAL NUMBER 201211
09/27/15

ARREST TIME: 00:00
SUBJECT NAME:
RADIO TEST
DOB: 11/11/11 SEX: M
STATE/D.L.: MO/1111
ARRESTING OFFICER:
OFFICER I.D. #:
OFFICER I.D. #:
OFFICER I.D. #:
PERMIT NUMBER:
EXPIRATION DATE:
MISCELLANEOUS DATA:

ARREST TIME: 00:00
SUBJECT NAME:
RADIO TEST
DOB: 11/11/11 SEX: M
STATE/D.L.: MO/1111
ARRESTING OFFICER:
OFFICER I.D. #:
OFFICER I.D. #:
OFFICER I.D. #:
PERMIT NUMBER:
EXPIRATION DATE:
MISCELLANEOUS DATA:

PLANK #1	2000	00:00
INTERNAL STORAGE	00:00	00:00
SUBJECT NAME	00:00	00:00
BLANK #1	00:00	00:00

Operator Signature



2208-02

Operator Signature



Face This Side Down - This Edge In First

BAC DataMaster Evidence Ticket

STATE OF MISSISSIPPI
BAC DATA MASTER SERIAL NUMBER 40111
#12707
#613

TOXICITY: OKAY
MODIFIED: OKAY
CHECKS:
INTERNAL: OKAY
EXTERNAL: OKAY
ALCOHOL: OKAY
HIGH DRESS: OKAY
DEFINITION: OKAY
FILTERS: OKAY
GLASS: OKAY
CALIBRATION: OKAY

RESISTER TEST

012345678910111213141516171819
ABCDEFGHIJKLMNOPQRSTUVWXYZ
abcdefghijklmnopqrstuvwxyz

Operator Signature

2200-02

Face This Side Down - This Edge In First

BAC DataMaster Evidence Ticket

STATE OF MISSISSIPPI
BAC DATA MASTER SERIAL NUMBER 40111
#12707
#613

TOXICITY: OKAY
MODIFIED: OKAY
CHECKS:
INTERNAL: OKAY
EXTERNAL: OKAY
ALCOHOL: OKAY
HIGH DRESS: OKAY
DEFINITION: OKAY
FILTERS: OKAY
GLASS: OKAY
CALIBRATION: OKAY

QUALITY CONTROL

Bottle Test	100	OK
INTERNAL STANDARD	100	OK
EXTERNAL STANDARD	100	OK
ALCOHOL TEST	100	OK
EXTERNAL STANDARD	100	OK
GLASS TEST	100	OK
EXTERNAL STANDARD	100	OK
BANK TEST	100	OK

RESISTER TEST
012345678910111213141516171819
ABCDEFGHIJKLMNOPQRSTUVWXYZ
abcdefghijklmnopqrstuvwxyz

Operator Signature

2200-



STATE OF MISSOURI
DEPARTMENT OF HEALTH AND SENIOR SERVICES
BREATH ALCOHOL PROGRAM

INSTRUMENT OPERATOR CARD

The named cardholder is authorized to operate an evidential breath alcohol instrument for the determination of the alcoholic content in breath form of expired air in Missouri.



Operator PALOMBO, PETER
Permit No 240177
Date Issued 4/22/2014 Date Expires 4/22/2016