1	4	4.
	8	$\langle \hat{\zeta} \rangle$
	Ŋ,	
A	ď.	

MISSOURI DEPARTMENT OF HEALTH AND SENIOR SERVICES STATE PUBLIC HEALTH LABORATORY

BREATH ALCOHOL PROGRAM RECEIVED INTOX DMT MAINTENANCE REPORT By Carol Day at 10:25 am, Sep 01, 2015 Complete this report at the time of the regular monthly preventive maintenance check (not to exceed 35 days); Complete this report whenever the instrument is serviced or repaired and whenever it is placed into service. Retain the original and send a copy within 15 days to the Breath Alcohol Program, DHSS. INTOX OMT SN DATE OF INSPECTION 500249 SIKESTON DPS 08/21/2015 OCATION OF INSTRUMENT (STREET AND CITY TIME OF INSPECTION 201 S. KINGSHIGHWAY SIKESTON,MO 63801 19:51:50 CHECKLIST: Place a mark in the box by each Item if found to be sallsfactory or is operating within established limits. (Write in observed values where determined). Unmarked items must be corrected before using instrument. ☑ DIAGNOSTIC RECORD DATE AND TIME 08/21/2015 19:51:52 ☑ DETECTOR ☑ PROGRAM FILTER 1 SAMPLE CHAMBER 49.0°C XI FILTER 2 BREATH TUBE 47.2°C XI FILTER 3 · Z PUMP INTERNAL STANDARD BREATH ANALYZER ACCURACY STANDARDS ☐ SIMULATOR STANDARD ☑ COMPRESSED ETHANOL-GAS MIXTURE STANDARD SUPPLIER INTOXIMATER LOT# AG422001 EXP. DATE 08/08/2016 SIMULATOR TEMP (34°C ± 0,2°C) SIMULATOR SN SIMULATOR EXP DATE CALIBRATION CHECK - (ONLY ONE STANDARD IS TO BE USED PER MAINTENANCE REPORT)
Run three tests using a standard. All three tests must be within ±5% of the standard value and must have a spread of .005 or less. Mark the box corresponding to the standard being used. ☑ 0.10% STANDARD - MUST READ BETWEEN 0.095% AND 0.105% INCLUSIVE ☐ 0.08% STANDARD - MUST READ BETWEEN 0.076% AND 0.084% INCLUSIVE ☐ 0.04% STANDARD - MUST READ BETWEEN 0.038% AND 0.042% INCLUSIVE TEST 1: 0,099 TEST 2: 0.099 TEST 3: 0,098 PERFORM R.F.I, TEST INDICATE THE NUMBER OF BREATH TESTS IN THE FOLLOWING RANGES SINCE THE LAST MAINTENANCE REPORT: REFUSALS: 0 .05 - .09: 1 .10-.14: 0 .15-.19: 0 OVER .19: 0 LIST ANY NEW PARTS AND DESCRIBE ANY ALTERATION OR MODIFICATION THAT WAS MADE TO RESTORE THE INSTRUMENT TO OPERATE SATISFACTORILY AND WITHIN ESTABLISHED LIMITS (USE OTHER SIDE IF NECESSARY) MAINTENANCE TEST INSPECTING OFFICER FRANKLIN C ADAMS TYPE A PERMIT NUMBER EXPIRATION DATE TELEPHONE NUMBER 240324 08/18/2016 573-471-6200 Breath Alcohol Program, MO Department of Health and Senior Services RETURN COMPLETED REPORT TO THE Southeast District Office

2875 James Bivd, Poplar Bluff, MO 63901



Alrgas USA LLC (LAB) 3500 Bernard Street St. Louis, Mo. 63103 Ph: (314) 533-3100 Fax: (314) 533-7328

Certificate of Analysis

Customer Name Intoximeters, Inc. 2081 Craig Road St. Louis, Mo 63146 Test Date: 11-Aug-2014

Lot# AG422001

Exp. Date 8-Aug-2016 Cyl. Type

Component Ethanol

Nitrogen

Certified Concentration 0.100 ± 2% BrAC (260 ppm)

Balance

Certification Traceable to N.I.S.T. RGM Ethanol Standards:

Serial No.	Concentration 391.8 ppm 269.8 ppm 209.0 ppm 103.7 ppm	Serial No.	Concentration
EB0010581		EB0010603	392.5 ppm
EB0010570		EB0010559	258.9 ppm
EB0010285		EB0010595	208.9 ppm
EB0010561		EB0010562	104.9 ppm
EB0010681	62.22 ppm	EB0010679	62.94 ppm

Analytical Method:

NDIR

Digitally signed by Quelity Control Date: 2014.08.11 16:03:33-05:00 Reason: Dry gas standard certification of analysis Location: Argas USA LLC (Leb)

Analyst:

Rod Marsala

ISO 17025:2005 A2LA accredited. Certificate Number 2989.01



CERTIFICATE OF ANALYSIS

Certified Alcohol Reference Solution for Simulator

Random Samples of Lot Number 14220 of Alcohol Reference Solution for Simulator were analyzed by gas chromatography on September 25, 2014, using a Perkin Elmer Gas Chromatograph Autosystem XL S/N: 610N9030209, and found to contain 0.1209% (w/vol) ethyl alcohol. The expiration date for this lot number is September 24, 2016 at 11:59 PM.

When used in a calibrated Simulator, operating at 34°C +/- .2°C, this solution will give a breath alcohol analysis instrument reading of 0.100 g/210L +/- 3%.

The alcohol and water used in this solution were free of test interfering substances.

Ted L. Pauley, President GUTH LABORATORIES, INC.

NIST Traceability:

Testing was conducted using Cerilliant Reference Standard lot number FN08051301 whose values are traceable to NIST.

All balances are calibrated annually by an outside agency using NIST traceable weights. Calibration verification is done prior to each use utilizing NIST traceable weights.