



MISSOURI DEPARTMENT OF HEALTH AND SENIOR SERVICES
 STATE PUBLIC HEALTH LABORATORY
 BREATH ALCOHOL PROGRAM
DATAMASTER MAINTENANCE REPORT

RECEIVED

REPORT #6

By Carol Day at 8:40 am, Sep 09, 2015

Complete this report at the time of the regular monthly preventive maintenance. Complete this report whenever the instrument is serviced or repaired and whenever it is placed into service. Retain the original and send a copy within 15 days to the Breath Alcohol Program, DHSS.

DATAMASTER SN 950090	NAME OF AGENCY Phelps Co. Sheriff's Dept.	DATE OF INSPECTION 8-31-15
LOCATION OF INSTRUMENT (STREET AND CITY) 500 W. 2ND, Rolla, Mo 65401		TIME OF INSPECTION 1721

CHECKLIST: Place a mark in the box by each item if found to be satisfactory or if operating within established limits. (Write in observed values where determined.) Unmarked items must be corrected before using instrument.

<input checked="" type="checkbox"/> DIAGNOSTIC CHECK (PRINTOUT ATTACHED)	DATE AND TIME (from printout) 8-31-15 1721
<input type="checkbox"/> COMPUTER	<input checked="" type="checkbox"/> DETECTOR
<input checked="" type="checkbox"/> PROGRAM	<input checked="" type="checkbox"/> FILTERS
<input checked="" type="checkbox"/> HEATERS SAMPLE CHAMBER 49 °C	<input checked="" type="checkbox"/> QUARTZ STANDARD
<input checked="" type="checkbox"/> FLOW DETECTOR	<input checked="" type="checkbox"/> CALIBRATION
<input checked="" type="checkbox"/> PUMP HIGH SPEED	<input checked="" type="checkbox"/> PRINTER

<input checked="" type="checkbox"/> INDICATOR LIGHTS
<input checked="" type="checkbox"/> SIMULATOR SOLUTION SUPPLIER Guth Labs LOT # 14220 EXP. DATE 9-24-16
<input checked="" type="checkbox"/> SIMULATOR TEMP (34°C ± 0.2°C) 34 °C SIMULATOR SN SD1313 EXP. DATE 1-21-15

CALIBRATION CHECK - (ONLY ONE STANDARD IS TO BE USED PER MAINTENANCE REPORT)

Run three tests using a standard solution. All three tests must be within ±5% of the standard value and must have a spread of .005 or less. Mark the box corresponding to the standard solution being used. (PRINTOUT ATTACHED)

<input checked="" type="checkbox"/> 0.100% STANDARD - MUST READ BETWEEN 0.095% AND 0.105% INCLUSIVE
<input type="checkbox"/> 0.080% STANDARD - MUST READ BETWEEN 0.076% AND 0.084% INCLUSIVE
<input type="checkbox"/> 0.040% STANDARD - MUST READ BETWEEN 0.038% AND 0.042% INCLUSIVE

TEST 1 .097	TEST 2 .097	TEST 3 .097
--------------------	--------------------	--------------------

PERFORM R.F.I. TEST (PRINTOUT ATTACHED)

INDICATE THE NUMBER OF BREATH TESTS IN THE FOLLOWING RANGES SINCE THE LAST MAINTENANCE REPORT: (DO NOT INCLUDE SELF-ADMINISTERED TESTS)

REFUSALS -	(0-.04) -	(.05-.09) 4	(.10-.14) 1	(.15-.19) -	OVER .19 -
------------	-----------	--------------------	--------------------	-------------	------------

LIST ANY NEW PARTS AND DESCRIBE ANY ALTERATION OR MODIFICATION THAT WAS MADE TO RESTORE THE INSTRUMENT TO OPERATE SATISFACTORILY AND WITHIN ESTABLISHED LIMITS (USE OTHER SIDE IF NECESSARY).

N/A

INSPECTING OFFICER	
SIGNATURE [Signature]	PRINT FULL NAME Paul J. Lambert
TYPE II PERMIT NUMBER/EXPIRATION DATE 240292 6-27-16	TELEPHONE NUMBER 573-426-3860

RETURN COMPLETED REPORT TO THE: **Breath Alcohol Program, MO Department of Health and Senior Services, Southeast District Office
 2875 James Blvd.
 Poplar Bluff, MO 63901**

**BAC DataMaster
Evidence Ticket**

STATE OF MISSISSIPPI
PHELPS COUNTY SHERIFFS DEPT

END DATA MASTER SERIAL NUMBER 9500000
09/03/15
17:01

----- DIAGNOSTIC CHECK -----

COMPUTER: OKAY

PROGRAM (04-07-2009): OKAY

HEATERS
SAMPLE CHAMBER: 49C

FLOW DETECTOR: OKAY

PUMP
HIGH SPEED: OKAY

DETECTOR: OKAY

FILTERS: OKAY

QUARTZ STANDARD: OKAY

CALIBRATION: OKAY

PRINTER TEST

!#%&'()*+,-./0123456789:;<=>@ABCDEFGHIJK
LMNOPQRSTUVWXYZ[\]^_`abdefghijklmnop
qrstuvwxyz{|}~

Operator Signature



2014

BAC DataMaster Evidence Ticket

STATE OF MISSOURI
OSAGE COUNTY SHERIFFS DEPT

BAC DATA MASTER SERIAL NUMBER 955000
06/20/15

TESTING OFFICER:
LAMBERT/PAUL/J
OFFICER I.D.# 204
PERMIT NUMBER: 240200
EXPIRATION DATE: 06/27/16
MISCELLANEOUS DATA:

----- SUPERVISOR MODE -----

BLANK TEST	.000	17:04
INTERNAL STANDARD	VERIFIED	17:04
EXTERNAL STANDARD	.007	17:04
BLANK TEST	.000	17:05
INTERNAL STANDARD	.007	17:05
BLANK TEST	.000	17:05
EXTERNAL STANDARD	.007	17:05
BLANK TEST	.000	17:06
EXTERNAL STANDARD	.007	17:06
BLANK TEST	.000	17:07

N = 9
SUM = .1
AVG = .007

Operator Signature 

**BAC DataMaster
Evidence Ticket**

STATE OF MISSOURI
PHELPS COUNTY SHERIFFS DEPT

390 DATAMASTER SERIAL NUMBER 989990
06/01/15

ARREST TIME: 16:00
SUBJECT NAME:
JONE/TURBLE
DOB: 01/01/50 SEX: F
STATE/I.D.: MO/90101010
ARRESTING OFFICER:
LAMBERT/PAUL/J
OFFICER I.D. #: 204
TESTING OFFICER:
LAMBERT/PAUL/J
OFFICER I.D. #: 204
PERMIT NUMBER: 240292
EXPIRATION DATE: 06/27/16
MISCELLANEOUS DATA:

--- BREATH ANALYSIS ---

BLANK TEST 200 17:00
INTERNAL STANDARD VERIFIED 17:00
MIDIC INTERFERENCE

Operator Signature





STATE OF MISSOURI
DEPARTMENT OF HEALTH AND SENIOR SERVICES
BREATH ALCOHOL PROGRAM

2

PERMIT
TYPE II
PAUL J LAMBERT

is hereby authorized to instruct and supervise operators, train instructors, inspect, calibrate, perform field service and repairs, and operate the following breath analyzer(s):

DATAMASTER, INTOX DMT

for the determination of the alcoholic content of blood from a sample of expired air. Permit issued under the provisions of sections 577.020 through 577.041, RSMo and 306.111 through 306.119 RSMo.

DATE **3/31/2015**

NUMBER **250072**

EXPIRES **3/31/2017**

MO 589 0771 (6-13)

DIRECTOR OF STATE PUBLIC HEALTH LABORATORY

DIRECTOR OF DEPARTMENT OF HEALTH AND SENIOR SERVICES

LAB-4 (16-10)



STATE OF MISSOURI
DEPARTMENT OF HEALTH AND SENIOR SERVICES
BREATH ALCOHOL PROGRAM

INSTRUMENT OPERATOR CARD

The named cardholder is authorized to operate an evidential breath alcohol instrument for the determination of the alcoholic content in breath form of expired air in Missouri.



Operator **LAMBERT, PAUL**
Permit No **250072**
Date Issued **3/31/2015** Date Expires **3/31/2017**



GUTH LABORATORIES, INC.

690 NORTH 67th STREET • HARRISBURG, PA 17111-4511 • TELEPHONE: 717-684-6470

CERTIFICATE OF ANALYSIS

Certified Alcohol Reference Solution for Simulator

Random Samples of Lot Number 14220 of Alcohol Reference Solution for Simulator were analyzed by gas chromatography on September 25, 2014, using a Perkin Elmer Gas Chromatograph Autosystem XL S/N: 610N9030209, and found to contain 0.1209% (w/vol) ethyl alcohol. The expiration date for this lot number is September 24, 2016 at 11:59 PM.

When used in a calibrated Simulator, operating at 34°C +/- .2°C, this solution will give a breath alcohol analysis instrument reading of 0.100 g/210L +/- 3%.

The alcohol and water used in this solution were free of test interfering substances.

Ted L. Pauley, President
GUTH LABORATORIES, INC.

NIST Traceability:

Testing was conducted using Cerilliant Reference Standard lot number FN08051301 whose values are traceable to NIST.

All balances are calibrated annually by an outside agency using NIST traceable weights. Calibration verification is done prior to each use utilizing NIST traceable weights.