MISSOURI DEPARTMENT OF HEALTH AND SENIOR SERVICES STATE PUBLIC HEALTH LABORATORY BREATH ALCOHOL PROGRAM

INTOX DMT MAINTENANCE REPORT

Complete this report at the time of the regular monthly preventive maintenance of RECEIVED days). Complete this report whenever the instrument is serviced or repaired and whene Retain the original and send a copy within 15 days to the Breath Alcohol Program

· ·		•	•			
OX DMT SN NAME OF AGENCY Stone County Sheriff's Office			<u> </u>	DATE OF INSPECTION 06/08/2016		
LOCATION OF INSTRUMENT (STREET AND CITY) 110 S Maple, Galena, MO				TIME OF INSPECTION 08:13:54		
CHECKLIST: Place a ma	ark in the box by each ite . Unmarked items must	em if found to be satisfact be corrected before usi	etory or is operating wing instrument.	ithin established limits. (Wri	te in observed	
☑ DIAGNOSTIC RECO	ORD					
DATE AND TIME 06/08/2016 08:13:56						
☑ PROGRAM]	☑ FILTER 1	FILTER 1		
☑ SAMPLE CHAMBER 48.8°C ☑ FILTER 2						
☑ BREATH TUBE 47.7°C ☑ FILTER 3						
⊠ PUMP			☑ INTERNAL STANDARD			
BREATH ANALYZER ACCURACY STANDARDS						
☑ SIMULATOR STANDARD			☐ COMPRESSED ETHANOL-GAS MIXTURE			
STANDARD SUPPLIER GUTH		LOT#_	15120	EXP. DATE <u>04/2</u>	29/2017	
SIMULATOR TEMP (34°C ± 0.2°C) 34.0	SIMULAT	OR SN <u>SD2282</u>	SIMULATOR EXP DATE	12/07/2016	
 □ CALIBRATION CHECK - (ONLY ONE STANDARD IS TO BE USED PER MAINTENANCE REPORT) Run three tests using a standard. All three tests must be within ±5% of the standard value and must have a spread of .005 or less. Mark the box corresponding to the standard being used. □ 0.10% STANDARD - MUST READ BETWEEN 0.095% AND 0.105% INCLUSIVE □ 0.08% STANDARD - MUST READ BETWEEN 0.076% AND 0.084% INCLUSIVE □ 0.04% STANDARD - MUST READ BETWEEN 0.038% AND 0.042% INCLUSIVE 						
		TEST 2: 0.100		TEST 3: 0,100		
☑ PERFORM R.F.I. TEST		12012.0.100		1.22.3.6.100		
INDICATE THE NUMBER OF BREATH TESTS IN THE FOLLOWING RANGES SINCE THE LAST MAINTENANCE REPORT:						
	004: 0	.0509: 0	.1014: 3	.1519: 0	OVER .19: 1	
	RIBE ANY ALTERATION OR MOD		II.	TO OPERATE SATISFACTORILY AND		
INCRECTING OFFICER						
INSPECTING OFFICER SIGNATURE			PRINT FULL NAME	PRINT FULL NAME DANIEL FLETCHER		
TYPE II PERMIT NUMBER	EXPIRATION DATE	TELEPHONE N	JMBER	<u> </u>		
240406		11/20/2016	417-357			
RETURN COMPLETED	S	Breath Alcohol Program, Southeast District Office 1875, James Blvd, Poplar		ealth and Senior Services		



STATE OF MISSOURI

DEPARTMENT OF HEALTH AND SENIOR SERVICES BREATH ALCOHOL PROGRAM

TYPE |



DANIEL D FLETCHER

is hereby authorized to instruct and supervise operators, train instructors, inspect, calibrate, perform field service and repairs, and operate the following breath analyzer(s):

ALCO-SENSOR IV WITH PRINTER, INTOX DMT

for the determination of the alcoholic content of blood from a sample of expired air. Permit issued under the provisions of sections 577.020 through 577.041, RSMo and 306.111 through 306.119 RSMo.

11/20/2014 DATE

NUMBER 240406

EXPIRES 11/20/2016

MO 580-0771 (6-1G)

Gal Veroland

DIRECTOR OF STATE PUBLIC HEALTH LABORATORY

2 2

DIRECTOR OF DEPARTMENT OF HEALTH AND SENIOR SERVICES



CERTIFICATE OF ANALYSIS

Certified Alcohol Reference Solution for Simulator

Random Samples of Lot Number 15120 of Alcohol Reference Solution for Simulator were analyzed by gas chromatography on May 4, 2015, using a Perkin Elmer Gas Chromatograph Autosystem XL S/N: 610N9030209, and found to contain 0.1209% (w/vol) ethyl alcohol. The expiration date for this lot number is April 29, 2017 at 11:59 PM.

When used in a calibrated Simulator, operating at 34°C +/- .2°C, this solution will give a breath alcohol analysis instrument reading of 0.100 g/210L +/- 3%.

The alcohol and water used in this solution were free of test interfering substances.

Ted L. Pauley, President GUTH LABORATORIES, INC.

NIST Traceability:

Testing was conducted using Cerilliant Reference Standard lot number FN08051301 whose values are traceable to NIST.

All balances are calibrated annually by an outside agency using NIST traceable weights. Calibration verification is done prior to each use utilizing NIST traceable weights.