

# MISSOURI DEPARTMENT OF HEALTH AND SENIOR SERVICES STATE PUBLIC HEALTH LABORATORY

### ALCO-SENSOR IV WITH PRINTER MAINTENANCE REPORT

REPORT #7

| A STATE OF THE STA |  |  |                                      |   |                        |
|--|--|--|--------------------------------------|---|------------------------|
| Complete this report in duplicate at the tin Send copy to Department of Health and S   | ne of the regular monthly p<br>enior Services; retain origir   | reventative maintena<br>al in department file.                                     | ance check, ar                       | nd whenever ins                         | strument is repaired.  |
| ALCO SENSOR IV SN  | PRINTER SN   |  |                                      | DATE OF INSPE                           | CTION                  |
| 107985   | 019.3586.8   | 320  |                                      | 07-20                                   | 1-21                   |
| LOCATION OF INSTRUMENT (STREET AND CITY)   | 1911   |  |                                      | TIME OF INSPE                           |                        |
| 200 Inch. Are South  | MO 65753   |  |                                      | 0708                                    |                        |
| CHECKLIST: Place a mark in the box by  | each item if found to be sati  | sfactory or if operatir  | ng within estab                      | lished limits. (V                       | /rite in observed val- |
| ues where determined.) Unmarked items r  | nust be corrected before us  | sing instrument.   |                                      |   |                        |
| DO DIGITAL READOUT (ALL ELEMENTS   |  |  |                                      |   |                        |
| TEMPERATURE OF ALCO SENSOR   | (10°C - 40°C)  |  |                                      |   |                        |
| PRINTER WORKING PROPERLY   |  |  |                                      |   |                        |
| TIME AND DATE DISPLAYING PROF  | PERLY  |  |                                      |   |                        |
| BREATH ALCOHOL ACCURACY STAND  | ARDS   |  |                                      |   |                        |
| SIMULATOR SOLUTION   |  | ☐ COMPRESSE  | D ETHANOL-(                          | SAS MIXTURE                             |                        |
| STANDARD SUPPLIER GUTH LAB   | L  | от# <u>Q(080</u>   | EXP. DAT                             | E_3/8/2                                 | <u> </u>               |
| Simulator Temperature (346C  | 2 0.2°C) <u>34°C</u> SIMI  | JLATOR SN MR 35  | 8 <b>4</b> simu                      | JLATOR EXP D                            | ATE 02/25/2022         |
| CALIBRATION CHECK - (ONLY ONE  Run three tests using a standard solution less. Check the box corresponding to 0.100% STANDARD - MUST RE.  0.080% STANDARD - MUST RE.  0.040% STANDARD - MUST RE.   | tion. All three tests must be<br>the standard solution being<br>AD BETWEEN 0.095% and<br>AD BETWEEN 0.076% and | e within ±5% of the s<br>used. (PRINTOUT A<br>0.105% INCLUSIVE<br>0.084% INCLUSIVE | tandard value<br>ATTACHED)<br>E<br>E | and must have                           | a spread of .005 cr    |
| TEST 1 . (02   | TEST 2 . (O(   |  | TEST 3 🖝                             | .10(                                    |                        |
| PRFI DETECTOR OPERATING  |  |  |                                      | *************************************** | - Lineary today        |
| INDICATE THE NUMBER OF BREATH T<br>(DO NOT INCLUDE SELF-ADMINISTERI  | ESTS IN THE FOLLOWING  | 3 RANGES SINCE   | THE LAST MA                          | INTENANCE F                             | REPORT:                |
| A 1 7  | (.0509)  | (.1014)  | (.1519)                              | Ø lov                                   | ER .19)                |
| REFUSALS (004)   |  | `  |                                      | /                                       |                        |
| List any new parts and describe any alter established limits (use other side if necess   | ation or modification that v<br>sary).   | as made to restore   | the mendinen                         | to operate sat                          | goldstorny and mon     |
| P  |  |  |                                      |   |                        |
| INSPECTING OFFICER   |  |  |                                      |   |                        |
| SIGNATURE SIGNATURE  |  |  | PRINT NAME                           |   |                        |
| 1 /2/14  |  |  | LOREN                                | NYSTRON                                 | 1                      |
| TYPE IL PERMIT NUMBER/EXPIRATION DATE  | 2/23/2023  | NIII.  | TELEPHONE NUME                       | 42-5511                                 |                        |
| Return completed report to the: Breatl   | n Alcohol∤Program, MO Del  | partment of Health a   | nd Senior Serv                       | vices, Southeas                         | t District Office      |
|  | James Boulevard<br>r Bluff, MO 63901   |  |                                      |   |                        |

07/20/21 07:10 .000 Calibration Check: 22 07/20/21 07:10 .102 Subject Name 1 + St Subject I.D. Operator Name, I.D. Stran #6 Location AS IV Serial no: 107985 Version no: 532B TEST RECORD 00989 Date Temp VOID: RFI 12 07/20/21 07:15 Subject Name Subject I.D. Operator Name, I.D.

AS IV Serial no: 107985

TEST RECORD 00986

Date Time 210L

Version no: 532B

Teme

Air Blank:

AS IV Serial no: 107985
Version no: 532B

TEST RECORD 00987

Temp Date Time 210L

Air Blank:
07/20/21 07:12 .000
Calibration Check:
22 07/20/21 07:12 .101

Subject Name

Test
Subject I.D.

H.Z.

Operator Name, I.D.

Location

Subject Name

Location

AS IV Serial no: 107985
Version no: 532B

TEST RECORD 00988

Yes a constant of the constant of

Test Record 00990

Temp Date Time 210L

Air Blank:
 07/20/21 07:18 .000

Calibration Check:
 23 07/20/21 07:18 .000

Subject Name

Test
Subject I.D.

Blank
Operator Name, I.D.

Location

Data P. D. Sobrin

AS IV Serial no: 107985



### STATE OF MISSOURI

DEPARTMENT OF HEALTH AND SENIOR SERVICES
BREATH ALCOHOL PROGRAM



# PERMIT TYPE II

# LOREN NYSTROM

is hereby authorized to instruct and supervise operators, train instructors, inspect, calibrate, perform field service and repairs, and operate the following breath analyzer(s):


## ALCO-SENSOR IV WITH PRINTER

for the determination of the alcoholic content of blood from a sample of expired air. Permit issued under the provisions of sections 577.020 through 577.041, RSMo and 306.111 through 306.119 RSMo.

| DATE    | 2/23/2021 | we when   |  |  |  |
|---------|-----------|---|--|--|--|
|         |           | DIRECTOR OF STATE PUBLIC HEALTH LABORATORY  |  |  |  |
| NUMBER  | 210032    |   |  |  |  |
| EXPIRES | 2/23/2023 | DIRECTOR OF DEPARTMENT OF HEALTH AND SERVICES   |  |  |  |
|         |           | DIRECTOR OF DEFACING OF FIEADITIANS GENERAL SECTION OF THE ACT OF |  |  |  |

MO 580-0771 (6-10)

LAB-4 (R6-10)



The named cardholder is authorized to operate an evidential breath alcohol instrument for the determination of the alcoholic content in breath form of expired air In Missouri.

Operator NYSTROM, LOREN

Permit No 210032

Date Issued 2/23/2021 Date Expires 2/23/2023

