

MISSOURI DEPARTMENT OF HEALTH AND SENIOR SERVICES STATE PUBLIC HEALTH LABORATORY

RECEIVED

By Tracy Crews at 8:07 am, Dec 16, 2020

ALCO-SENSOR IV WITH PRINTER MAINTENANCE REPORT

REPORT #7

waters.			
Complete this report in duplicate at the time of Send copy to Department of Health and Senio	of the regular monthly preventative mainten or Services; retain original in department file	nance check, and whenever instrument is repaired.	
ALCO SENSOR IV SN 030810	PRINTER SN 84.9324.143	DATE OF INSPECTION 12/15/2020	
LOCATION OF INSTRUMENT (STREET AND CITY) 1207 South Main Street El Dorado Spring	gs	TIME OF INSPECTION ICI 141 PM	
CHECKLIST: Place a mark in the box by each	item if found to be satisfactory or if operati	ng within established limits. (Write in observed val-	
ues where determined.) Unmarked items mus	t be corrected before using instrument.		
DIGITAL READOUT (ALL ELEMENTS OF	PERATIONAL)		
▼ TEMPERATURE OF ALCO SENSOR (10°C - 40°C)			
✓ PRINTER WORKING PROPERLY			
TIME AND DATE DISPLAYING PROPERLY			
BREATH ALCOHOL ACCURACY STANDAR	DS		
☑ SIMULATOR SOLUTION ☐ COMPRESSED ETHANOL-GAS MIXTURE			
STANDARD SUPPLIER Guth Laborator	ries Inc LOT # 19370	EXP. DATE 12/09/2021	
SIMULATOR TEMPERATURE (34°C ± 0.	2°C)34.00 _ SIMULATOR SNMF	22946 SIMULATOR EXP DATE 02/13/2020	
less. Check the box corresponding to the s 0.100% STANDARD - MUST READ I 0.080% STANDARD - MUST READ I	standard solution being used. (PRINTOUT) BETWEEN 0.095% and 0.105% INCLUSIV BETWEEN 0.076% and 0.084% INCLUSIV BETWEEN 0.038% and 0.042% INCLUSIV	E E	
TEST 1 • 197	EST 2 . 097	TEST 3 - 696	
☑ RFI DETECTOR OPERATING			
INDICATE THE NUMBER OF BREATH TESTS IN THE FOLLOWING RANGES SINCE THE LAST MAINTENANCE REPORT: (DO NOT INCLUDE SELF-ADMINISTERED TESTS)			
REFUSALS (004)	(.0509) (.1014)	(.1519) O (OVER .19)	
List any new parts and describe any alteration or modification that was made to restore the instrument to operate satisfactorily and within established limits (use other side if necessary).			
NA			
	ž		
WORLD CHILD			
INSPECTING OFFICER SIGNATURE		PRINT NAME	
· 0 /		Jarrod Schiereck	
TYPE II PERMIT NUMBER/EXPIRATION DATE 200225/08/13/2022		TELEPHONE NUMBER (417) 876-2313	
2875 Jame	cohol Program, MO Department of Health at es Boulevard off, MO 63901	nd Senior Services, Southeast District Office	

AS IV Serial no: 030810 AS IV Serial no: 030810 Version not 532B Version no: 532B AS IV Serial no: 030810 Version no: 532B TEST RECORD 00333 TEST RECORD 00334 9/ TEST RECORD 00335 Time 2101. Time 210L 9/ Air Blank: Time 210L Air Blank: 12/15/20 14:41 .000 12/15/20 14:43 .000 Subject Test: Auto Air Blanks Calibration Check: 22 12/15/20 14:41 .000 12/15/20 14:45 .000 23 12/15/20 14:43 .097 Calibration Checks 23 12/15/20 14:45 .097 Subject Name Operator Hame, I.D. 109 Operator Name, I.D. 109 JARROD Schiered Operator Name, I.D. 109 1267 SOUTH MAIN 1207 SOUTH BARTH 1207 South MAIN Elbrada Springs Elbonale Jangs El Donado Grugo

> AS IV Serial no: 030810 Varsion no: 5328

> > TEST RECORD 00336

Time 210L

Air Blank: 12/15/20 14:46 .000 Calibration Check: 24 12/15/20 14:46 .096

Subject I.D.

Operator Name, I.D. 109 1207 South main LE Donalo for 1745 AS IV Serial no: \$30810 Version no: 532B

TEST RECORD 00337

Time 210L

12 12/15/20 14:48

ject Name

Subject I.D.

Operator Hame, I.I.109

Elborado Springs



CERTIFICATE OF ANALYSIS

Certified Alcohol Reference Solution for Simulator

Random Samples of Lot Number 19370 of Alcohol Reference Solution for Simulator were analyzed by gas chromatography on December 10, 2019, using a Perkin Elmer Gas Chromatograph Autosystem XL S/N: 610N9030209, and found to contain 0.1199% (w/vol) ethyl alcohol. The expiration date for this lot number is December 9, 2021 at 11:59 PM.

When used in a calibrated Simulator, operating at 34°C +/- .2°C, this solution will give a breath alcohol analysis instrument reading of 0.100 g/210L +/- 3%.

The alcohol and water used in this solution were free of test interfering substances.

Ted L. Pauley, President GUTH LABORATORIES, INC.

NIST Traceability:

Testing was conducted using Cerilliant Reference Standard lot number FN04271602 whose values are traceable to NIST.

All balances are calibrated annually by an outside agency using NIST traceable weights. Calibration verification is done prior to each use utilizing NIST traceable weights.



STATE OF MISSOURI

DEPARTMENT OF HEALTH AND SENIOR SERVICES
BREATH ALCOHOL PROGRAM



PERMIT TYPE II

JARROD SCHIERECK

is hereby authorized to instruct and supervise operators, train instructors, inspect, calibrate, perform field service and repairs, and operate the following breath analyzer(s):

ALCO-SENSOR IV WITH PRINTER

for the determination of the alcoholic content of blood from a sample of expired air. Permit issued under the provisions of sections 577.020 through 577.041, RSMo and 306.111 through 306.119 RSMo.

DATE8/13/2020	mun
	DIRECTOR OF STATE PUBLIC HEALTH LABORATORY
NUMBER 200225	
EXPIRES 8/13/2022	for I the literain
3-0771 (6-10)	DIRECTOR OF DEPARTMENT OF HEALTH AND SENIOR SERVICES
	LAB-4 (R6:10)

STATE OF MISSOURI
DEPARTMENT OF HEALTH AND SENIOR SERVICES
BREATH ALCOHOL PROGRAM

INSTRUMENT OPERATOR CARD

The named cardholder is authorized to operate an evidential breath alcohol instrument for the determination of the alcoholic content in breath form of expired at in Missouri.

Operator SC

SCHIERECK, JARROD 200225

Permit No 200225 Date Issued 8/13/2020

Date Expires 8/13/2022

