

### MISSOURI DEPARTMENT OF HEALTH AND SENIOR SERVICES STATE PUBLIC HEALTH LABORATORY

#### ALCO-SENSOR IV WITH PRINTER MAINTENANCE REPORT

REPORT #7

Complete this report in duplicate at the tim Send copy to Department of Health and Se	e of the regular monthly penior Services; retain origi	preventative mainten nal in department file	ance check, and	whenever instrume	ent is repaired.
ALCO SENSOR IV SN 111647	NAME OF AGENCY NORMANDY PD	)		DATE OF INSPECTION 04/28/2021	
LOCATION OF INSTRUMENT (STREET AND CITY) 7700 NATURAL BRIDGE, NORMAND	Y, MO, 63121			TIME OF INSPECTION 9:41 am	
CHECKLIST: Place a mark in the box by ea where determined.) Unmarked items must			within establishe	d limits. (Write in ol	served values
☐ DIGITAL READOUT (ALL ELEMENTS	OPERATIONAL)			:-	
☑ TEMPERATURE OF ALCO SENSOR	(10°C - 40°C)				
PRINTER WORKING PROPERLY			3 ,	,	
☑ TIME AND DATE DISPLAYING PROP	ERLY				
BREATH ALCOHOL ACCURACY STANDA	ARDS				
☑ SIMULATOR SOLUTION		☐ COMPRESSE	D ETHANOL-GA	S MIXTURE	
☑ STANDARD SUPPLIER GUTH	L	OT # <u>21080</u>	EXP. DATE	03/08/2023	_
☑ SIMULATOR TEMPERATURE (34°C ±	: 0.2°C)34.0 SII	M. SNSD161	5 SIM. NI	IST EXP DATE 07	/24/2021
Run three tests using a standard solutiless. Check the box corresponding to the 0.100% STANDARD - MUST REA 0.080% STANDARD - MUST REA 0.040% STANDARD - MUST REA	ne standard solution being .D BETWEEN 0.095% and .D BETWEEN 0.076% and	g used. (PRINTOUT / d 0.105% INCLUSIVI d 0.084% INCLUSIVI	ATTACHED) E E		
TEST 1 • .099	TEST 2099		TEST 3 ☞ .099	)	
RFI DETECTOR OPERATING				,*	
INDICATE THE NUMBER OF BREATH TE (DO NOT INCLUDE SELF-ADMINISTERE		G RANGES SINCE	THE LAST MAIN	TENANCE REPOR	RT:
REFUSALS 0 (004) 0	(.0509)	(.1014)	(.1519)	0 (OVER .19	9) 0
List any new parts and describe any altera established limits (use other side if necessa NONE AT THIS TIME	tion or modification that v			operate satisfacto	rily and within
INSPECTING OFFICER					
SIGNATURE Buch	7.66		PRINT NAME BRANDON W.	. LANGE	,
TYPE II PERMIT NUMBER/EXPIRATION DATE 290166 08/01/2021			TELEPHONE NUMBER (314) 385-330	0	
	Alcohol Program, MO Dep , fax, or email.	partment of Health ar	nd Senior Service	es, Southeast Distri	ct Office

AS IV Serial no: 111647 Version no: 532B TEST RECORD 00271 Temp Date Time 210L Air Blank: 04/28/21 09:41 .000 Calibration Check: 24 04/28/21 09:41 :099 Subject Name I twipm Subject I.D. AIG Orecator Name: I.D. 206 Location

AS IV Serial no: 111647 Version no: 532B TEST RECORD 00272 9/ Temp' Date Time 210L Air Blank: 04/28/21 09:43 .000 Calibration Check: 24 04/28/21 09:43 .099 Subject Name Maint 2 Subject I.D. NIX Operator Name, I.D. 2 amou Location

Normandy

AS IV Serial no: 111647 Version no: 532B TEST RECORD 00273 Date Time 210L Air Blank: 04/28/21 09:45 .000 Calibration Check: 25 04/28/21 09:45 .099 Subject Name Maint 3 Subject I.D. NA Oremator Name: I.D. Dann 206 Location Normandy

AS IV Serial no: 111647
Version no: 532B

TEST RECORD 00274

Temp Date Time 210L

VOID: RFI
12 04/28/21 09:46

Subject Name

Subject I.D.

Operator Name: I.D.

Location

Normany Ro

AS IV Serial no: 111647 Version no: 532B TEST RECORD 00275 Temp Date Time 210L Air Blank: 04/28/21 09:47 .000 Subject Test: Auto 26 04/28/21 09:47 .000 Subject Name ) anox Subject I.D. 766 Operator Name, I.D. 266 Location Normondy &D



## HEATH ALGERIAL PROGRAM

# PERMIT TYPE II



क्ष्मां कार्यका महावादी के मुन्ति कार्य अण्डापंड क्ष्मां कार्य के क्ष्मां के क्ष्मां के कार्य क

and operate the tolknown breath analyz	SENSOR IV WITH PRI	NTER	=
	and a faction of the	Trefles to endervoy all repursit	F
Figure de la company de la com	[AB	Fra January	2
DAE STRAIG	ने में में में में में में में में में मे	शुम्रह्मणब्द्धामध्यामा स्वतिष्या एवत्	
William Wille			<del>-</del> .
EXERES STITATE	DIRECTOR DE DEPÀF		<b>b</b> )
如马勒而功臣初			
,		•	
	•	,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,	
	DEFAULT OF MISSOURI		
Tieston Instan	uned racificided is authorized in operata मा स्पर्धे कार्ने breath alcohol कार्रोंक the determination मंत्रिक बदेन हैं जो दिन्दा माने कार्यकार में महर्मीए जेगा में स्पूर्णेल में		
元Mis Clust	ooot ofor lange brandon fino 290186 (seqed)112019 Dafe Equips: 8012021		•



#### CERTIFICATE OF ANALYSIS

Certified Alcohol Reference Solution for Simulator

Random Samples of Lot Number 21080 of Alcohol Reference Solution for Simulator were analyzed by gas chromatography on March 10, 2021, using a Perkin Elmer Gas Chromatograph Autosystem XL S/N: 610N9030209, and found to contain 0.1214% (w/vol) ethyl alcohol. The expiration date for this lot number is March 8, 2023 at 11:59 PM.

When used in a calibrated Simulator, operating at 34°C +/- .2°C, this solution will give a breath alcohol analysis instrument reading of 0.100 g/210L +/- 3%.

The alcohol and water used in this solution were free of test interfering substances.

Ted L. Pauley, President GUTH LABORATORIES, INC.

NIST Traceability:

Testing was conducted using Cerilliant Reference Standard lot number FN03052002 whose values are traceable to NIST.

All balances are calibrated annually by an outside agency using NIST traceable weights. Calibration verification is done prior to each use utilizing NIST traceable weights.