RECEIVED

By Tracy Crews at 10:15 am, Aug 19, 2024



MISSOURI DEPARTMENT OF HEALTH AND SENIOR SERVICE STATE PUBLIC HEALTH LABORATORY

ALCO-SENSOR IV WITH PRINTER MAINTENANCE REPORT

REPORT #7

C. ASHAMA.				
Complete this report in duplicate at the time o Send copy to Department of Health and Senío	f the regular monthly p r Services; retain origin	reventative maintenar al in deparlment file.		······································
ALCO SENSOR IV SN 087968	NAME OF AGENCY Missouri State Park Rangers			INSPECTION 2024
LOCATION OF INSTRUMENT (STREET AND CITY) 2390 S. Jefferson Ave, Lebanon, MO 65538			4:00 p	
CHECKLIST: Place a mark in the box by each	item if found to be satisf	actory or if operating v	vithin established limits	. (Write in observed values
where determined.) Unmarked items must be	corrected before using	instrument.		
DIGITAL READOUT (ALL ELEMENTS OF	PERATIONAL)			and a second
☑ TEMPERATURE OF ALCO SENSOR (10	°C - 40°C)			4)
PRINTER WORKING PROPERLY				
☑ TIME AND DATE DISPLAYING PROPER				
BREATH ALCOHOL ACCURACY STANDAR	DS			
SIMULATOR SOLUTION COMPRESSED ETHANOL-GAS MIXTURE				
STANDARD SUPPLIER Guth Laborator	ries L	OT # 23180	EXP. DATE 05/17/	2025
☑ SIMULATOR TEMPERATURE (34°C ± 0.	2°C) 34 SH	и. sn <u>MP6020</u>	SłM. NIST EX	P DATE 09/06/2024
0.100% STANDARD - MUST READ 0.080% STANDARD - MUST READ 0.040% STANDARD - MUST READ	BETWEEN 0.076% and	0.084% INCLUSIVE		
TEST 1 - ,101	TEST 2 ★ .101		TEST 3 ★ .101	
☑ RFI DETECTOR OPERATING				
INDICATE THE NUMBER OF BREATH TEST (DO NOT INCLUDE SELF-ADMINISTERED		G RANGES SINCE T	HE LAST MAINTENA	NCE REPORT:
REFUSALS 0 (0-,04) 0	(.0509)	(.1014)	(.1519)	(OVER +19)
List any new parts and describe any alteration established limits (use other side if necessary Instrument operates within guidelines.	n or modification that v		he instrument to opera	ate satisfactorily and within
INSPECTING OFFICER SIGNATURE			PRINT NAME	
TYPE II PERIOR NUMBURUEX PRAYTON DATE			David Miller	N ₂ U ₃
230220			(573) 508-2295	
Return completed report to the: Breath Ale	cohol Program, MO De ax, or email.	partment of Health an	d Senior Services, Sou	utheast District Office



CERTIFICATE OF ANALYSIS

Certified Alcohol Reference Solution for Simulator

Random Samples of Lot Number 23180 of
Alcohol Reference Solution for Simulator were analyzed by
gas chromatography on May 18, 2023, using a Perkin Elmer Gas Chromatograph
Autosystem XL S/N: 610N9030209, and found to contain 0.1220% (w/vol)
ethyl alcohol. The expiration date for this lot
number is May 17, 2025 at 11:59 PM.

When used in a calibrated Simulator, operating at 34°C +/- .2°C, this solution will give a breath alcohol analysis instrument reading of 0.100 g/210L +/- 3%.

The alcohol and water used in this solution were free of test interfering substances.

Ted L. Pauley, President GUTH LABORATORIES, INC.

NIST Traceability: Testing was conducted using Cerilliant Reference Standard lot number FN11172002 whose values are traceable to NIST. All balances are calibrated annually by an outside agency using NIST traceable weights.

Calibration verification is done prior to each use utilizing NIST traceable weights.



STATE OF MISSOURI

DEPARTMENT OF HEALTH AND SENIOR SERVICES **BREATH ALCOHOL PROGRAM**



PERMIT TYPE II

DAVID MILLER

is hereby authorized to instruct and supervise operators, train instructors, inspect, calibrate, perform field service and repairs, and operate the following breath analyzer(s):

ALCO-SENSOR IV WITH PRINTER

	from a sample of expired air. Permit issued under the provisions of sections
577.020 through 577.041, RSMo and 306.111 throug	
30	Mile Massur
DATE 10/11/2023	/ lake / lassing
	DIRECTOR OF STATE PUBLIC HEALTH LABORATORY
NUMBER 230220	
	Davla I. Nichelson
EXPIRES 10/11/2025	
	DIRECTOR OF DEPARTMENT OF HEALTH AND SENIOR SERVICES
MO 680-0771 (6-10)	LAU-4 (136-40)



The named cardholder is authorized to operate an evidential breath alcohol instrument for the determination of the alcoholic content in breath form of expired a

Operator MILLER, DAVID Permit No 230220

Date Issued 10/11/2023 Date Expires 10/11/2025

