



MISSOURI DEPARTMENT OF HEALTH AND SENIOR SERVICES STATE PUBLIC HEALTH LABORATORY BREATH ALCOHOL PROGRAM

INTOX DMT MAINTENANCE REPORT

REPORT #1

Complete this report at the time of the regular monthly preventive maintenance check (not to exceed 35 days). Complete this report whenever the instrument is serviced or repaired and whenever it is placed into service. Retain the original and send a copy within 15 days to the Breath Alcohol Program, DHSS. INTOX DMT SN 500194 LOCATION OF INSTRUMENT (STREET AND CITY) 108 N 5th St, Canton, MO CHECKLIST: Place a mark in the box by each item if found to be satisfactory or is operating within established limits. (Write in observed values where determined). Unmarked items must be corrected before using instrument. DIAGNOSTIC RECORD DATE AND TIME 04/03/2024 15:41:34 DETECTOR PROGRAM FILTER 1 SAMPLE CHAMBER 48.7°C FILTER 2 FILTER 3 PUMP INTERNAL STANDARD COMPRESSED ETHANOL-GAS MIXTURE
Canton Police Department Canton Police Department O4/03/2024
108 N 5th St, Canton, MO CHECKLIST: Place a mark in the box by each item if found to be satisfactory or is operating within established limits. (Write in observed values where determined). Unmarked items must be corrected before using instrument. DIAGNOSTIC RECORD DATE AND TIME 04/03/2024 15:41:34
☑ DIAGNOSTIC RECORD DATE AND TIME04/03/2024_15:41:34
☑ DIAGNOSTIC RECORD DATE AND TIME04/03/2024_15:41:34
☑ PROGRAM ☑ FILTER 1 ☑ SAMPLE CHAMBER 48.7°C ☑ FILTER 2 ☑ BREATH TUBE 47.6°C ☑ FILTER 3 ☑ PUMP ☑ INTERNAL STANDARD BREATH ANALYZER ACCURACY STANDARDS ☑ COMPRESSED ETHANOL-GAS MIXTURE
☑ SAMPLE CHAMBER 48.7°C ☑ FILTER 2 ☑ BREATH TUBE 47.6°C ☑ FILTER 3 ☑ PUMP ☑ INTERNAL STANDARD BREATH ANALYZER ACCURACY STANDARDS ☑ COMPRESSED ETHANOL-GAS MIXTURE
☑ BREATH TUBE 47.6°C ☑ FILTER 3 ☑ PUMP ☑ INTERNAL STANDARD BREATH ANALYZER ACCURACY STANDARDS ☑ COMPRESSED ETHANOL-GAS MIXTURE
 ☑ PUMP ☑ INTERNAL STANDARD BREATH ANALYZER ACCURACY STANDARDS ☑ SIMULATOR STANDARD ☑ COMPRESSED ETHANOL-GAS MIXTURE
BREATH ANALYZER ACCURACY STANDARDS SIMULATOR STANDARD COMPRESSED ETHANOL-GAS MIXTURE
☐ SIMULATOR STANDARD ☐ COMPRESSED ETHANOL-GAS MIXTURE
STANDARD SUPPLIER INTOXIMETERS LOT # AG220102 EXP. DATE 07/20/2024
SIMULATOR TEMP (34°C ± 0.2°C) SIM. SN SIM. NIST EXP DATE
 □ CALIBRATION CHECK - (ONLY ONE STANDARD IS TO BE USED PER MAINTENANCE REPORT) Run three tests using a standard. All three tests must be within ±5% of the standard value and must have a spread of .005 or less. Mark the box corresponding to the standard being used. □ 0.10% STANDARD - MUST READ BETWEEN 0.095% AND 0.105% INCLUSIVE □ 0.08% STANDARD - MUST READ BETWEEN 0.076% AND 0.084% INCLUSIVE □ 0.04% STANDARD - MUST READ BETWEEN 0.038% AND 0.042% INCLUSIVE
TEST 1: 0.098 TEST 2: 0.098 TEST 3: 0.097
☑ PERFORM R.F.I. TEST
INDICATE THE NUMBER OF BREATH TESTS IN THE FOLLOWING RANGES SINCE THE LAST MAINTENANCE REPORT:
REFUSALS: 0 004: 0 .0509: 0 .1014: 0 .1519: 0 OVER .19: 0
LIST ANY NEW PARTS AND DESCRIBE ANY ALTERATION OR MODIFICATION THAT WAS MADE TO RESTORE THE INSTRUMENT TO OPERATE SATISFACTORILY AND WITHIN ESTABLISHED LIMITS (USE OTHER SIDE IF NECESSARY)
INSPECTING OFFICER
SIGNATURE PRINT FULL NAME ROY T SMITH
TYPE II PERMIT NUMBER 240040 EXPIRATION DATE 02/05/2026 1TELEPHONE NUMBER 660-385-2132
RETURN COMPLETED REPORT TO THE Breath Alcohol Program, Missouri Department of Health and Senior Services by mail, fax, or email