#### RECEIVED

By Tracy Crews at 7:11 am, Dec 02, 2022



MISSOURI DEPARTMENT OF HEALTH AND SENIOR SERVICES STATE PUBLIC HEALTH LABORATORY BREATH ALCOHOL PROGRAM

## INTOX DMT MAINTENANCE REPORT

REPORT #1

500186	<sup>AME OF AGENCY</sup> <b>Missouri State Highwa</b>	ay Patrol	12/01/2022	
COCATION OF INSTRUMENT (STREET AND CIT 350 E High Street, Jefferson	City		TIME OF INSPECTION 19:25:54	
CHECKLIST: Place a mark in the values where determined). Unmar	box by each item if foun	d to be satisfactory or is o		rite in observed
☑ DIAGNOSTIC RECORD	ned hems must be comed	ned before using instrume	ent.	
DATE AND TIME	22 19:25:58	☑ DETEC	TOR	
☑ PROGRAM		☑ FILTER		
☑ SAMPLE CHAMBER 48	.7°C	☑ FILTER	NO.	
☐ BREATH TUBE 47.2°C		☑ FILTER		No.
☑ PUMP			VAL STANDARD	
BREATH ANALYZER ACCURAC			INC O INNUNIO	
☐ SIMULATOR STANDARD		☑ COMPR	RESSED ETHANOL-GAS MIXTURE	
STANDARD SUPPLIER INTO		LOT# <u>AG125101</u>		
The state of the second of the second				
☐ SIMULATOR TEMP (34°C ± 0.  ☐ CALIBRATION CHECK - (ON Run three tests using a standa of .005 or less. Mark the box of .010% STANDARD - M.	NLY ONE STANDARD I			
☑ CALIBRATION CHECK - (ON Run three tests using a standa	NLY ONE STANDARD IN THE PROPERTY OF STANDARD IN THE PROPERTY OF STANDARD IN THE STANDARD IN TH	S TO BE USED PER M/e within ±5% of the standard being used.  I 0.095% AND 0.105% IN I 0.076% AND 0.084% IN	AINTENANCE REPORT) ard value and must have a spread ICLUSIVE ICLUSIVE	
□ CALIBRATION CHECK - (ON Run three tests using a standa of .005 or less. Mark the box of □ 0.10% STANDARD - N □ 0.08% STANDARD - N □ 0.04% STANDARD - N TEST 1: 0.099	NLY ONE STANDARD IN THE PROPERTY OF STANDARD IN THE PROPERTY OF STANDARD IN THE STANDARD IN TH	S TO BE USED PER M/e within ±5% of the standard being used. I 0.095% AND 0.105% IN I 0.076% AND 0.084% IN I 0.038% AND 0.042% IN	AINTENANCE REPORT) ard value and must have a spread NCLUSIVE NCLUSIVE	
□ CALIBRATION CHECK - (ON Run three tests using a standa of .005 or less. Mark the box of .005 or less. Mark the box of .010% STANDARD - M.     □ 0.08% STANDARD - M.     □ 0.04% STANDARD - M.  TEST 1: 0.099	NLY ONE STANDARD IS Ird. All three tests must be corresponding to the stan MUST READ BETWEEN MUST READ BETWEEN MUST READ BETWEEN TEST 2: 0	S TO BE USED PER M/e within ±5% of the standard being used. N 0.095% AND 0.105% IN 0.076% AND 0.084% IN 0.038% AND 0.042% IN 0.100	AINTENANCE REPORT) ard value and must have a spread ICLUSIVE ICLUSIVE ICLUSIVE TEST 3: 0.099	
CALIBRATION CHECK - (ON Run three tests using a standa of .005 or less. Mark the box or .0.10% STANDARD - M	NLY ONE STANDARD IS Ird. All three tests must be corresponding to the stan MUST READ BETWEEN MUST READ BETWEEN MUST READ BETWEEN TEST 2: 0	S TO BE USED PER M/e within ±5% of the standard being used. N 0.095% AND 0.105% IN 0.076% AND 0.084% IN 0.038% AND 0.042% IN 0.100	AINTENANCE REPORT) ard value and must have a spread ICLUSIVE ICLUSIVE ICLUSIVE TEST 3: 0.099	E REPORT:
CALIBRATION CHECK - (ON Run three tests using a standar of .005 or less. Mark the box of .005 or less. Mark the box of .0.10% STANDARD - №      □ 0.08% STANDARD - №      □ 0.04% STANDARD - №  TEST 1: 0.099      区 PERFORM R.F.I. TEST  INDICATE THE NUMBER OF BR  REFUSALS: 0 004: 0	NLY ONE STANDARD II Ird. All three tests must be corresponding to the stan MUST READ BETWEEN MUST READ BETWEEN TEST 2: 0	S TO BE USED PER M/e within ±5% of the standard being used. N 0.095% AND 0.105% IN 0.076% AND 0.084% IN 0.038% AND 0.042% IN 0.100	AINTENANCE REPORT) ard value and must have a spread  ICLUSIVE ICLUSIVE  TEST 3: 0.099  SINCE THE LAST MAINTENANC	
CALIBRATION CHECK - (ON Run three tests using a standa of .005 or less. Mark the box of .005 or less. Mark the box of .005 or less. Mark the box of .0.10% STANDARD - MODE .0.08% STANDARD - MODE .0.04% OF LIST ANY NEW PARTS AND DESCRIBE ANY ALTER STAND .0.04% OTHER SIDE IF NECESTABLISHED LIMITS (USE OTHER SID	NLY ONE STANDARD II  Ird. All three tests must be corresponding to the stan MUST READ BETWEEN MUST READ BETWEEN TEST 2: 0  REATH TESTS IN THE II  .0509: 6  ERATION OR MODIFICATION THAT ESSARY)	S TO BE USED PER MA e within ±5% of the standard being used. I 0.095% AND 0.105% IN I 0.076% AND 0.084% IN I 0.038% AND 0.042% IN I.100  FOLLOWING RANGES  .1014: 6	AINTENANCE REPORT) ard value and must have a spread  ICLUSIVE ICLUSIVE TEST 3: 0.099  SINCE THE LAST MAINTENANC	
CALIBRATION CHECK - (ON Run three tests using a standar of .005 or less. Mark the box of .005 or less. Mark the box of .0.10% STANDARD - №      □ 0.08% STANDARD - №      □ 0.04% STANDARD - №  TEST 1: 0.099      区 PERFORM R.F.I. TEST  INDICATE THE NUMBER OF BR  REFUSALS: 0 004: 0	NLY ONE STANDARD II  Ird. All three tests must be corresponding to the stan MUST READ BETWEEN MUST READ BETWEEN TEST 2: 0  REATH TESTS IN THE II  .0509: 6  ERATION OR MODIFICATION THAT ESSARY)	S TO BE USED PER MA e within ±5% of the standard being used. I 0.095% AND 0.105% IN I 0.076% AND 0.084% IN I 0.038% AND 0.042% IN I.100  FOLLOWING RANGES  .1014: 6	AINTENANCE REPORT) ard value and must have a spread  ICLUSIVE ICLUSIVE  TEST 3: 0.099  SINCE THE LAST MAINTENANC	
CALIBRATION CHECK - (ON Run three tests using a standa of .005 or less. Mark the box of .005 or less. Mark the box of .005 or less. Mark the box of .0.10% STANDARD - MODEL .0.08% STANDARD - MODEL .0.04% STANDARD - MODEL .0	NLY ONE STANDARD II  Ird. All three tests must be corresponding to the stan MUST READ BETWEEN MUST READ BETWEEN TEST 2: 0  REATH TESTS IN THE II  .0509: 6  ERATION OR MODIFICATION THAT ESSARY)	S TO BE USED PER MA e within ±5% of the standard being used. I 0.095% AND 0.105% IN I 0.076% AND 0.084% IN I 0.038% AND 0.042% IN I.100  FOLLOWING RANGES  .1014: 6	AINTENANCE REPORT) ard value and must have a spread  ICLUSIVE ICLUSIVE  TEST 3: 0.099  SINCE THE LAST MAINTENANC	
CALIBRATION CHECK - (ON Run three tests using a standa of .005 or less. Mark the box of .005 or less. Mark the box of .005 or less. Mark the box of .0.10% STANDARD - MODE .0.08% STANDARD - MODE .0.04% OF LIST ANY NEW PARTS AND DESCRIBE ANY ALTER STAND .0.04% OTHER SIDE IF NECESTABLISHED LIMITS (USE OTHER SID	NLY ONE STANDARD II  Ird. All three tests must be corresponding to the stan MUST READ BETWEEN MUST READ BETWEEN TEST 2: 0  REATH TESTS IN THE II  .0509: 6  ERATION OR MODIFICATION THAT ESSARY)	S TO BE USED PER MA e within ±5% of the standard being used. I 0.095% AND 0.105% IN I 0.076% AND 0.084% IN I 0.038% AND 0.042% IN I.100  FOLLOWING RANGES  .1014: 6	AINTENANCE REPORT) ard value and must have a spread  ICLUSIVE ICLUSIVE  TEST 3: 0.099  SINCE THE LAST MAINTENANC	
CALIBRATION CHECK - (ON Run three tests using a standa of .005 or less. Mark the box of .005 or less. Mark the box of .005 or less. Mark the box of .0.10% STANDARD - MODEL .0.08% STANDARD - MODEL .0.04% STANDARD - MODEL .0	NLY ONE STANDARD II  Ird. All three tests must be corresponding to the stan MUST READ BETWEEN MUST READ BETWEEN TEST 2: 0  REATH TESTS IN THE II  .0509: 6  ERATION OR MODIFICATION THAT ESSARY)	S TO BE USED PER MA e within ±5% of the standard being used. I 0.095% AND 0.105% IN I 0.076% AND 0.084% IN I 0.038% AND 0.042% IN I.100  FOLLOWING RANGES  .1014: 6	AINTENANCE REPORT) ard value and must have a spread  ICLUSIVE ICLUSIVE  TEST 3: 0.099  SINCE THE LAST MAINTENANC	
CALIBRATION CHECK - (ON Run three tests using a standa of .005 or less. Mark the box of .010% STANDARD - M	NLY ONE STANDARD II  Ird. All three tests must be corresponding to the stan MUST READ BETWEEN MUST READ BETWEEN TEST 2: 0  REATH TESTS IN THE II  .0509: 6  ERATION OR MODIFICATION THAT ESSARY)	S TO BE USED PER MA e within ±5% of the standard being used. I 0.095% AND 0.105% IN I 0.076% AND 0.084% IN I 0.038% AND 0.042% IN I.100  FOLLOWING RANGES  .1014: 6	AINTENANCE REPORT) ard value and must have a spread  ICLUSIVE ICLUSIVE  TEST 3: 0.099  SINCE THE LAST MAINTENANC	
□ CALIBRATION CHECK - (ON Run three tests using a standar of .005 or less. Mark the box of	NLY ONE STANDARD II  Ird. All three tests must be corresponding to the stan MUST READ BETWEEN MUST READ BETWEEN TEST 2: 0  REATH TESTS IN THE II  .0509: 6  ERATION OR MODIFICATION THAT ESSARY)	S TO BE USED PER MA e within ±5% of the standard being used. I 0.095% AND 0.105% IN I 0.076% AND 0.084% IN I 0.038% AND 0.042% IN I.100  FOLLOWING RANGES  .1014: 6 T WAS MADE TO RESTORE THE INS	AINTENANCE REPORT) ard value and must have a spread  ICLUSIVE ICLUSIVE  TEST 3: 0.099  SINCE THE LAST MAINTENANC  .1519: 14  STRUMENT TO OPERATE SATISFACTORILY AND W	
CALIBRATION CHECK - (ON Run three tests using a standa of .005 or less. Mark the box of .005 or less. Mark the box of .005 or less. Mark the box of .0.10% STANDARD - MODE .0.08% STANDARD - MODE .0.04% OF LIST ANY NEW PARTS AND DESCRIBE ANY ALTER STAND .0.04% OTHER SIDE IF NECESTABLISHED LIMITS (USE OTHER SID	NLY ONE STANDARD IS ITAL All three tests must be corresponding to the stan MUST READ BETWEEN MUST READ BETWEEN TEST 2: 0  REATH TESTS IN THE I  .0509: 6  ERATION OR MODIFICATION THAT ESSARY)  ARDS. TIME CALIBRATED.	S TO BE USED PER MA e within ±5% of the standard being used. I 0.095% AND 0.105% IN I 0.076% AND 0.084% IN I 0.038% AND 0.042% IN I.100  FOLLOWING RANGES  1014: 6 T WAS MADE TO RESTORE THE INS	AINTENANCE REPORT) ard value and must have a spread  ICLUSIVE ICLUSIVE  TEST 3: 0.099  SINCE THE LAST MAINTENANC  .1519: 14  STRUMENT TO OPERATE SATISFACTORILY AND W	



Airgas USA LLC (LAB) 3500 Bernard Street St. Louis, Mo. 63103 Ph: (314) 533-3100 Fax: (314) 533-7328

## **Certificate of Analysis**

**Customer Name** 

Exclusive Supplier Intoximeters, Inc. 2081 Craig Road St. Louis, Mo 63146 Test Date: 8-Sep-2021

Lot # AG125101 Model 108

**Exp Date** 8-Sep-2023 Cyl. Type 108

Component Ethanol

Certified Concentration 0.100 ± 2% BrAC (260 ppm)

Nitrogen

Certification Traceable to N.I.S.T. RGM and to CRM Ethanol Standards:

RGM Serial No.	Concentration	RGM Serial No.	Concentration
EB0010581	392.1 ppm	EB0010603	393.0 ppm
EB0010570	259.8 ppm	EB0010559	258.2 ppm
EB0010285	208.0 ppm	EB0010595	208.3 ppm
EB0010561	103.6 ppm	EB0010562	104.2 ppm
EB0010681	52.12 ppm	EB0010579	52.81 ppm

CRM Serial No. CC434668

Concentration 800.0 ppm 253.0 ppm

CRM Serial No. 0056649

0056662

Concentration 390.1 ppm 150.2 ppm

Analytical Method: NDIR

CC234503

Digitally signed by Quality Control Reason Dry gas standard certification of analysis Location Airgas USA LLC (Lab) Date:09.09.2021 14:28

Approved for Release:

ISO 17025:2017 A2LA accredited. Certificate Number 3082.06 ISO 17034:2016 A2LA accredited. Certificate Number 3082.07



## STATE OF MISSOURI

DEPARTMENT OF HEALTH AND SENIOR SERVICES
BREATH ALCOHOL PROGRAM



# PERMIT TYPE II

# JACOB M. McKINNEY

is hereby authorized to instruct and supervise operators, train instructors, inspect, calibrate, perform field service and repairs, and operate the following breath analyzer(s):

### **INTOX DMT**

for the determination of the alcoholic content of blood from a sample of expired air. Permit issued under the provisions of sections 577.020 through 577.041. RSMo and 306.111 through 306.119 RSMo.

577.525 tillough 577.541, Holvio and 500.111 tillough 500.115 Holvio	1
DATE 8/12/2021	Laura a Nay
	DIRECTOR OF STATE PUBLIC HEALTH LABORATORY
NUMBER 210178	THE STATE OF THE S
EXPIRES 8/12/2023	DIRECTOR OF DEPARTMENT OF HEALTH AND SENIOR SERVICES

MO 580-0771 (6-10)

LAB-4 (R6-10)



The named cardholder is authorized to operate an evidential breath alcohol instrument for the determination of the alcoholic content in breath form of expired ali in Missouri.

Operator McKINNEY, JACOB Permit No 210178

