

MISSOURI DEPARTMENT OF HEALTH AND SENIOR SERVICES STATE PUBLIC HEALTH LABORATORY BREATH ALCOHOL PROGRAM

	MAINTENANCE REPORT		REPORT #3	
Complete this report at the time of the regular monthly preventive maintenance check (not to exceed 35 days). Complete this report whenever the instrument is serviced or repaired and whenever it is placed				
days). Complete this report whenev	er the instrument is serviced	or repaired and when	ever it is placed	
into service. Retain the original INTOX EC/IR II SN				
12707	NAME OF AGENCY	DATE OF INSPI		
LOCATION OF INSTRUMENT (STREET AND CLTS	Blue Springs Police Dept			
1100 SW. Smith St. Blue Spring		TIME OF INSPI	ECTION	
		14:25 CST		
CHECKLIST: Place a mark in the box established limits. (Write in obse	by each item if found to be	satisfactory or is op	erating within	
before using instrument.	ived values where determined)	. Unmarked Items mus	t be corrected	
X DIAGNOSTIC RECORD				
X BLANK CHECK	X CO2 C	HECK		
X FC 1 TEMP				
The state of the s				
X SRC TEMP X FCB CHECK				
X DET TEMP	X CRC C	OMP CHECK		
X BT TEMP X CRC CAL CHECK				
X STD 2 TEMP X PRINT TEST				
X ETH CHECK				
BREATH ANALYZER ACCURACY STAND	ARDS			
SIMULATOR SOLUTION		EGGED PERAMOT CAS A	AT VIII ID D	
A WASHINGTON OF THE PARTY OF TH				
SIMULATOR TEMP (34°C +0.2°C)			EXP. DATE 03/18/2026	
LISTMODIATOR TEMP (54 C ±0.2 C)	SIM. SN	SIM. NIST	EXP DATE	
X CALIBRATION CHECK - (ONLY ON	E STANDARD IS TO BE USED P	ER MAINTENANCE REPO	ORT)	
Run three tests using a stan	dard solution. All three	tests must be withi	n +5% of the standard value	
Run three tests using a stan and must have a spread of .0	dard solution. All three	tests must be withi	n +5% of the standard value	
Run three tests using a stan and must have a spread of .0 used.	dard solution. All three 05 or less. Mark the box	tests must be withicorresponding to th	n +5% of the standard value	
Run three tests using a stan and must have a spread of .0 used. X 0.10% STANDARD - MUST READ	dard solution. All three 05 or less. Mark the box BETWEEN 0.095% AND 0.105%	tests must be withicorresponding to the	n +5% of the standard value	
Run three tests using a stan and must have a spread of .0 used. X 0.10% STANDARD - MUST READ 0.08% STANDARD - MUST READ	dard solution. All three 05 or less. Mark the box BETWEEN 0.095% AND 0.105% BETWEEN 0.076% AND 0.084%	tests must be withicorresponding to the INCLUSIVE INCLUSIVE	n +5% of the standard value	
Run three tests using a stan and must have a spread of .0 used. X 0.10% STANDARD - MUST READ	dard solution. All three 05 or less. Mark the box BETWEEN 0.095% AND 0.105% BETWEEN 0.076% AND 0.084%	tests must be withicorresponding to the INCLUSIVE INCLUSIVE	n +5% of the standard value	
Run three tests using a stan and must have a spread of .0 used. X 0.10% STANDARD - MUST READ 0.08% STANDARD - MUST READ 0.04% STANDARD - MUST READ	dard solution. All three 05 or less. Mark the box BETWEEN 0.095% AND 0.105% BETWEEN 0.076% AND 0.084% BETWEEN 0.038% AND 0.042%	tests must be within corresponding to the INCLUSIVE INCLUSIVE INCLUSIVE	In $\pm 5\%$ of the standard value are standard solution being	
Run three tests using a stan and must have a spread of .0 used. X 0.10% STANDARD - MUST READ 0.08% STANDARD - MUST READ 0.04% STANDARD - MUST READ TEST 1 5 0.100 g/210L	dard solution. All three 05 or less. Mark the box BETWEEN 0.095% AND 0.105% BETWEEN 0.076% AND 0.084% BETWEEN 0.038% AND 0.042% TEST 2 P 0.100 g/210L	tests must be within corresponding to the INCLUSIVE INCLUSIVE INCLUSIVE	n ±5% of the standard value ne standard solution being	
Run three tests using a stan and must have a spread of .0 used. X 0.10% STANDARD - MUST READ 0.08% STANDARD - MUST READ 0.04% STANDARD - MUST READ	dard solution. All three 05 or less. Mark the box BETWEEN 0.095% AND 0.105% BETWEEN 0.076% AND 0.084% BETWEEN 0.038% AND 0.042% TEST 2 P 0.100 g/210L	tests must be within corresponding to the INCLUSIVE INCLUSIVE INCLUSIVE	n ±5% of the standard value ne standard solution being	
Run three tests using a stan and must have a spread of .0 used. X 0.10% STANDARD - MUST READ 0.08% STANDARD - MUST READ 0.04% STANDARD - MUST READ TEST 1 © 0.100 g/210L INDICATE THE NUMBER OF BREATH	dard solution. All three 05 or less. Mark the box BETWEEN 0.095% AND 0.105% BETWEEN 0.076% AND 0.084% BETWEEN 0.038% AND 0.042% TEST 2 12 0.100 g/210L TESTS IN THE FOLLOWING RAN	tests must be within corresponding to the INCLUSIVE INCLUSIVE INCLUSIVE TEST 3 FEST SINCE THE LAST	n ±5% of the standard value me standard solution being 0.100 g/210L MAINTENANCE REPORT:	
Run three tests using a stan and must have a spread of .0 used. X 0.10% STANDARD - MUST READ 0.08% STANDARD - MUST READ 0.04% STANDARD - MUST READ TEST 1 © 0.100 g/210L INDICATE THE NUMBER OF BREATH	dard solution. All three 05 or less. Mark the box BETWEEN 0.095% AND 0.105% BETWEEN 0.076% AND 0.084% BETWEEN 0.038% AND 0.042% TEST 2 *** 0.100 g/210L TESTS IN THE FOLLOWING RAN	tests must be within corresponding to the INCLUSIVE INCLUSIVE INCLUSIVE TEST 3 F GES SINCE THE LAST	n ±5% of the standard value the standard solution being 0.100 g/210L MAINTENANCE REPORT: 0 OVER .19 0	
Run three tests using a stan and must have a spread of .0 used. X 0.10% STANDARD - MUST READ 0.08% STANDARD - MUST READ 0.04% STANDARD - MUST READ TEST 1 2 0.100 g/210L INDICATE THE NUMBER OF BREATH STANDARD 1 004 0 LIET ANY NEW PARTS AND DESCRIBE ANY ALT	dard solution. All three 05 or less. Mark the box BETWEEN 0.095% AND 0.105% BETWEEN 0.076% AND 0.084% BETWEEN 0.038% AND 0.042% TEST 2 18 0.100 g/210L TESTS IN THE FOLLOWING RAN .0509 0 .1014	tests must be within corresponding to the INCLUSIVE INCLUSIVE INCLUSIVE TEST 3 F GES SINCE THE LAST 0 .1519	n ±5% of the standard value the standard solution being 0.100 g/210L MAINTENANCE REPORT: 0 OVER .19 0	
Run three tests using a stan and must have a spread of .0 used. X 0.10% STANDARD - MUST READ 0.08% STANDARD - MUST READ 0.04% STANDARD - MUST READ TEST 1 © 0.100 g/210L INDICATE THE NUMBER OF BREATH	dard solution. All three 05 or less. Mark the box BETWEEN 0.095% AND 0.105% BETWEEN 0.076% AND 0.084% BETWEEN 0.038% AND 0.042% TEST 2 18 0.100 g/210L TESTS IN THE FOLLOWING RAN .0509 0 .1014	tests must be within corresponding to the INCLUSIVE INCLUSIVE INCLUSIVE TEST 3 F GES SINCE THE LAST 0 .1519	n ±5% of the standard value the standard solution being 0.100 g/210L MAINTENANCE REPORT: 0 OVER .19 0	
Run three tests using a stan and must have a spread of .0 used. X 0.10% STANDARD - MUST READ 0.08% STANDARD - MUST READ 0.04% STANDARD - MUST READ TEST 1 2 0.100 g/210L INDICATE THE NUMBER OF BREATH STANDARD 1 004 0 LIET ANY NEW PARTS AND DESCRIBE ANY ALT	dard solution. All three 05 or less. Mark the box BETWEEN 0.095% AND 0.105% BETWEEN 0.076% AND 0.084% BETWEEN 0.038% AND 0.042% TEST 2 18 0.100 g/210L TESTS IN THE FOLLOWING RAN .0509 0 .1014	tests must be within corresponding to the INCLUSIVE INCLUSIVE INCLUSIVE TEST 3 F GES SINCE THE LAST 0 .1519	n ±5% of the standard value the standard solution being 0.100 g/210L MAINTENANCE REPORT: 0 OVER .19 0	
Run three tests using a stan and must have a spread of .0 used. X 0.10% STANDARD - MUST READ 0.08% STANDARD - MUST READ 0.04% STANDARD - MUST READ TEST 1 . 0.100 g/210L INDICATE THE NUMBER OF BREATH SEFUSALS 1 004 0 LIET ANY NEW PARTS AND DESCRIBE ANY ALT SATISFACTORILY AND WITHIN ESTABLISHED L	dard solution. All three 05 or less. Mark the box BETWEEN 0.095% AND 0.105% BETWEEN 0.076% AND 0.084% BETWEEN 0.038% AND 0.042% TEST 2 18 0.100 g/210L TESTS IN THE FOLLOWING RAN .0509 0 .1014	tests must be within corresponding to the INCLUSIVE INCLUSIVE INCLUSIVE TEST 3 F GES SINCE THE LAST 0 .1519	n ±5% of the standard value the standard solution being 0.100 g/210L MAINTENANCE REPORT: 0 OVER .19 0	
Run three tests using a stan and must have a spread of .0 used. X 0.10% STANDARD - MUST READ 0.08% STANDARD - MUST READ 0.04% STANDARD - MUST READ TEST 1 . 0.100 g/210L INDICATE THE NUMBER OF BREATH SEFUSALS 1 004 0 LIET ANY NEW PARTS AND DESCRIBE ANY ALT SATISFACTORILY AND WITHIN ESTABLISHED L	dard solution. All three 05 or less. Mark the box BETWEEN 0.095% AND 0.105% BETWEEN 0.076% AND 0.084% BETWEEN 0.038% AND 0.042% TEST 2 18 0.100 g/210L TESTS IN THE FOLLOWING RAN .0509 0 .1014	tests must be within corresponding to the INCLUSIVE INCLUSIVE INCLUSIVE TEST 3 F GES SINCE THE LAST 0 .1519	n ±5% of the standard value the standard solution being 0.100 g/210L MAINTENANCE REPORT: 0 OVER .19 0	
Run three tests using a stan and must have a spread of .0 used. X 0.10% STANDARD - MUST READ 0.08% STANDARD - MUST READ 0.04% STANDARD - MUST READ TEST 1 2 0.100 g/210L INDICATE THE NUMBER OF BREATH STANDARD 1 004 0 LIST ANY NEW PARTS AND DESCRIBE ANY ALT SATISFACTORILY AND WITHIN ESTABLISHED L NOV 2024 MAINTENANCE	dard solution. All three 05 or less. Mark the box BETWEEN 0.095% AND 0.105% BETWEEN 0.076% AND 0.084% BETWEEN 0.038% AND 0.042% TEST 2 18 0.100 g/210L TESTS IN THE FOLLOWING RAN .0509 0 .1014	tests must be within corresponding to the INCLUSIVE INCLUSIVE INCLUSIVE TEST 3 F GES SINCE THE LAST 0 .1519	n ±5% of the standard value the standard solution being 0.100 g/210L MAINTENANCE REPORT: 0 OVER .19 0	
Run three tests using a stan and must have a spread of .0 used. X 0.10% STANDARD - MUST READ 0.08% STANDARD - MUST READ 0.04% STANDARD - MUST READ TEST 1 2 0.100 g/210L INDICATE THE NUMBER OF BREATH SATISFACTORILY AND WITHIN ESTABLISHED L NOV 2024 MAINTENANCE	dard solution. All three 05 or less. Mark the box BETWEEN 0.095% AND 0.105% BETWEEN 0.076% AND 0.084% BETWEEN 0.038% AND 0.042% TEST 2 10 0.100 g/210L TESTS IN THE FOLLOWING RAN .0509 0 .1014 BRATION OR MODIFICATION THAT WAS IMITS (USE OTHER SIDE IF NECESSAR	tests must be within corresponding to the INCLUSIVE INCLUSIVE INCLUSIVE TEST 3 F GES SINCE THE LAST 0 .1519	n ±5% of the standard value the standard solution being 0.100 g/210L MAINTENANCE REPORT: 0 OVER .19 0	
Run three tests using a stan and must have a spread of .0 used. X 0.10% STANDARD - MUST READ 0.08% STANDARD - MUST READ 0.04% STANDARD - MUST READ TEST 1 2 0.100 g/210L INDICATE THE NUMBER OF BREATH SATISFACTORILY AND WITHIN ESTABLISHED L NOV 2024 MAINTENANCE	dard solution. All three 05 or less. Mark the box BETWEEN 0.095% AND 0.105% BETWEEN 0.076% AND 0.084% BETWEEN 0.038% AND 0.042% TEST 2 10 0.100 g/210L TESTS IN THE FOLLOWING RAN .0509 0 .1014 ERATION OR MODIFICATION THAT WAS IMITS (USE OTHER SIDE IF NECESSAR	tests must be within corresponding to the INCLUSIVE INCLUSIVE INCLUSIVE TEST 3 F GES SINCE THE LAST O .1519 MADE TO RESTORE THE INST ().	n ±5% of the standard value the standard solution being 0.100 g/210L MAINTENANCE REPORT: 0 OVER .19 0	
Run three tests using a stan and must have a spread of .0 used. X 0.10% STANDARD - MUST READ 0.08% STANDARD - MUST READ 0.04% STANDARD - MUST READ TEST 1 . 0.100 g/210L INDICATE THE NUMBER OF BREATH STANDARD 1 004 0 LIET ANY NEW PARTS AND DESCRIBE ANY ALT SATISFACTORILY AND WITHIN ESTABLISHED L NOV 2024 MAINTENANCE INSERCTING OFFICER SIGNATURE	dard solution. All three 05 or less. Mark the box BETWEEN 0.095% AND 0.105% BETWEEN 0.076% AND 0.084% BETWEEN 0.038% AND 0.042% TEST 2 10 0.100 g/210L TESTS IN THE FOLLOWING RAN .0509 0 .1014 BRATION OR MODIFICATION THAT WAS IMITS (USE OTHER SIDE IF NECESSAR	tests must be within corresponding to the inclusive incl	n ±5% of the standard value the standard solution being 0.100 g/210L MAINTENANCE REPORT: 0 OVER .19 0	
Run three tests using a stan and must have a spread of .0 used. X 0.10% STANDARD - MUST READ 0.08% STANDARD - MUST READ 0.04% STANDARD - MUST READ TEST 1 2 0.100 g/210L INDICATE THE NUMBER OF BREATH STANDARD 1 004 0 LIET ANY NEW PARTS AND DESCRIBE ANY ALT SATISFACTORILY AND WITHIN ESTABLISHED L NOV 2024 MAINTENANCE INSPECTING OFFICER SIGNATURE TYPE IT FERMIT NUMBER EXPIRE	dard solution. All three 05 or less. Mark the box BETWEEN 0.095% AND 0.105% BETWEEN 0.076% AND 0.084% BETWEEN 0.038% AND 0.042% TEST 2 10 0.100 g/210L TESTS IN THE FOLLOWING RAN .0509 0 .1014 ERATION OR MODIFICATION THAT WAS IMITS (USE OTHER SIDE IF NECESSAR PRINT FOLLOWING CONTROL OF C. L. TELEPHONE	tests must be within corresponding to the inclusive incl	n ±5% of the standard value the standard solution being 0.100 g/210L MAINTENANCE REPORT: 0 OVER .19 0	
Run three tests using a stan and must have a spread of .0 used. X 0.10% STANDARD - MUST READ 0.08% STANDARD - MUST READ 0.04% STANDARD - MUST READ 1.00 G/210L TEST 1 © 0.100 G/210L INDICATE THE NUMBER OF BREATH 1.00 OF BREATH 1.	dard solution. All three 05 or less. Mark the box BETWEEN 0.095% AND 0.105% BETWEEN 0.076% AND 0.084% BETWEEN 0.038% AND 0.042% TEST 2 10 0.100 g/210L TESTS IN THE FOLLOWING RAN .0509 0 .1014 ERATION OR MODIFICATION THAT WAS IMITS (USE OTHER SIDE IF NECESSAR PRINT BULL OFC. L ATTON DATE TELEPHONE (816)	tests must be within corresponding to the inclusive incl	n ±5% of the standard value the standard solution being 0.100 g/210L MAINTENANCE REPORT: 0 OVER .19 0	
Run three tests using a stan and must have a spread of .0 used. X 0.10% STANDARD - MUST READ 0.08% STANDARD - MUST READ 0.04% STANDARD - MUST READ 1.000 G/210L INDICATE THE NUMBER OF BREATH 1.000 G/210L REFUSALS 1 004 0 LIET ANY NEW PARTS AND DESCRIBE ANY ALT SATISFACTORILY AND WITHIN ESTABLISHED L. NOV 2024 MAINTENANCE INSPECTING OFFICER SIGNATURE TYPE IT FERMIT RUMBER 230155 08/0	dard solution. All three 05 or less. Mark the box BETWEEN 0.095% AND 0.105% BETWEEN 0.076% AND 0.084% BETWEEN 0.038% AND 0.042% TEST 2 ** 0.100 g/210L TESTS IN THE FOLLOWING RAN .0509 0 .1014 ERATION OF MUDIFICATION THAT WAS IMITS (USE OTHER SIDE IF NECESSAR PRINT FULL OFC. L ATION DATE TELEPHONE 1/2025 (816)	tests must be within corresponding to the INCLUSIVE INCLUSIVE INCLUSIVE TEST 3 F GES SINCE THE LAST O .1519 MADE TO RESTORE THE INST ().	n ±5% of the standard value are standard solution being 0.100 g/210L MAINTENANCE REPORT: 0 OVER .19 0 RUMENT TO OPERATE	
Run three tests using a stan and must have a spread of .0 used. X 0.10% STANDARD - MUST READ 0.08% STANDARD - MUST READ 0.04% STANDARD - MUST READ 1.000 G/210L INDICATE THE NUMBER OF BREATH 1.000 G/210L REFUSALS 1 004 0 LIET ANY NEW PARTS AND DESCRIBE ANY ALT SATISFACTORILY AND WITHIN ESTABLISHED L. NOV 2024 MAINTENANCE INSPECTING OFFICER SIGNATURE TYPE IT FERMIT RUMBER 230155 08/0	dard solution. All three 05 or less. Mark the box BETWEEN 0.095% AND 0.105% BETWEEN 0.076% AND 0.084% BETWEEN 0.038% AND 0.042% TEST 2 ** 0.100 g/210L TESTS IN THE FOLLOWING RAN .0509 0 .1014 ERATION OF MUDIFICATION THAT WAS IMITS (USE OTHER SIDE IF NECESSAR PRINT FULL OFC. L ATION DATE TELEPHONE 1/2025 (816)	tests must be within corresponding to the INCLUSIVE INCLUSIVE INCLUSIVE TEST 3 F GES SINCE THE LAST O .1519 MADE TO RESTORE THE INST ().	n ±5% of the standard value are standard solution being 0.100 g/210L MAINTENANCE REPORT: 0 OVER .19 0 RUMENT TO OPERATE	
Run three tests using a stan and must have a spread of .0 used. X 0.10% STANDARD - MUST READ 0.08% STANDARD - MUST READ 0.04% STANDARD - MUST READ 1.00% STANDARD ST	dard solution. All three 05 or less. Mark the box BETWEEN 0.095% AND 0.105% BETWEEN 0.076% AND 0.084% BETWEEN 0.038% AND 0.042% TEST 2 ** 0.100 g/210L TESTS IN THE FOLLOWING RAN .0509 0 .1014 ERATION OF MUDIFICATION THAT WAS IMITS (USE OTHER SIDE IF NECESSAR PRINT FULL OFC. L ATION DATE TELEPHONE 1/2025 (816)	tests must be within corresponding to the INCLUSIVE INCLUSIVE INCLUSIVE TEST 3 F GES SINCE THE LAST O .1519 MADE TO RESTORE THE INST ().	n ±5% of the standard value are standard solution being 0.100 g/210L MAINTENANCE REPORT: 0 OVER .19 0 RUMENT TO OPERATE	



Airgas USA LLC (LAB) 3500 Bernard Street St. Louis, Mo. 63103 Ph: (314) 533-3100 Fax: (314) 533-7328

Certificate of Analysis

Customer Name Exclusive Supplier Intoximeters, Inc. 2081 Craig Road St. Louis, Mo 63146

Test Date: 19-Mar-2024

Lot # AG407801 **Model** 108

Exp Date

Cyl. Type

Component

Certified Concentration

18-Mar-2026

108

Ethanol Nitrogen 0.100 ± 2% BrAC (272 ppm)

Certification Traceable to N.I.S.T. RGM and to CRM Ethanol Standards:

RGM Serial No. Concentration EB0010581 391.8 ppm EB0010570 259.8 ppm EB0010285 209.0 ppm EB0010561 103.7 ppm EB0010681 52.22 ppm

RGM Serial No. Concentration EB0010603 392.5 ppm EB0010559 258.9 ppm EB0010562 104.2 ppm EB0010579 52.94 ppm

CRM Serial No. CC727481 CC727496

Concentration 799.4 ppm 253.4 ppm

CRM Serial No. CC727493 CC727498

Concentration 389.8 ppm 150.2 ppm

Analytical Method: NDIR

Digitally signed by Quality Control Reason:Dry gas standard certification of analysis Location:Argas USA LLC (Lab) Date:03.22.2024 07:50

Approved for Release:

Yusef Woods

ISO 17025:2017 A2LA accredited. Certificate Number 3082.06 ISO 17034:2016 A2LA accredited. Certificate Number 3082.07



STATE OF MISSOURI

DEPARTMENT OF HEALTH AND SENIOR SERVICES **BREATH ALCOHOL PROGRAM**



PERMIT TYPE II JORDAN LITZ

is hereby authorized to instruct and supervise operators, train instructors, inspect, calibrate, perform field service and repairs, and operate the following breath analyzer(s):

INTOX EC/IR II

for the determination of the alcoholic content of blood from a sample of expired air. Permit issued under the provisions of sections 577.020 through 577.041, RSMo and 306.111 through 306.119 RSMo.

DATE 8/1/2023	Mike Masson
5/112	DIRECTOR OF STATE PUBLIC HEALTH LABORATORY
NUMBER 230155	
	Davla I. Nichelson
EXPIRES 8/1/2025	- Company
	DIRECTOR OF DEPARTMENT OF HEALTH AND SENIOR SERVICES
MO 580-0771 (6-10)	LAR 4 (DE 10

LAB-4 (R6-10)



The named cardholder is authorized to operate an evidential breath alcohol instrument for the determination of the alcoholic content in breath form of expired ail in Missouri.

Operator LITZ, JORDAN

Permit No 230155 Date Issued 8/1/2023 Date Expires 8/1/2025

