

MISSOURI DEPARTMENT OF HEALTH AND SENIOR SERVICES STATE PUBLIC HEALTH LABORATORY

BREATH ALCOHOL PROGRAM

INTOX EC/IR II MAINTENANCE REPORT

REPORT #3

Complete this report at the time of the regular monthly preventive maintenance check (not to exceed 35 days). Complete this report whenever the instrument is serviced or repaired and whenever it is placed into service. Retain the original and send a copy within 15 days to the Breath Alcohol Program, DHSS.

INTOX EC/IR II SN	NAME OF AGENCY	DATE OF INSPECTION	
12818	ST. ANN PD	12/08/2024	
LOCATION OF INSTRUMENT (STREET AND CIT	()	TIME OF INSPECTION	
10405 St. Charles Rock R St. A	17:45 CST		
CHECKLIST: Place a mark in the box	by each item if found to be satisfact	tory or is operating within	
	erved values where determined). Unmark	ked items must be corrected	
before using instrument.			
X DIAGNOSTIC RECORD	The second of the second s		

X BLANK CHECK	X BLANK CHECK			X CO2 CHECK						
XFC 1 TEMP			X FLOW CHECK							
X SRC TEMP			X FCB CHEC	K						
X DET TEMP			X CRC COMP CHECK							
X BT TEMP			X CRC CAL	CHECK						
X STD 2 TEMP			X PRINT TE	ST						
X ETH CHECK	and the states of the		have the second			al de lasse lande				
REATH ANALYZER A	CCURACY ST	ANDARDS			1 States				and a second second	
SIMULATOR SOL	UTION	一、调查算机,非常自由的		X COMPRESS	ED ET	HANOL-GAS	MIXTU	JRE		
STANDARD SUPPLI	ER IN	TOXIMETERS		LOT# AG309	502	Stand Station .	EXP.	DATE 04/	05/2	025
SIMULATOR TEMP			SIM.			SIM. NIS		DIIID		
Run three tests and must have a used. X 0.10% STANDAR 0.08% STANDAR	using a s spread of D - MUST R D - MUST R	tandard solut .005 or less EAD BETWEEN 0 EAD BETWEEN 0	ion. A Mark . 095% A .076% A	the box cor: ND 0.105% IN ND 0.084% IN	ts mus respon CLUSI	st be with nding to VE VE	hin <u>+</u> 5			
Run three tests and must have a used. X 0.10% STANDAR 0.08% STANDAR 0.04% STANDAR	using a s spread of D - MUST R D - MUST R D - MUST R	tandard solut .005 or less EAD BETWEEN 0 EAD BETWEEN 0	ion. A Mark .095% A .076% A .038% A	11 three test the box cor: ND 0.105% IN ND 0.084% IN ND 0.042% IN	ts mus respon CLUSI	st be with nding to VE VE VE	hin <u>+</u> 5 the st			
Run three tests and must have a used. X 0.10% STANDAR 0.08% STANDAR 0.04% STANDAR EST 1 0.100 g/	using a s spread of D - MUST R D - MUST R D - MUST R	tandard solut .005 or less EAD BETWEEN 0 EAD BETWEEN 0 TEST 2	ion. A Mark .095% A .076% A .038% A	ll three test the box cor: ND 0.105% IN ND 0.084% IN ND 0.042% IN g/210L	ts musicular CLUSI CLUSI CLUSI	st be with nding to VE VE VE TEST 3	hin <u>+</u> 5 the st	andard sol	utio	on being
Run three tests and must have a used. 0.10% STANDAR 0.08% STANDAR 0.04% STANDAR EST 1 0.100 g/	using a s spread of D - MUST R D - MUST R D - MUST R /210L ER OF BREA	tandard solut .005 or less EAD BETWEEN 0 EAD BETWEEN 0 TEST 2	ion. A Mark .095% A .076% A .038% A .038% A	ll three test the box cor: ND 0.105% IN ND 0.084% IN ND 0.042% IN g/210L	ts musicus CLUSI CLUSI CLUSI SINC	st be with nding to VE VE VE TEST 3	hin <u>+</u> 5 the st	andard sol	SPOR	on being
Run three tests and must have a used. 0.10% STANDAR 0.08% STANDAR 0.04% STANDAR EST 1 0.100 g/ NDICATE THE NUMB EFUSALS 1	using a s spread of D - MUST R D - MUST R D - MUST R /210L ER OF BREA 004 0 DESCRIBE ANY	tandard solut .005 or less EAD BETWEEN 0 EAD BETWEEN 0 EAD BETWEEN 0 TEST 2 4 TH TESTS IN T .0509 ALTERATION OR MO	ion. A Mark . 095% A . 076% A . 038% A . 038% A . 0.100 HE FOLLO	ll three test the box cor: ND 0.105% IN ND 0.084% IN ND 0.042% IN g/210L owing ranges .1014	ts musicus CLUSI CLUSI CLUSI SINC	st be with ding to VE VE VE TEST 3 E THE LAS .1519	hin <u>+</u> 5 the st 0.10 T MAIN 0	OVER .1	SPOR	on being
Run three tests and must have a used. 0.10% STANDAR 0.08% STANDAR 0.04% STANDAR EST 1 0.100 g/ NDICATE THE NUMB EFUSALS 1	using a s spread of D - MUST R D - MUST R D - MUST R /210L ER OF BREA 004 0 DESCRIBE ANY	tandard solut .005 or less EAD BETWEEN 0 EAD BETWEEN 0 EAD BETWEEN 0 TEST 2 4 TH TESTS IN T .0509 ALTERATION OR MO	ion. A Mark . 095% A . 076% A . 038% A . 038% A . 0.100 HE FOLLO	ll three test the box cor: ND 0.105% IN ND 0.084% IN ND 0.042% IN g/210L owing ranges .1014	ts musicus CLUSI CLUSI CLUSI SINC	st be with ding to VE VE VE TEST 3 E THE LAS .1519	hin <u>+</u> 5 the st 0.10 T MAIN 0	OVER .1	SPOR	on being
Run three tests and must have a used. 0.10% STANDAR 0.08% STANDAR 0.04% STANDAR EST 1 0.100 g/ NDICATE THE NUMB EFUSALS 1 ST ANY NEW PARTS AND TISFACTORILY AND WITH	using a s spread of D - MUST R D - MUST R D - MUST R /210L ER OF BREA 004 0 DESCRIBE ANY HIN ESTABLISH	tandard solut .005 or less EAD BETWEEN 0 EAD BETWEEN 0 EAD BETWEEN 0 TEST 2 4 TH TESTS IN T .0509 ALTERATION OR MO	ion. A Mark . 095% A . 076% A . 038% A . 038% A . 0.100 HE FOLLO	ll three test the box cor: ND 0.105% IN ND 0.084% IN ND 0.042% IN g/210L owing ranges .1014	ts musicus CLUSI CLUSI CLUSI SINC	st be with ding to VE VE VE TEST 3 E THE LAS .1519	hin <u>+</u> 5 the st 0.10 T MAIN 0	OVER .1	SPOR	on being
Run three tests and must have a used. 0.10% STANDAR 0.08% STANDAR 0.04% STANDAR EST 1 0.100 g/ NDICATE THE NUMB EFUSALS 1 ST ANY NEW PARTS AND TISFACTORILY AND WITH	using a s spread of D - MUST R D - MUST R D - MUST R /210L ER OF BREA 004 0 DESCRIBE ANY HIN ESTABLISH	tandard solut .005 or less EAD BETWEEN 0 EAD BETWEEN 0 EAD BETWEEN 0 TEST 2 4 TH TESTS IN T .0509 ALTERATION OR MO	ion. A Mark . 095% A . 076% A . 038% A . 038% A . 0.100 HE FOLLO	ll three test the box cor: ND 0.105% IN ND 0.084% IN ND 0.042% IN g/210L owing ranges .1014	ts musicus CLUSI CLUSI CLUSI SINC	st be with ding to VE VE VE TEST 3 E THE LAS .1519	hin <u>+</u> 5 the st 0.10 T MAIN 0	OVER .1	SPOR	on being
Run three tests and must have a used. 0.10% STANDAR 0.08% STANDAR 0.04% STANDAR EST 1 0.100 g/ NDICATE THE NUMB EFUSALS 1 ST ANY NEW PARTS AND ATISFACTORILY AND WITH REALS 1 ST ANY NEW PARTS AND ATISFACTORILY AND WITH	using a s spread of D - MUST R D - MUST R D - MUST R /210L ER OF BREA 004 0 DESCRIBE ANY HIN ESTABLISH	tandard solut .005 or less EAD BETWEEN 0 EAD BETWEEN 0 EAD BETWEEN 0 TEST 2 4 TH TESTS IN T .0509 ALTERATION OR MO	ion. A Mark . 095% A . 076% A . 038% A . 038% A . 0.100 HE FOLLO	ll three test the box cor: ND 0.105% IN ND 0.084% IN ND 0.042% IN g/210L owing ranges .1014	ts musicus CLUSI CLUSI CLUSI SINC	st be with ding to VE VE VE TEST 3 E THE LAS .1519	hin <u>+</u> 5 the st 0.10 T MAIN 0	OVER .1	SPOR	on being
Run three tests and must have a used. 0.10% STANDAR 0.08% STANDAR 0.04% STANDAR EST 1 0.100 g/ NDICATE THE NUMB EFUSALS 1 ST ANY NEW PARTS AND ATISFACTORILY AND WITH NSPECTING OFFICER	using a s spread of D - MUST R D - MUST R D - MUST R /210L ER OF BREA 004 0 DESCRIBE ANY HIN ESTABLISH	tandard solut .005 or less EAD BETWEEN 0 EAD BETWEEN 0 EAD BETWEEN 0 TEST 2 4 TH TESTS IN T .0509 ALTERATION OR MO	ion. A Mark . 095% A . 076% A . 038% A . 038% A . 0.100 HE FOLLO	<pre>11 three test the box cor: ND 0.105% ING ND 0.084% ING ND 0.042% ING g/210L OWING RANGES .1014 N THAT WAS MADE IF NECESSARY).</pre>	ts musicespoint CLUSI CLUSI CLUSI SINC	st be with ding to VE VE VE TEST 3 E THE LAS .1519	hin <u>+</u> 5 the st 0.10 T MAIN 0	OVER .1	SPOR	on being
Run three tests and must have a used. 0.10% STANDAR 0.08% STANDAR 0.04% STANDAR EST 1 0.100 g/ NDICATE THE NUMB EFUSALS 1 ST ANY NEW PARTS AND ATISFACTORILY AND WITH ST ANY NEW PARTS AND ATISFACTORILY AND WITH INSPECTING OFFICER	using a s spread of D - MUST R D - MUST R D - MUST R /210L ER OF BREA 004 0 DESCRIBE ANY HIN ESTABLISH	tandard solut .005 or less EAD BETWEEN 0 EAD BETWEEN 0 EAD BETWEEN 0 TEST 2 4 TH TESTS IN T .0509 ALTERATION OR MO	ion. A Mark . 095% A . 076% A . 038% A . 038% A . 0.100 HE FOLLO	<pre>11 three test the box cor: ND 0.105% ING ND 0.084% ING ND 0.042% ING g/210L OWING RANGES .1014 N THAT WAS MADE IF NECESSARY).</pre>	ts musicespoint CLUSI CLUSI CLUSI SINC	st be with ding to VE VE VE TEST 3 .1519 TORE THE IN	hin <u>+</u> 5 the st 0.10 T MAIN 0	OVER .1	SPOR	on being
Run three tests and must have a used. 0.10% STANDAR 0.08% STANDAR 0.04% STANDAR EST 1 0.100 g/ NDICATE THE NUMB EFUSALS 1 IST ANY NEW PARTS AND ATISFACTORILY AND WITH	LUSING A S Spread of D - MUST R D - MUST R D - MUST R /210L ER OF BREA 004 0 DESCRIBE ANY HIN ESTABLISH	tandard solut .005 or less EAD BETWEEN 0 EAD BETWEEN 0 EAD BETWEEN 0 TEST 2 4 TH TESTS IN T .0509 ALTERATION OR MO	ion. A Mark . 095% A . 076% A . 038% A . 038% A . 0.100 HE FOLLO	<pre>11 three test the box cor: ND 0.105% ING ND 0.084% ING ND 0.042% ING g/210L OWING RANGES .1014 N THAT WAS MADE IF NECESSARY).</pre>	ts musicespoint CLUSI CLUSI CLUSI SINCE 0 TO RES	st be with ding to VE VE VE TEST 3 .1519 TORE THE IN	hin <u>+</u> 5 the st 0.10 T MAIN 0	OVER .1	SPOR	on being

Breath Alcohol Program, Missouri Department of Health and Senior Services, by mail, fax, or e-mail LAB 163 AN EQUAL OPPORTUNITY/AFFIRMATIVE ACTION EMPLOYER MO 580-2899(5-19) services provided on a nondiscriminatory basis



Airgas USA LLC (LAB) 3500 Bernard Street St. Louis, Mo. 63103 Ph: (314) 533-3100 Fax: (314) 533-7328

Certificate of Analysis

Customer Name Exclusive Supplier Intoximeters, Inc. 2081 Craig Road St. Louis, Mo 63146

pre-

Test Date: 5-Apr-2023

Lot # AG309502 Model 108

Exp DateCyl. TypeComponentCertifie5-Apr-2025108Ethanol0.100 ±Nitrogen

Certified Concentration 0.100 ± 2% BrAC (272 ppm)

Certification Traceable to N.I.S.T. RGM and to CRM Ethanol Standards:

RGM Serial No.	Concentration	RGM Serial No.	Concentration
EB0010581	391.8 ppm	EB0010603	and the second
EB0010570	259.8 ppm	EB0010559	392.5 ppm
EB0010285	209.0 ppm	EB0010562	258.9 ppm 104.2 ppm
EB0010561	103.7 ppm	EB0010579	52.94 ppm
EB0010681	52.22 ppm		Jz.J4 ppm
CRM Serial No.	Concentration	CRM Serial No.	Concentration
CC727481	800.0 ppm	CC727493	390.0 ppm
CC727496	253.0 ppm	CC727498	150.0 ppm
Analytical Method: NDIR			

Digitally signed by:Quality Control Reason:Dry gas standard certification of analysis Location:Airgas USA LLC (Lab) Date:04.05.2023 17:55

Approved for Release:

Nor Marsala

Rod Marsala

ISO 17025:2017 A2LA accredited. Certificate Number 3082.06 ISO 17034:2016 A2LA accredited. Certificate Number 3082.07

Page 1 of 1

