

MISSOURI DEPARTMENT OF HEALTH AND SENIOR SERVICES STATE PUBLIC HEALTH LABORATORY BREATH ALCOHOL PROGRAM

INTOX EC/IR II MAINTE			EPORT #3
Complete this report at the time of the regu			
days). Complete this report whenever the inst			
into service. Retain the original and send a INTOX EC/IR II SN NAME OF A		DATE OF INSPECTION	
INTOX EC/IR II SN NAME OF AU 12854 ST JOHN		05/15/2024	
LOCATION OF INSTRUMENT (STREET AND CITY)	FD	TIME OF INSPECTION	
8944 ST. CHARLES ROCK RD ST. JOHN		11:51 CDT	
CHECKLIST: Place a mark in the box by each it	tem if found to be satisfa-		
established limits. (Write in observed value)			
before using instrument.			
X DIAGNOSTIC RECORD			
X BLANK CHECK	X CO2 CHECK		
X FC 1 TEMP	X FLOW CHECK		w
X SRC TEMP	X FCB CHECK		
l territoria	X CRC COMP CH	POV	
X DET TEMP			
X BT TEMP	X CRC CAL CHE	CK	
X STD 2 TEMP	X PRINT TEST		
X ETH CHECK			
BREATH ANALYZER ACCURACY STANDARDS			***************************************
SIMULATOR SOLUTION	X COMPRESSED	ETHANOL-GAS MIXTURE	******
X STANDARD SUPPLIER INTOXIMETERS	LOT# AG329101	EXP. DATE 10/18/2025	
SIMULATOR TEMP (34°C +0.2°C)	SIM. SN	SIM. NIST EXP DATE	
	ודגאו משת תשסוו שם מס אודו	TWENTIME DEDORT	
X CALIBRATION CHECK - (ONLY ONE STANDAR			_
Run three tests using a standard solu	tion. All three tests i	must be within $\pm 5\%$ of the standard v	alue
Run three tests using a standard solution and must have a spread of .005 or less	tion. All three tests i	must be within $\pm 5\%$ of the standard v	alue ng
Run three tests using a standard solute and must have a spread of .005 or lest used.	tion. All three tests ness. Mark the box corresp	must be within $\pm 5\%$ of the standard vocation to the standard solution bei	alue ng
Run three tests using a standard solute and must have a spread of .005 or less used. 0.10% STANDARD - MUST READ BETWEEN	tion. All three tests rs. Mark the box corresponded on the corresponded to the corresp	must be within $\pm 5\%$ of the standard vocation to the standard solution bei	alue ng
Run three tests using a standard solute and must have a spread of .005 or less used. 0.10% STANDARD - MUST READ BETWEEN 0.08% STANDARD - MUST READ STANDARD - MUST	tion. All three tests not so that the box corresponds of the correspon	must be within $\pm 5\%$ of the standard vocation to the standard solution beingstones of the standard solution beingstones. The standard solution being the standard soluti	alue ng
Run three tests using a standard solute and must have a spread of .005 or less used. 0.10% STANDARD - MUST READ BETWEEN	tion. All three tests not so that the box corresponds of the correspon	must be within $\pm 5\%$ of the standard vocation to the standard solution beingstones of the standard solution beingstones. The standard solution being the standard soluti	alue ng
Run three tests using a standard solution and must have a spread of .005 or less used. 0.10% STANDARD - MUST READ BETWEEN 0.08% STANDARD - MUST READ BETWEEN 0.04% STANDARD - MUST READ STANDARD - MU	tion. All three tests not so that the box corresponds of the correspon	must be within $\pm 5\%$ of the standard vocation to the standard solution beingstones of the standard solution beingstones. The standard solution being the standard soluti	alue ng
Run three tests using a standard solution and must have a spread of .005 or less used. 0.10% STANDARD - MUST READ BETWEEN OF	tion. All three tests is s. Mark the box corresponded on the corre	must be within ±5% of the standard voonding to the standard solution beingstreet the standard votes and solution beingstreet the standard solution being	alue ng
Run three tests using a standard solution and must have a spread of .005 or less used. 0.10% STANDARD - MUST READ BETWEEN 0.08% STANDARD - MUST READ BETWEEN 0.04% STANDARD - MUST READ STANDARD - MU	tion. All three tests is s. Mark the box corresponded on the corre	must be within ±5% of the standard voonding to the standard solution beingstreet the standard votes and solution beingstreet the standard solution being	alue
Run three tests using a standard solution and must have a spread of .005 or less used. 0.10% STANDARD - MUST READ BETWEEN OF	tion. All three tests is a mark the box corresponds. Mark the box corresponds on the corr	must be within ±5% of the standard voonding to the standard solution beingstreet the standard votes and solution beingstreet the standard solution being	alue
Run three tests using a standard solution and must have a spread of .005 or less used. 0.10% STANDARD - MUST READ BETWEEN 0.08% STANDARD - MUST READ BETWEEN 0.04% STANDARD - MUST READ BETWEEN 0.04% STANDARD - MUST READ BETWEEN 0.04% STANDARD - MUST READ BETWEEN 0.05 TEST 1 0.080 g/210L TEST 2 TEST 1 0.080 g/210L TEST 2 TEST 1 0.080 g/210L TEST 2 TEST 1 TEST 1 0.080 g/210L TEST 2 TEST 2 TEST 1 0.080 g/210L TEST 2 TEST 2 TEST 1 TEST 2 TEST 3 TEST	tion. All three tests is a mark the box corresponds. Mark the box corresponds on the following control on the following ranges since the following ranges si	must be within ±5% of the standard vectoring to the standard solution being sive SIVE SIVE TEST 3 *** 0.080 g/210L NCE THE LAST MAINTENANCE REPORT: .1519 0 OVER .19 0	alue
Run three tests using a standard solution and must have a spread of .005 or less used. 0.10% STANDARD - MUST READ BETWEEN 0.08% STANDARD - MUST READ BETWEEN 0.04% STANDARD - MUST READ BETWEEN 0.04% STANDARD - MUST READ BETWEEN 0.004% STANDARD - MUST READ BETWEEN 0.000	tion. All three tests is a mark the box corresponds. Mark the box corresponds on the following control on the following ranges since the following ranges si	must be within ±5% of the standard vectoring to the standard solution being sive SIVE SIVE TEST 3 *** 0.080 g/210L NCE THE LAST MAINTENANCE REPORT: .1519 0 OVER .19 0	alue
Run three tests using a standard solution and must have a spread of .005 or less used. 0.10% STANDARD - MUST READ BETWEEN 0.08% STANDARD - MUST READ BETWEEN 0.04% STANDARD - MUST READ BETWEEN 0.04% STANDARD - MUST READ BETWEEN 0.04% STANDARD - MUST READ BETWEEN 0.05% STANDARD - MUST READ BETWEEN 0.05% STANDARD - MUST READ BETWEEN 0.05% STANDARD OF BREATH TESTS IN 1.00% STANDAR	tion. All three tests is a mark the box corresponds. Mark the box corresponds on the following control on the following ranges since the following ranges si	must be within ±5% of the standard vectoring to the standard solution being sive SIVE SIVE TEST 3 *** 0.080 g/210L NCE THE LAST MAINTENANCE REPORT: .1519 0 OVER .19 0	alue
Run three tests using a standard solution and must have a spread of .005 or less used. 0.10% STANDARD - MUST READ BETWEEN 0.08% STANDARD - MUST READ BETWEEN 0.04% STANDARD - MUST READ BETWEEN 0.04% STANDARD - MUST READ BETWEEN 0.04% STANDARD - MUST READ BETWEEN 0.05 TEST 1 0.080 g/210L TEST 2 TEST 1 0.080 g/210L TEST 2 TEST 1 0.080 g/210L TEST 2 TEST 1 TEST 1 0.080 g/210L TEST 2 TEST 2 TEST 1 0.080 g/210L TEST 2 TEST 2 TEST 1 TEST 2 TEST 3 TEST	tion. All three tests is a mark the box corresponds. Mark the box corresponds on the following control on the following ranges since the following ranges si	must be within ±5% of the standard vectoring to the standard solution being sive SIVE SIVE TEST 3 *** 0.080 g/210L NCE THE LAST MAINTENANCE REPORT: .1519 0 OVER .19 0	alue
Run three tests using a standard solution and must have a spread of .005 or less used. 0.10% STANDARD - MUST READ BETWEEN 0.08% STANDARD - MUST READ BETWEEN 0.04% STANDARD - MUST READ BETWEEN 0.04% STANDARD - MUST READ BETWEEN 0.04% STANDARD - MUST READ BETWEEN 0.05% STANDARD - MUST READ BETWEEN 0.05% STANDARD - MUST READ BETWEEN 0.05% STANDARD OF BREATH TESTS IN 1.00% STANDAR	tion. All three tests is a mark the box corresponds. Mark the box corresponds on the following control on the following ranges since the following ranges si	must be within ±5% of the standard vectoring to the standard solution being sive SIVE SIVE TEST 3 *** 0.080 g/210L NCE THE LAST MAINTENANCE REPORT: .1519 0 OVER .19 0	alue
Run three tests using a standard solution and must have a spread of .005 or less used. 0.10% STANDARD - MUST READ BETWEEN 0.08% STANDARD - MUST READ BETWEEN 0.04% STANDARD - MUST READ BETWEEN 0.04% STANDARD - MUST READ BETWEEN 0.04% STANDARD - MUST READ BETWEEN 0.05% STANDARD - MUST READ BETWEEN 0.05% STANDARD - MUST READ BETWEEN 0.05% STANDARD OF BREATH TESTS IN 1.00% STANDAR	tion. All three tests is a mark the box corresponds. Mark the box corresponds on the following control on the following ranges since the following ranges si	must be within ±5% of the standard vectoring to the standard solution being sive SIVE SIVE TEST 3 *** 0.080 g/210L NCE THE LAST MAINTENANCE REPORT: .1519 0 OVER .19 0	alue
Run three tests using a standard solution and must have a spread of .005 or less used. 0.10% STANDARD - MUST READ BETWEEN 0.08% STANDARD - MUST READ BETWEEN 0.04% STANDARD - MUST READ BETWEEN 0.04% STANDARD - MUST READ BETWEEN 0.04% STANDARD - MUST READ BETWEEN 0.05 TEST 1 0.080 g/210L TEST 2 INDICATE THE NUMBER OF BREATH TESTS IN TEST 2 INDICATE THE NUMBER OF BREATH TEST 2 INDICATE T	tion. All three tests is a mark the box corresponds. Mark the box corresponds on the following control on the following ranges since the following ranges si	must be within ±5% of the standard vectoring to the standard solution being sive SIVE SIVE TEST 3 *** 0.080 g/210L NCE THE LAST MAINTENANCE REPORT: .1519 0 OVER .19 0	alue
Run three tests using a standard solution and must have a spread of .005 or less used. 0.10% STANDARD - MUST READ BETWEEN 0.08% STANDARD - MUST READ BETWEEN 0.04% STANDARD - MUST READ BETWEEN 0.04% STANDARD - MUST READ BETWEEN 0.04% STANDARD - MUST READ BETWEEN 0.05% STANDARD - MUST READ BETWEEN 0.05% STANDARD - MUST READ BETWEEN 0.05% STANDARD OF BREATH TESTS IN 1.00% STANDAR	tion. All three tests is a mark the box corresponds. Mark the box corresponds on the following control on the following ranges since the following ranges si	must be within ±5% of the standard vectoring to the standard solution being sive SIVE SIVE TEST 3 *** 0.080 g/210L NCE THE LAST MAINTENANCE REPORT: .1519 0 OVER .19 0	alue
Run three tests using a standard solution and must have a spread of .005 or less used. 0.10% STANDARD - MUST READ BETWEEN 0.008% STANDARD - MUST READ BETWEEN 0.04% STANDARD - MUST READ BETWEEN 0.04% STANDARD - MUST READ BETWEEN 0.04% STANDARD - MUST READ BETWEEN 0.004% STANDARD - MUST READ BETWEEN 0.004% STANDARD - MUST READ BETWEEN 0.004% STANDARD OF BREATH TESTS IN 1.000 STANDARD NOT SENT FOR REPAIRS INSPECTING OFFICER SIGNATURE 1.000 STANDARD	tion. All three tests in s. Mark the box corresponds. Mark the box corresponds of the cor	must be within ±5% of the standard vectoring to the standard solution being sive SIVE TEST 3 ** 0.080 g/210L NCE THE LAST MAINTENANCE REPORT: .1519 0 OVER .19 0 RESTORE THE INSTRUMENT TO OPERATE	alue
Run three tests using a standard solution and must have a spread of .005 or less used. 0.10% STANDARD - MUST READ BETWEEN 0.08% STANDARD - MUST READ BETWEEN 0.04% STANDARD - MUST READ BETWEEN 0.05 TEST 1 0.080 g/210L TEST 2 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1	tion. All three tests in s. Mark the box corresponds. Mark the box corresponds of the corresponding control of the corresponding corresponds of the corresponding corresponds of the corresponding corresponds of the corresponding corresponds of the corresponds of the corresponding corresponding corresponds of the corresponding corresponding corresponds of the corresponding corresponds of the corresponding corresponds of the corresponding corresponds of the corresponding corresponds on the corresponding corresponding corresponds on the corresponds on the corresponding corresponds on the corresponds on the corresponding corresponds on the corresponding cor	must be within ±5% of the standard vectoring to the standard solution being sive SIVE TEST 3 ** 0.080 g/210L NCE THE LAST MAINTENANCE REPORT: .1519 0 OVER .19 0 RESTORE THE INSTRUMENT TO OPERATE	alue
Run three tests using a standard solution and must have a spread of .005 or less used. 0.10% STANDARD - MUST READ BETWEEN 0.008% STANDARD - MUST READ BETWEEN 0.04% STANDARD - MUST READ BETWEEN 0.05 OF BREATH TESTS IN STANDARD	tion. All three tests in s. Mark the box corresponds. Mark the box corresponds of the corresponding control of the corresponding corresponds of the corresponding corresponds of the corresponding corresponds of the corresponding corresponds of the corresponds of the corresponding corresponds of the corresponds of the corresponding corresponds of the corresponds of the corresponding corresponds of the corresponds of the corresponding corresponds of the corresponds of the corresponding corresponds of the corresponding corresponds of the corresponding corresponds of the corresponding corresp	must be within ±5% of the standard vectoring to the standard solution being sive SIVE TEST 3 ** 0.080 g/210L NCE THE LAST MAINTENANCE REPORT: .1519 0 OVER .19 0 RESTORE THE INSTRUMENT TO OPERATE	alue
Run three tests using a standard solution and must have a spread of .005 or less used. 0.10% STANDARD - MUST READ BETWEEN 0.08% STANDARD - MUST READ BETWEEN 0.04% STANDARD - MUST READ BETWEEN 0.05 TEST 1 0.080 g/210L TEST 2 INDICATE THE NUMBER OF BREATH TESTS IN 1.1 TEST 2 INDICATE THE NUMBER OF BREATH TESTS IN 1.2 TEST 2 INDICATE THE NUMBER OF BREATH TESTS IN 1.3 TEST 2 INDICATE THE NUMBER OF BREATH TEST 3 INDICATE THE NUMBER OF	tion. All three tests in s. Mark the box corresponds. Mark the box corresponds of the corresponding control of the corresponding corresponds of the corresponding corresponds of the corresponding corresponds of the corresponding corresponds of the corresponds of the corresponding corresponding corresponds of the corresponding corresponding corresponds of the corresponding corresponds of the corresponding corresponds of the corresponding corresponds of the corresponding corresponds on the corresponding corresponding corresponds on the corresponds on the corresponding corresponds on the corresponds on the corresponding corresponds on the corresponding cor	must be within ±5% of the standard vectoring to the standard solution being sive SIVE TEST 3 ** 0.080 g/210L NCE THE LAST MAINTENANCE REPORT: .1519 0 OVER .19 0 RESTORE THE INSTRUMENT TO OPERATE	alue
Run three tests using a standard solution and must have a spread of .005 or less used. 0.10% STANDARD - MUST READ BETWEEN 0.08% STANDARD - MUST READ BETWEEN 0.04% STANDARD - MUST READ BETWEEN 0.05 TEST 1 0.080 g/210L TEST 2 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1	print full name Talbott, Tin Print full name Talbott, Tin Print full name Talbott, Tin Telephone number (314) 427-87	must be within ±5% of the standard vectoring to the standard solution being sive SIVE TEST 3 ** 0.080 g/210L NCE THE LAST MAINTENANCE REPORT: .1519 0 OVER .19 0 RESTORE THE INSTRUMENT TO OPERATE	alue



Airgas USA LLC (LAB) 3500 Bernard Street St. Louis, Mo. 63103 Ph: (314) 533-3100

Fax: (314) 533-3100 Fax: (314) 533-7328

Certificate of Analysis

Customer Name Exclusive Supplier Intoximeters, Inc. 2081 Craig Road St. Louis, Mo 63146 Test Date: 18-Oct-2023

Lot # AG329101 Model 108

Exp Date 18-Oct-2025

Cyl. Type 108 Component

Certified Concentration

Ethanol

 0.080 ± 0.002 BrAC (208 ppm)

Nitrogen

Certification Traceable to N.I.S.T. RGM and to CRM Ethanol Standards:

RGM Serial No.	Concentration	RGM Serial No.	Concentration
EB0010581	391.8 ppm	EB0010603	392.5 ppm
EB0010570	259.8 ppm	EB0010559	258.9 ppm
EB0010285	209.0 ppm	EB0010562	104.2 ppm
EB0010561	103.7 ppm	EB0010579	52.94 ppm
EB0010681	52.22 ppm		

 CRM Serial No.
 Concentration
 CRM Serial No.
 Concentration

 CC727481
 800.0 ppm
 CC727493
 390.0 ppm

 CC727496
 253.0 ppm
 CC727498
 150.0 ppm

Analytical Method: NDIR

Digitally signed by:Quality Control Reason:Dry gas standard certification of analysis Location:Airgas USA LLC (Lab) Date:10.19.2023 17:39

Approved for Release:

Yusef Woods

ISO 17025:2017 A2LA accredited. Certificate Number 3082.06 ISO 17034:2016 A2LA accredited. Certificate Number 3082.07



STATE OF MISSOURI

DEPARTMENT OF HEALTH AND SENIOR SERVICES
BREATH ALCOHOL PROGRAM



PERMIT TYPE II

TINA TALBOTT

is hereby authorized to instruct and supervise operators, train instructors, inspect, calibrate, perform field service and repairs, and operate the following breath analyzer(s):

INTOX EC/IR II				
for the determination of the alcoholic content of blood from a sample of expired air. Permit issued under the provisions of sections 577.020 through 577.041, RSMo and 306.111 through 306.119 RSMo.				
DATE4/2/2024	DIRECTOR OF STATE PUBLIC HEALTH LABORATORY			
NUMBER 240086				
EXPIRES 4/2/2026	DIRECTOR OF DEPARTMENT OF HEALTH AND SERVICES			
IO 600 0774 (6 40)	I AB-4 (R6-10)			

